



North Carolina Department of Health and Human Services
Division of Health Service Regulation
Office of the Director

2701 Mail Service Center • Raleigh, North Carolina 27699-2701
<http://www.ncdhhs.gov/dhst/>

Pat McCrory
Governor

Aldona Z. Wos, M.D.
Ambassador (Ret.)
Secretary DHHS

Drexdal Pratt, Director
Phone: 919-855-3750
Fax: 919-733-2757

February 28, 2013

Catharine W. Cummer, Regulatory Counsel, Strategic Planning
Duke University Health System
3100 Tower Blvd, Suite 1300
Durham NC 27707

Exempt from Review - Replacement Equipment

Facility: Duke University Health System d/b/a Duke University Hospital
Project Description: Replace MRI scanner
County: Durham
FID #: 943138

Dear Ms. Cummer:

In response to your letter of February 14, 2013, the above referenced proposal is exempt from certificate of need review in accordance with N.C.G.S 131E-184(a)(7). Therefore, you may proceed to acquire, without a certificate of need, the Siemens Magnetom Aera 1.5T MRI scanner to replace the existing Siemens Avanto 1.5T MRI scanner [Serial # 25445]. This determination is based on your representations that the existing unit will be removed from North Carolina and will not be used again in the State without first obtaining a certificate of need. Further please be advised that as soon as the replacement equipment is acquired, you must provide the CON Section and the Medical Facilities Planning Section with the serial number of the new equipment to update the inventory, if not already provided.

It should be noted that this Agency's position is based solely on the facts represented by you and that any change in facts as represented would require further consideration by this Agency and a separate determination. If you have any questions concerning this matter, please feel free to contact this office.

Sincerely,

Michael J. McKillip
Project Analyst

Craig R. Smith, Chief
Certificate of Need Section





Catharine W. Cummer
Regulatory Counsel, Strategic Planning

Mike

February 14, 2013

Via Electronic Mail

Michael J. McKillip, Project Analyst
Certificate of Need Section
Division of Health Service Regulation
2704 Mail Service Center
Raleigh, NC 27699-2704

Re: Equipment Replacement Project at Duke University Hospital – MRI

Dear Mr. McKillip:

I am writing to provide you with prior written notice of the planned acquisition of a replacement MRI scanner at Duke University Hospital. The existing and replacement equipment both can provide MRI scanning services. A completed equipment comparison form and a quote for the new equipment are enclosed. The total extended selling price reflecting the fair market value of the equipment is \$1,675,381, plus a \$300,000 trade-in allowance. All installation costs for the replacement equipment in the new DUH Duke Medical Pavilion building (for which a CON was already issued) are included in the vendor quote. Duke is not currently intending on purchasing any of the training options (nor incurring any other costs to make this equipment operational), and understands that the total project cost must remain below \$2 million in order for the project to be exempt from certificate of need review. As a result, the total cost plus trade-in allowance of this equipment will be \$1,975,381. The existing equipment is currently in use at Duke, but upon placement of the replacement equipment into service it will be removed from service in the state unless its use elsewhere is first approved by the state.

It is our understanding that with this prior notice, this acquisition is exempt from certificate of need review as the acquisition of replacement equipment pursuant to N.C.G.S. Section 131E-184, and we accordingly intend to proceed immediately with this project. If you have any questions, please let me know as soon as possible.

Mr. Mike McKillip
February 14, 2013

Thank you for your attention to this matter. Should you have any questions, please let me know.

Very truly yours,

A handwritten signature in cursive script that reads "Catharine W. Cummer".

Catharine W. Cummer

Enclosures

cc: Monte Brown
Britt Crewse

EQUIPMENT COMPARISON

| | EXISTING EQUIPMENT | REPLACEMENT EQUIPMENT |
|--|--------------------|---|
| Type of Equipment (List Each Component) | MRI Scanner | MRI Scanner |
| Manufacturer of Equipment | Siemens | Siemens |
| Tesla Rating for MRIs | 1.5 | 1.5T |
| Model Number | Avanto 1.5 Tesla | MAGNETOM Aera |
| Serial Number | 25445 | N/A |
| Provider's Method of Identifying Equipment | MR 1 | N/A |
| Specify if Mobile or Fixed | Fixed | Fixed |
| Mobile Trailer Serial Number/VIN # | | |
| Mobile Tractor Serial Number/VIN # | | |
| Date of Acquisition of Each Component | 2005 | 2013 |
| Does Provider Hold Title to Equipment or Have a Capital Lease? | Title | Title |
| Specify if Equipment Was/Is New or Used When Acquired | New | New |
| Total Capital Cost of Project (Including Construction, etc.) | | \$1,675,381 |
| Total Cost of Equipment | \$1,330,000 | \$ 1,975,381 (including trade-in allowance) |
| Fair Market Value of Equipment | | \$1,975,381 |
| Net Purchase Price of Equipment | \$1,330,000 | \$1,675,381 |
| Locations Where Operated | DUH – Duke North | DUH – DMP |
| Number Days In Use/To be Used in N.C. Per Year | 365 | 365 |
| Percent of Change in Patient Charges (by Procedure) | NA | 0 |
| Percent of Change in Per Procedure Operating Expenses (by Procedure) | NA | 0 |
| Type of Procedures Currently Performed on Existing Equipment | MR Scanning | NA |
| Type of Procedures New Equipment is Capable of Performing | NA | MR Scanning |

SIEMENS

Siemens Medical Solutions USA, Inc.
51 Valley Stream Parkway, Malvern, PA 19355
Fax: (866) 309-6967

SIEMENS REPRESENTATIVE
Stuart Wadley - (919) 605-9227

Customer Number: 0000005901

Date: 2/14/2013

DUKE UNIVERSITY HEALTH SYSTEM-DUKE MEDICAL CENTER
Erwin Road
Durham, NC 27710

Siemens Medical Solutions USA, Inc. is pleased to submit the following quotation for the products and services described herein at the stated prices and terms, subject to your acceptance of the terms and conditions on the face and back hereof, and on any attachment hereto.

| <u>Table of Contents</u> | <u>Page</u> |
|---|-------------------|
| MAGNETOM Aera - USA | 3 |
| General Terms and Conditions..... | 11 |
| Warranty Information..... | 18 |
| Detailed Technical Specifications | 19 |
| Cut Sheets..... | following page 42 |

Proposal valid until 3/31/2013

For the items that are shipped early - The warranty for the Products will be twelve (12) months commencing on the earlier of first patient use or sixty (60) days from installation completion. Pricing is contingent on Customer purchasing POS Service Agreement at time of purchase.

This proposal includes the trade-in of equipment referenced in Trade Sheet Project # 2013-099

This Quotation includes the trade-in equipment described herein and referenced by the Project number identified herein, as further described in the associated Trade-In Sheet which is incorporated herein by reference. Purchaser certifies that the description of the trade-in equipment as set forth on the Trade-In Sheet is a true and accurate representation of the equipment, and that the equipment is in good working condition unless otherwise noted on the Trade-In Sheet. The trade-in equipment must be available for inspection within two weeks of the scheduled de-installation date.

The trade-in equipment must be made available for removal no later than turn over of the new equipment. In the event that access to the trade-in equipment is denied past 14 days post-turnover, then Purchaser shall pay to Seller a rental fee in the amount 10% of the total trade-in value including Elevate discount (no less than \$1000) for each month, or part thereof, that access is denied. If the de-installation of the trade-in equipment is delayed by Purchaser for reasons other than a force majeure event, or if upon inspection by Seller it is determined that the equipment does not meet the manufacturer's operating specifications, or if any items listed as included on the Trade-In Sheet are not made available at the time of de-installation, then trade-in value will be re-evaluated and any loss in value or additional costs incurred by Seller shall be deducted from the established trade-in value and the pricing set forth on this Quotation will be adjusted by change order.

Purchaser further acknowledges and agrees that (i) the trade-in equipment will be free and clear of all liens and encumbrances including, but not limited to, unpaid leases and loans, and that upon request, it will execute a bill of sale or other documents reasonably satisfactory to Siemens to transfer title and ownership of the equipment to Seller, (ii) it is Purchaser's sole responsibility to delete all protected health information and any other confidential information from the equipment prior to de-installation, without damaging or cannibalizing the equipment or otherwise affecting the operation of the equipment in accordance with its specifications, (iii) the equipment, including all updates, upgrades, modifications, enhancements, revisions, software, S/W disks and manuals, shall

SIEMENS

Siemens Medical Solutions USA, Inc.
51 Valley Stream Parkway, Malvern, PA 19355
Fax: (866) 309-6967

SIEMENS REPRESENTATIVE
Stuart Wadley - (919) 605-9227

be returned to Siemens in good operating condition, reasonable wear and tear excepted, and (iv) to the extent not prohibited by applicable law, Purchaser shall indemnify and hold Seller harmless from and against any and all claims, demands, causes of action, damages, liability, costs and expenses (including reasonable attorney's fees) resulting or arising from Purchaser's failure to comply with its obligations with respect to the trade-in equipment. Purchaser must provide a clear path for the removal of the trade-in equipment. Any additional costs due to the need to use a larger rig (other than a standard 80 ton rig), as well as any construction activities, street closings, permits, etc., required to deinstall/remove the equipment are out-of-scope costs and will be the responsibility of Purchaser. FOR MR SYSTEMS: cryogen levels must be least 75% upon time of deinstallation. FOR MOBILE SYSTEMS: system must be road worthy and a state issued title transferring ownership to Seller must be received by Seller prior to the removal of the mobile system. Purchaser must vacate the room of all items not listed on the Trade-In Sheet, or otherwise clearly identify all items listed on the Trade-In Sheet, prior to the start of the de-installation. If this is not done, Seller will have no liability for items which are subsequently removed or scrapped. In addition, if the purchase and installation of the new equipment covered by this Quotation is not completed, then Seller shall invoice Purchaser for all costs and expenses incurred by Seller in connection with the deinstallation and removal of the trade-in equipment, including but not limited to labor, materials, rigging out, and transportation, which costs shall be paid by Purchaser within thirty (30) days of the invoice date.

Accepted and Agreed to by:

Siemens Medical Solutions USA, Inc.

By (sign): _____
Name: Stuart Wadley
Title: Account Executive
Date: _____

**DUKE UNIVERSITY HEALTH SYSTEM-DUKE
MEDICAL CENTER**

By (sign): _____
Name: _____
Title: _____
Date: _____

All pages of the signed proposal must be returned to Siemens to process the order - Thank you.

SIEMENS

Siemens Medical Solutions USA, Inc.
51 Valley Stream Parkway, Malvern, PA 19355
Fax: (866) 309-6967

SIEMENS REPRESENTATIVE
Stuart Waddey - (919) 605-9227

Quote Nr: 1-4MQDWY Rev. 1

Terms of Payment: 10% Down, 80% Delivery, 10% Installation
Free On Board: Shipping Point

Purchasing Agreement: Not Applicable

MAGNETOM Aera - USA

All items listed below are included for this system: (See Detailed Technical Specifications at end of Proposal.)

| Qty | Part No. | Item Description | Extended Price |
|-----|----------|--|----------------|
| 1 | 14416900 | MAGNETOM Aera - System MAGNETOM Aera - 1.5T Tim+Dot system - The integration of the next generation Tim - "Tim 4G" and the Siemens unique Dot Engines (Day optimizing throughput Engine). Short and open appearance (145 cm system length with 70 cm Open Bore Design). Tim 4G's redesigned RF system and all-new coil architecture. - Siemens unique DirectRF(tm) technology enable Tim's new all digital-in/ digital-out design - All-new coil architecture including Dual-Density Signal Transfer Technology - Whole-body superconductive Zero Helium Boil-Off 1.5T magnet - TrueForm Magnet and Gradient Design - Actively Shielded water-cooled Siemens gradient system - Head/Neck 20 DirectConnect, Spine 32 DirectConnect, Body 18, Flex Large/Small 4 Dot offers patient personalization, user guidance and process automation that result in consistent examination results. - Brain Dot Engine is designed to simplify general brain examinations through personalized, guided and automated workflows. - Dot Display and Dot Control Centers - efficient patient preparation. Additional features include: -Tim Application Suite including Neuro, Angio, Cardiac, Body, Onco, Breast, Ortho, Pediatric and Scientific Suite - syngo MR software including 1D/2D PACE, syngo BLADE, iPAT ² , Phoenix, Inline Technologies. - High performance host computer and measurement and reconstruction system. The system (magnet, electronics and control room) can be installed in 30sqm space. For system cooling either the Eco Chiller options or the Separator is required. | \$1,163,730 |
| 1 | 14416902 | Tim [204x48] XQ Gradients #Ae Tim [204x48] XQ-gradients performance level Tim 4G with it's newly designed RF system and innovative coil architecture enables high resolution imaging and increased throughput. Up to 204 simultaneously connected coil elements in combination with the standard 48 independent RF channels, allow for more flexible parallel imaging. Maximum SNR through the new Tim 4G matrix coil technology. This option includes also Advanced High Order Shim. XQ - gradients The XQ- gradients are designed combining high performance and linearity to support clinical whole body imaging at 1.5T. The force compensated gradient system minimizes vibration levels and acoustic noise. The XQ gradients combine 45 mT/m peak amplitude with a slew rate of 200 T/m/s. | \$77,748 |
| 1 | 14416906 | Tim Dockable Table #Ae The new Tim Dockable Table is designed for maximum patient comfort and smooth patient preparation. Tim Dockable Table can support up to 250 kg (550 lbs) patients without restricting the vertical or horizontal movement. The one step docking mechanism and the innovative multi-directional navigation wheel ensure easy maneuvering and handling. Critically ill or immobile patients can now be prepared outside the examination room for maximum patient care, flexibility and speed. | \$46,649 |
| 1 | 08464872 | PC Keyboard US english #Tim Standard PC keyboard with 101 keys. | \$1 |

SIEMENS

Siemens Medical Solutions USA, Inc.
51 Valley Stream Parkway, Malvern, PA 19355
Fax: (866) 309-6967

SIEMENS REPRESENTATIVE
Stuart Waddey - (919) 605-9227

| Qty | Part No. | Item Description | Extended Price |
|-----|---------------------|---|----------------|
| 1 | 14416914 | Pure White Design #T+D The MAGNETOM Aera / MAGNETOM Skyra design is available in different light and appealing variants which perfectly integrates into the different environments. The color of the main face plate cover of the Pure White Design Variant with the integrated Dot Control Centers and the unique Dot Display is brilliant white surrounded by a brilliant silver trim. The asymmetrical deco area on the left side is colored white matte and also with a brilliant surrounding silver trim. The table cover is presented also in the same color and material selection. | \$1 |
| 1 | MR_TRADE_IN_ALLOW | MR Trade-in-Allowance Avanto Serial Nbr 25445; Proj Nbr 2013-099; De-Install Dte 9/2013; Expiration Dte 5/22/2013 (\$300,000.00) | -\$300,000 |
| 1 | 08857828 | UPS Cable #Tim Power cable for connecting the UPS Powerware PW 9130-3000i (14413662) to the ACC of MAGNETOM Tim and MAGNETOM Tim+Dot systems for backing up the computer. Standard cable length: 9 m. | \$1,500 |
| 1 | 14413662 | UPS Powerware PW9130G-3000T-XLEU UPS system Eaton PW9130G-3000T-XLEU for MAGNETOM Tim, MAGNETOM Tim+Dot and MAGNETOM Symphony systems for safeguarding computers. Power output: 3.0 kVA / 2.7 kW Bridge time: 5 min full load / 14 min half load Input voltage: 230 VAC | \$2,700 |
| 1 | 14413663 | UPS Battery module UPS battery module Eaton PW 9130N-3000T-EBM for all MAGNETOM Tim, MAGNETOM Tim+Dot and MAGNETOM Symphony systems for safeguarding computers. Extension for: PW9130I-3000T Battery type: Closed, maintenance-free Extension of the bridge time to: 24 minutes with a module Dimensions (H x W x D): Battery module: 346 x 214 x 412 mm incl. bracket set Weight: approx. 50 kg | \$1,000 |
| 1 | 14407258 | MR Workplace Table 1.2m Table suited for syngo Acquisition Workplace and syngo MR Workplace based on syngo Hardware. | \$995 |
| 1 | 14407261 | MR Workplace Container, 50cm 50 cm wide extra case for the syngo host computer with sliding front door to allow change of storage media (CD/DVD/USB). | \$1,244 |
| 1 | 4MR5142869 | Armrest #MR | \$240 |
| 1 | CHILINST_AVT | Chiller Start-up and Warranty for TIM | \$3,750 |
| 1 | MR_ADDL_RIGGING | Additional Rigging MR \$18,125 | \$18,125 |
| 1 | MRWSE | MR Wall sign -English | \$41 |
| 1 | MR_PREINST_DOCK | T+D Preinstall kit for dockable table | \$550 |
| 1 | MR_STD_RIGGING_INST | MR Standard Rigging and Installation MR Standard Rigging and Installation This quotation includes standard rigging and installation of your new MAGNETOM system Standard rigging into a room on ground floor level of the building during standard working hours (Mon. - Fri./ 8 a.m. to 5 p.m.) It remains the responsibility of the Customer to prepare the room in accordance with the SIEMENS planning documents Any rigging requiring a crane over 80 tons and/or special site requirements (e.g. removal of existing systems, etc.) is an incremental cost and the responsibility of the Customer. All other "out of scope" charges (not covered by the standard rigging and installation) will be identified during the site assessment and remain the responsibility of the Customer. | \$0 |
| 1 | MR_BTL_INST_ALL | MR Standard Rigging & Install | \$15,000 |
| 1 | MR_CRYO | Standard Cryogens | \$8,000 |

SIEMENS

Siemens Medical Solutions USA, Inc.
51 Valley Stream Parkway, Malvern, PA 19355
Fax: (866) 309-6967

SIEMENS REPRESENTATIVE
Stuart Wadley - (919) 605-9227

| Qty | Part No. | Item Description | Extended Price |
|-----|----------|--|----------------|
| 1 | MR_PM | MR Project Management A Siemens Project Manager (PM) will be the single point of contact for the implementation of your Siemens equipment. The assigned PM will work with the customer's facilities management, architect or building contractor to assist you in ensuring that your site is ready for installation. Your PM will provide initial and final drawings and will coordinate the scheduling of the equipment, installation, and rigging, as well as the initiation of on-site clinical education. | \$0 |
| 1 | 14402592 | Inline Composing syngo #Tim Automatic anatomical or angiographic composing of multiple adjacent coronal or sagittal images for presentation and further evaluation. Composed images can be automatically loaded into Graphical Slice Positioning for scan planning purposes. | \$15,550 |
| 1 | 14402593 | Tim Planning Suite With the Tim Planning Suite, multiple regions in the entire body can be examined in a minimum of time through measurement planning on a single FoV of any desired size. | \$18,659 |
| 1 | 14416923 | Abdomen Dot Engine #T+D The Abdomen Dot Engine: Personalized Exam Strategies - Guidance - Automatic sequence scaling - Auto Navigator - Auto-FoV - Timeline setup and monitoring - Automatic Voice Commands - Auto Bolus Detection - Inline radial range calculation for MRCP - Inline Subtraction - Inline Registration | \$37,319 |
| 1 | 08464740 | Flow Quantification #Tim Special sequences for quantitative assessment of flow. | \$12,440 |
| 1 | 07365419 | Argus Flow | \$9,330 |
| 1 | 14418563 | Neuro Perfusion Evaluation,USA #T+D Neuro Perfusion Evaluation syngo provides a task card for detailed post-processing of brain perfusion data sets. Color display of the relative Mean Transit Time (relMTT), relative Cerebral Blood Volume (relCBV), corrected rel CBV, and relative Cerebral Blood Flow (relCBF) is supported. Flexible selection of the Arterial Input Function (AIF). Furthermore a calculation of maps using the pre-selected local Arterial Input Functions (AIF) is provided. The detailed evaluation of brain perfusion data sets generates parameter maps for TTP and PBP and for the hemodynamic parameters relMTT, relCBV, rel CBVcor and relCBF. | \$9,330 |
| 1 | 14416944 | DTI Package #T+D The DTI Package is a bundle of: - Diffusion Tensor Imaging - DTI Evaluation and - DTI Tractography syngo The bundle comprehends all acquisition and postprocessing tools for comprehensive DTI exams. | \$38,376 |
| 1 | 14409110 | Arterial Spin Labeling 2D ASL is a non contrast enhanced brain perfusion technique. EPI sequence enhanced for PASL (Pulsed Arterial Spin Labeling) with preparation module (inversion pulse, saturation pulses) and selectable prospective motion correction. Perfusion-weighted color maps and relative cerebral blood flow (relCBF) color maps are calculated with Inline technology. | \$34,209 |
| 1 | 14402527 | SWI #Tim Susceptibility Weighted Imaging is a high-resolution 3D imaging technique for the brain with ultra-high sensitivity for microscopic magnetic field inhomogeneities caused by deoxygenated blood, products of blood decomposition and microscopic iron deposits. Among other things, the method allows for the highly sensitive proof of cerebral hemorrhages and the high-resolution display of venous cerebral blood vessels. | \$15,550 |
| 1 | 14416941 | Spectroscopy Package #T+D The Spectroscopy Package is a comprehensive software package which bundles Single Voxel Spectroscopy, 2D Chemical Shift Imaging, 3D Chemical Shift Imaging and syngo Spectroscopy Evaluation. Sequences and protocols for proton spectroscopy, 2D and 3D proton chemical shift imaging (2D CSI and 3D CSI) to examine metabolic changes in the brain (e.g. in tumors and degenerative diseases) and in the prostate are included. Furthermore included is the comprehensive syngo Spectroscopy Evaluation Software which enables fast evaluation of spectroscopy data on the syngo Acquisition Workplace. | \$31,099 |

SIEMENS

Siemens Medical Solutions USA, Inc.
51 Valley Stream Parkway, Malvern, PA 19355
Fax: (866) 309-6967

SIEMENS REPRESENTATIVE
Stuart Waddey - (919) 605-9227

| Qty | Part No. | Item Description | Extended Price |
|-----|----------|--|----------------|
| 1 | 14405328 | TWIST syngo #Tim This package contains a Siemens unique sequence and protocols for time-resolved (4D) MR angiographic and dynamic imaging in general with high spatial and temporal resolution. syngo TWIST supports comprehensive dynamic MR angio exams in all body regions. It offers temporal information of vessel filling in addition to conventional static MR angiography, which can be beneficial in detecting or evaluating malformations such as shunts. In case of general dynamic imaging, for example an increase in spatial resolution by a factor of up to 2 at 60 seconds temporal resolution (compared to conventional dynamic imaging) is possible due to intelligent k-space sampling strategies. Alternatively, increased temporal resolution at constant spatial resolution is possible. | \$18,659 |
| 1 | 14416955 | Body 18 #Ae The new Tim 4G coil technology with Dual Density Signal Transfer and SlideConnect Technology combines key imaging benefits: excellent image quality, high patient comfort, and unmatched flexibility: - 18 channels (inherent) or up to 30 (in combination with the Spine 32) - Dual Density Signal Transfer - Ultra light-weight - SlideConnect Technology The Body 18 is part of the standard configuration. The 18-channel coil with its 18 intergrated pre-amplifiers ensures excellent signal-to-noise ratio. The 18 coil elements provide extensive coverage in all directions. The single SlideConnect plug allows for fast and easy patient preparation. The light-weight coil ensures highest patient comfort. The Body 18 Coil features: - 18-element design with 18 integrated preamplifiers (3 clusters of 6 elements each) - Operates in an integrated fashion with the Spine 32 as an 30 channel body coil - Can be combined with further Body 18 coils for larger coverage - Can be positioned in different orientations (0°, 90°, 180°, 270°) for patient specific adaptations - No c oil tuning - iPAT compatible in all directions The highly flexible design enables a wide variety of applications including: - Thorax (incl. heart) - Abdomen - Pelvis - Hip Typically combined with: - Head / Neck 20 - Spine 32 - Additional Body 18 coil(s) (optional) - Peripheral Angio 36 (optional) - Flex Large 4 - Flex Small 4 - Loop coils (optional) - Endorectal coil (optional) | \$46,649 |
| 1 | 14405341 | Mapit syngo #Tim Based on the T1, T2 or T2* properties of the cartilage syngo ParametricMap allows the early detection of osteoarthritic break down of cartilage structures even before morphological changes occur. The method supports therapeutic decisions in individual patients and can be used to control treatments non-invasively, replacing surgeries or biopsies. The assessment of T1, T2 and T2* properties of tissues in other body regions is also possible. syngo ParametricMap provides very fast 2D and 3D high resolution imaging sequences and the Inline calculation of parametric maps for the T1, T2 and T2* properties of the imaged tissue. | \$7,277 |
| 1 | 07820082 | Inline Perfusion #Tim Automatic real-time calculation of Global Bolus Plot (GBP), Percentage of Baseline at Peak map (BBP), and Time-to-Peak map (TTP) with Inline technology. | \$9,330 |
| 1 | 14416958 | Peripheral Angio 36 #Ae The new Tim 4G coil technology with Dual Density Signal Transfer and SlideConnect Technology combines key imaging benefits: excellent image quality, high patient comfort, and unmatched flexibility: - 36 channels - Dual Density Signal Transfer - Ultra light-weight - SlideConnect Technology The 36-channel coil includes 36 integrated pre-amplifiers for excellent signal-to-noise ratio. The single SlideConect Plug allows for fast and easy patient preparation. The Peripheral Angio 36 features: - 36-element design with 36 integrated preamplifiers, distributed over 6 planes with 6 elements each - Operates in an integrated fashion with Body 18 coils and with the Spine 32 . For Whole-Body examinations also with the Head/ Neck 20 - Automatic table feed and active coil switch - Can be utilized head and feet first - Both legs are independently covered with coil elements, maximizing the coil filling factor and the signal-to-noise ratio - No coil tuning - iPAT-compatible - Dual-Density Signal Transfer enables ultra-high density coil designs by integrating key RF components into the local coil - SlideConnect technology for easy coil set up - One cable only for easy handling - Includes special non-ferromagnetic coil cart for safe, user-friendly storage Applications: - High-resolution angiography of both legs incl. Pelvis (by additional use of the Body 18) with highest signal-to-noise ratio - Visualization of the iliac arteries and aorta in combinaton with Body 18 - Bilateral examinations of long bones of the legs Typically combined with: Head/ Neck 20, Body 18, Spine 32, and all flexible coils such as Flex Large 4 or Flex Small 4 | \$37,319 |

SIEMENS

Siemens Medical Solutions USA, Inc.
51 Valley Stream Parkway, Malvern, PA 19355
Fax: (866) 309-6967

SIEMENS REPRESENTATIVE
Stuart Wadley - (919) 605-9227

| Qty | Part No. | Item Description | Extended Price |
|-----|----------|--|----------------|
| 1 | 14416961 | Hand/Wrist 16 #Ae The new Tim 4G coil technology with Dual Density Signal Transfer and SlideConnect Technology combines key imaging benefits: excellent image quality, high patient comfort, and unmatched flexibility. Hand/Wrist 16 for examinations of the left or right hand and wrist region consists of a base plate and an iPAT compatible 16-channel coil and allows high resolution imaging of the wrist and the hand within one examination. Hand/Wrist 16 will be connected via a SlideConnect plug for fast and easy patient preparation. | \$40,429 |
| 1 | 14416962 | Foot/Ankle 16 #Ae The new Tim 4G coil technology with Dual Density Signal Transfer and DirectConnect Technology combines key imaging benefits: excellent image quality, high patient comfort, and unmatched flexibility. Foot/Ankle 16 for examinations of the left or right foot and ankle region consists of a base plate and an iPAT compatible 16-channel coil and allows high resolution imaging of the foot and ankle within one examination. Foot/Ankle 16 is a cable-less coil and will be connected via DirectConnect for fast and easy patient preparation. | \$46,649 |
| 1 | 14407269 | Loop Coil, large #Tim No-tune loop receive coil for examinations of the upper or lower extremities (e.g. shoulder). | \$2,488 |
| 1 | 14407270 | Loop Coil, medium #Tim No-tune loop receive coil for examinations of small structures near the surface (e.g. hands, fingers, skin). | \$2,488 |
| 1 | 14407271 | Loop Coil, small #Tim No-tune loop receive coil for examinations of small structures near the surface (e.g. hands, fingers, skin). | \$2,488 |
| 1 | 14405224 | Composing syngo #Tim This application provides dedicated evaluation software for creation of full-format images from overlapping MR volume data sets and MIPs (starting from syngo MR B13) acquired at multiple stages. | \$9,330 |
| 1 | 14409198 | Native syngo #Tim Integrated software package with sequences and protocols for non-contrast enhanced 3D MRA with high spatial resolution. syngo NATIVE particularly enables imaging of abdominal and peripheral vessels and is an alternative to MR angiography techniques with contrast medium, especially for patients with severe renal insufficiency. | \$31,099 |
| 1 | 14430403 | Tx/Rx 15-channel Knee Coil DDST #Ae New 15-channel transmitter/receiver coil for joint examinations in the area of the lower extremities. Main features : - 15-element design (3x5 coil elements) with 15 integrated preamplifiers, - iPAT-compatible - SlideConnect Technology | \$44,783 |
| 1 | 14416968 | CP Extremity Coil #Ae Circularly Polarized no-tune transmit/receive coil for joint examinations in the region of the lower extremities. | \$18,659 |
| 1 | 14416960 | Shoulder 16 Coil Kit #Ae The new Tim 4G coil technology with Dual Density Signal Transfer and SlideConnect Technology combines key imaging benefits: excellent image quality, high patient comfort, and unmatched flexibility. The Shoulder 16 Coil Kit for examinations of the left or right shoulder consists of a base plate and two different sized iPAT compatible 16 channel coils (Shoulder Large 16 and Shoulder Small 16). These will be attached and can be relocated on the base plate. The 16-element coils with 16 integrated pre-amplifiers ensure maximum signal-to-noise ratio. Shoulder Large 16 and Shoulder Small 16 will be connected via a SlideConnect plug for fast and easy coil set-up and patient preparation. | \$37,319 |
| 1 | 14416952 | Coil Storage Cart #T+D Specially designed non-ferromagnetic cart for easy storage of the most commonly used coils and accessories. | \$3,110 |

SIEMENS

Siemens Medical Solutions USA, Inc.
 51 Valley Stream Parkway, Malvern, PA 19355
 Fax: (866) 309-6967

SIEMENS REPRESENTATIVE
 Stuart Waddey - (919) 605-9227

| Qty | Part No. | Item Description | Extended Price |
|----------------------|----------|---|--------------------|
| 1 | 14416975 | Eco Chiller 60kW The KKT ECO 133 -L chiller is a dedicated 20°C cooling system for MAGNETOM Aera and MAGNETOM Skyra which automatically adapts to the different cooling requirements (e.g. system in operation, standby, ..) to reduce the energy consumption for cooling. The cooling system must be used in combination with the IFP (Interface Panel), if there is no on-site chilled water supply at all. The IFP is included in the scope of supply. | \$38,000 |
| 1 | 14418521 | syngo Expert-i #T+D This software application enables remote access to the system (connected via local area network) for planning and processing. | \$6,842 |
| 1 | 14416969 | Flex Coil Interface 1.5T Interface with integrated preamplifiers for the connection of the following flexible coils to the system: - Flex Large 4 - Flex Small 4 - Loop Coil, large - Loop Coil, medium - Loop Coil, small - Endorectal Coil | \$9,330 |
| System Total: | | | \$1,675,381 |

OPTIONS:

| Qty | Part No. | Item Description | Extended Price | Initial to Accept |
|-----|---------------------|---|----------------|-------------------|
| 1 | MR_INITIAL_32 | Initial onsite training 32 hrs MR_INITIAL_32 Up to (32) hours of on-site clinical education training, scheduled consecutively (Monday - Friday) during standard business hours for a maximum of (4) imaging professionals. Training will cover agenda items on the ASRT approved checklist. Uptime Clinical Education phone support is provided during the warranty period for specified posted hours. This educational offering must be completed (12) months from install end date. If training is not completed within the applicable time period, Siemens obligation to provide the training will expire without refund. | + \$7,400 | X _____ |
| 1 | MR_FOLLOWU P_32 | Follow-up training 32 hrs Up to (32) hours of follow-up on-site clinical education training, scheduled consecutively (Monday - Friday) during standard business hours for a maximum of (4) imaging professionals. Uptime Clinical Education phone support is provided during the warranty period for specified posted hours. This educational offering must be completed (12) months from install end date. If training is not completed within the applicable time period, Siemens obligation to provide the training will expire without refund. | + \$7,500 | X _____ |
| 1 | MR_INT_DOT_ BCLS | MR Dot Training Class Tuition for (1) imaging professional to attend Classroom Course at Siemens Training Center. The objectives of this class are to introduce the user interface of the common syngo platform, including Dot, and instructions on building protocols, demonstration of software functions, and hands-on sessions. This class includes lunch, economy airfare, and lodging for (1) imaging professional. All arrangements must be arranged through Siemens designated travel agency. This educational offering must be completed (12) months from install end date. If training is not completed within the applicable time period, Siemens obligation to provide the training will expire without refund. | + \$4,200 | X _____ |
| 1 | MR_ADD_32 | Additional onsite training 32 hours Up to (32) hours of on-site clinical education training, scheduled consecutively (Monday - Friday) during standard business hours for a maximum of (4) imaging professionals. Training will cover agenda items on the ASRT approved checklist if applicable. This educational offering must be completed (12) months from install end date. If training is not completed within the applicable time period, Siemens obligation to provide the training will expire without refund. | + \$7,400 | X _____ |

SIEMENS

Siemens Medical Solutions USA, Inc.
51 Valley Stream Parkway, Malvern, PA 19355
Fax: (866) 309-6967

SIEMENS REPRESENTATIVE
Stuart Waddey - (919) 605-9227

| Qty | Part No. | Item Description | Extended Price | Initial to Accept |
|-----|------------------|--|----------------|-------------------|
| 1 | M3SSMREPIC BC | Spectris Solaris EP Injector ICBC Includes Spectris Solaris EP injector and Integrated Continuous Battery Charger (ICBC). - Optimized color touch screen with few keystrokes. - Six user-programmable phases for added flexibility. - Independent Keep Vein Open (KVO) allows more time to focus on patient. - Large 115 mL saline syringe allows for longer KVO and multiple flushes. - Design of low pressure tubing eliminates dead space in the "T" connection that can waste contrast. - The clear barrel design with molded FluidDots help detect the presence of air in a syringe. - Pressure Limit Setting control software enables user to select from one to six preset maximum pressure limits, ranging from 100-300 psi, and to view current pressure during injection next to the pre-selected maximum value on the Solaris display. Installation, applications and one year warranty provided by Medrad. Not for mobile use, refer to Siemens part number M3SSMR300EPM for the Solaris injector used in a mobile environment. This product has been tested and verified for compatibility with the following Siemens' products: MAGNETOM Trio, Espree, Essenza, Verio, Avanto, Symphony, Aera, Skyra and Biograph mMR. Compatibility with other products cannot be guaranteed and use with any other products may void service contracts and/or system warranties. | + \$55,646 | X _____ |
| 1 | 14413612 | Tissue 4D syngo #Tim Tissue 4D is an application for visualizing and post-processing dynamic contrast-enhanced 3D datasets. This card provides two evaluation options: - Standard curve evaluation - Curve evaluation according to a pharmacokinetic model | + \$35,000 | X _____ |
| 1 | 14416929 | Advanced Cardiac Package #T+D This package contains special sequences and protocols for advanced cardiac imaging including 3D and 4D syngo BEAT functionalities. It supports advanced techniques for ventricular function imaging, dynamic imaging, tissue characterization, coronary imaging, and more. | + \$50,000 | X _____ |
| 1 | 14407334 | Argus 4D Ventr.Function syngo #Tim syngo Argus 4D Ventricular Function software processes MR cine images of the heart and generates quantitative results for physicians in the diagnostic process. | + \$27,000 | X _____ |
| 1 | 14416908 | Tim Whole Body Suite #T+D Tim Whole Body Suite puts it all together. This suite enables table movement for imaging of up to 205 cm (6' 9") FoV without compromise. In combination with Tim's newly designed ultra highdensity array higher spatial and temporal resolution can be achieved along with unmatched flexibility of any coverage up to Whole Body. For faster exams and greater diagnostic confidence. | + \$35,000 | X _____ |
| 1 | 14426310 | Angio Dot Engine The timing of contrast injection and scan is widely considered the most challenging part of an angiographic exam. Angio Dot guides the user through angiographic single or multi station examinations by providing semi-automatic detection of arterial and venous timing windows using a test bolus technique. This information is fed back into the next planning steps automatically adapting scan parameters to the individual patient and patient's condition. Where needed, AutoVoiceCommands support the communication with the patient and ensure optimal timing of breathing, scanning and contrast media. All steps of contrast injection are presented in a simple, automated graphic on the monitor. syngo Inline Composing and the Tim Planning Suite are included. | + \$55,000 | X _____ |
| 1 | 07365484 | Image Fusion syngo This application provides a dedicated evaluation software for spatial alignment (matching) and visualization of image data either from different modalities (CT,MR,NM,PET) or from the same modality but from multiple examinations of the same patient. It supports optimal diagnostic outcome (fusion of morphological and functional information) and therapy planning. | + \$25,000 | X _____ |

SIEMENS

Siemens Medical Solutions USA, Inc.

51 Valley Stream Parkway, Malvern, PA 19355

Fax: (866) 309-6967

SIEMENS REPRESENTATIVE

Stuart Wadley - (919) 605-9227

FINANCING: The equipment listed above may be financed through Siemens. Ask us about our full range of financial products that can be tailored to meet your business and cash flow requirements. For further information, please contact your local Sales Representative.

ACCESSORIES: Don't forget to ask us about our line of OEM imaging accessories to complete your purchase. All accessories can be purchased or financed as part of this order. To purchase accessories directly or to receive our accessories catalog, please call us directly at 1-888-222-9944 ext. 7 or contact your local Sales Representative.

COMPLIANCE: Compliance with legal and internal regulations is an integral part of all business processes at Siemens. Possible infringements can be reported to our Helpdesk "Tell us" function at www.siemens.com/tell-us.

TERMS AND CONDITIONS OF SALE

1. GENERAL

1.1 Contract Terms. These terms and conditions constitute an integral part of any contract between the Seller identified on the first page hereof to sell products ("Products") and Purchaser and shall govern the sale of the Products. Seller shall not be bound by, and specifically objects to, any terms, conditions or other provisions which are different from or in addition to the provisions of this Agreement (whether or not it would materially alter this Agreement) which is proffered by Purchaser in any purchase order, receipt, acceptance, confirmation, correspondence or otherwise (even if provided to Seller concurrently with this Agreement), unless Seller specifically agrees to any such provision in a writing signed by Seller. Neither Seller's lack of objection to any such terms, nor delivery of the Products or provision of any services hereunder, shall constitute the agreement of Seller to any such terms. Purchaser acknowledges that this is a commercial and not a consumer transaction.

1.2 Acceptance. Purchaser shall be deemed to have assented to, and waived any objection to, this Agreement upon the earliest to occur of any of the following: Purchaser's completion or execution of this Agreement; Purchaser's acceptance of all or any part of the Products subject to this Agreement; Purchaser's issuance of a purchase order for any Products identified on Seller's quotation or proposal; or delivery of the Products to the common carrier for shipment pursuant hereto.

1.3 Refurbished/Used Products. For Products identified on Seller's Quotation as used or refurbished Products, these Products have been previously owned and used. When delivered to Purchaser, the Products may have received mechanical, electrical and/or cosmetic reconditioning, as needed, and will comply with the manufacturer's specifications. Since pre-owned Products may be offered simultaneously to several customers, the sale of such Products to Purchaser cannot be guaranteed and is subject to continuing availability at the time Purchaser accepts Seller's offer to sell the Products. If the Products are no longer available, Seller will use its best efforts to identify other products in its inventory that may be suitable for purchase by Purchaser, and if substitute products are not acceptable to Purchaser, then Seller will cancel the order and refund to Purchaser any deposits previously paid. The warranty period for any used or refurbished Products will be separately stated on the Quotation.

1.4 Third Party Products. If this Agreement includes the sale of third party products not manufactured by Seller as indicated on Seller's Quotation, then Purchaser agrees and acknowledges that (a) Purchaser has made the selection of these products on its own, (b) the products are being acquired by Seller solely at the request of and for the benefit of Purchaser, in order to eliminate the need for Purchaser to issue a separate purchase order to the manufacturer of the products, (c) no representation, warranty or guarantee has been made by Seller with respect to the products, (d) the obligation of Purchaser to pay Seller for the products is absolute and unconditional, (f) Purchaser will assert no claim whatsoever against the Seller with respect to the products, and will look solely to the manufacturer

regarding any such claims, (g) Purchaser will indemnify and hold Seller harmless from and against any and all claims, regardless of the form of action, related to, resulting from or caused by the products or any work or service provided by the manufacturer of the products or any other party, (h) use of the products may be subject to the Purchaser's agreement to comply with any software licensing terms imposed by the manufacturer, as well as any applicable laws, rule and regulations; and (i) the manufacturer, and not Seller, is solely responsible for any required installation, testing, validation, tracking, product recall, warranty service, maintenance, support, and complaint handling, as well as any other applicable FDA regulatory requirements.

2. PRICES

2.1 Quotations. Unless otherwise agreed to in writing or set forth in the quotation, all prices quoted by Seller are based on U.S. dollars, and include standard and customary packaging. F.O.B. terms are set forth in Section 6.2 hereof. Domestic prices apply only to purchasers located in, and who will use the Products in, the U.S. International prices apply to all purchasers located outside of, or who will use or ship or facilitate shipment of the Products outside of, the U.S. Unless otherwise stated, the quotation shall only be valid for forty-five (45) days from the date of the quotation.

2.2 Delay in Acceptance of Delivery. Should the agreed delivery date be postponed by Purchaser after the Seller's last call date after which the manufacture of the Products commences, Seller shall have the right to deliver to storage at Purchaser's risk and expense, and payments due upon delivery shall become due when Seller is ready to deliver.

3. TAXES

3.1 Seller agrees to honor any valid exemption certificate provided by Purchaser. Purchaser has a tax exemption certificate which it will provide to Seller prior to execution of this Agreement.

4. TERMS OF PAYMENT; DEFAULT

4.1 Payments; Due Date. Unless otherwise set forth in the quotation, Seller's payment terms are as follows: an initial deposit of 10% of the purchase price for each Product is due upon submission of the purchase order, an additional 80% of the purchase price is due upon delivery of each Product, and the final 10% of the purchase price is due upon completion of installation or when the Products are available for first patient use, whichever occurs first. Unless otherwise agreed, all undisputed payments other than the initial deposit are due thirty (30) days from receipt of invoice. All amounts payable pursuant to this Agreement are denominated in United States dollars, and Purchaser shall pay all such amount in lawful money of the United States. In the event that Purchaser makes any payments hereunder by credit card, Seller has the right to charge the Purchaser any credit card fees imposed on the Seller by the financial institution.

4.2 Late Payment. A service charge of 1½% per month, not to exceed the maximum rate allowed by law, shall be made on any portion of Purchaser's outstanding balance which is not paid within forty-five (45) days after invoice receipt date, which charge shall be determined and

compounded on a daily basis from the due date until the date paid. Payment of such service charge shall not excuse or cure Purchaser's breach or default for late payment. In addition, in the event that Purchaser fails to make any payment to Seller within this thirty 45 day period, including but not limited to any payment under any service contract, promissory note or other agreement with Seller, then Seller shall have no obligation to continue performance under any agreement with Purchaser.

4.3 Payment of Lesser Amount. If Purchaser pays, or Seller otherwise receives, a lesser amount than the full amount provided for under this Agreement, such payment or receipt shall not constitute or be construed other than as on account of the earliest amount due Seller. Seller may accept any check or payment in any amount without prejudice to Seller's right to recover the balance of the amount due or to pursue any other right or remedy. No endorsement or statement on any check or payment or in any letter accompanying a check or payment or elsewhere shall constitute or be construed as an accord or satisfaction.

4.4 Where Payment Due Upon Installation or Completion. Should any terms of payment provide for either full or partial payment upon installation or completion of installation or thereafter, and the installation or completion is delayed for any reason for which Seller is not responsible, then the Products shall be deemed installed upon delivery and, if no other terms were agreed upon in writing signed by the parties, the balance of payments shall be due no later than thirty (30) days from delivery regardless of the actual installation date.

4.5 Default; Termination. Each of the following shall constitute an event of default under this Agreement: (i) a failure by Purchaser to make any payment due Seller within ten (10) days of receipt of notice of non-payment from Seller; (ii) a failure by Purchaser to perform any other obligation under this Agreement within thirty (30) days of receipt of notice from Seller; (iii) a default by Purchaser or any affiliate of Purchaser under any other obligation to or agreement with Seller, Siemens Financial Services, Inc. or Siemens Medical Solutions Health Services Corporation, or any assignee of the foregoing (including, but not limited to, a promissory note, lease, rental agreement, license agreement or purchase contract); or (iv) the commencement of any insolvency, bankruptcy or similar proceedings by or against the Purchaser (including any assignment by Purchaser for the benefit of creditors). Upon the occurrence of any event of default, at Seller's election: (a) the entire amount of any indebtedness and obligation due Seller under this Agreement and interest thereon shall become immediately due and payable without notice, demand, or period of grace; (b) Seller may suspend the performance of any of Seller's obligations hereunder, including, but not limited to, obligations relating to delivery, installation and warranty services; (c) Purchaser shall put Seller in possession of the Products upon demand; (d) Seller may enter any premises where the Products are located and take possession of the Products without notice or demand and without legal proceedings; (e) at the request of Seller, Purchaser shall assemble the Products and make

TERMS AND CONDITIONS OF SALE

them available to Seller at a place designated by Seller which is reasonable and convenient to all parties; (f) Seller may sell or otherwise dispose of all or any part of the Products and apply the proceeds thereof against any indebtedness or obligation of Purchaser under this Agreement (Purchaser agrees that a period of 10 days from the time notice is sent to Purchaser shall be a reasonable period of notification of sale or other disposition of the Products by or for Seller); (g) if this Agreement or any indebtedness or obligation of Purchaser under this Agreement is referred to an attorney for collection or realization, Purchaser shall pay to Seller all costs of collection and realization (including, without limitation, a reasonable sum for attorneys' fees, expenses of title search, all court costs and other legal expenses) incurred thereby; and (h) Purchaser shall pay any deficiency remaining after collection of or realization by Seller on the Products. In addition, Seller may terminate this Agreement upon written notice to Purchaser in the event that Purchaser is not approved for credit or upon the occurrence of any material adverse change in the financial condition or business operations of Purchaser.

4.6 Financing. Notwithstanding any arrangement that Purchaser may make for the financing of the purchase price of the Products, the parties agree that any such financing arrangement shall have no effect on the Purchaser's payment obligations under this Agreement, including but not limited to Sections 4.1 and 4.2 above.

5. EXPORT TERMS

5.1 Unless other arrangements have been made, payment on export orders shall be made by irrevocable confirmed letter of credit, payable in U.S. dollars against Seller's invoice and standard shipping documents. Such letter of credit shall be in an amount equal to the full purchase price of the Products and shall be established in a U.S. bank acceptable to Seller. Purchaser shall procure all necessary permits and licenses for shipment and compliance with any governmental regulations concerning control of final destination of Products.

5.2 Purchaser shall not, directly or indirectly, violate any U.S. law, regulation or treaty, or any other international treaty or agreement, relating to the export or reexport of any Product or associated technical data, to which the U.S. adheres or with which the U.S. complies. Purchaser shall defend, indemnify and hold Seller harmless from any claim, damage, liability or expense (including but not limited to reasonable attorney's fees) arising out of or in connection with any violation of the preceding sentence. If Purchaser purchases a Product at the domestic price and exports such Product, or transfers such Product to a third party for export, outside of the U.S., Purchaser shall pay to Seller the difference between the domestic price and the international retail price of such Product pursuant to the payment terms set forth herein. Purchaser shall deliver to Seller, upon Seller's request, written assurance regarding compliance with this section in form and content acceptable to Seller.

6. DELIVERY, RISK OF LOSS

6.1 Delivery Date. Delivery and completion schedules are approximate only and are based on conditions at the time of acceptance of Purchaser's order by Seller. Seller shall make every

reasonable effort to meet the delivery date(s) quoted or acknowledged, but shall not be liable for any failure to meet such date(s). Partial shipments may be made.

6.2 Risk of Loss; Title Transfer. Unless otherwise agreed to in writing, the following shall apply:

(a) For Products that do not require installation by Seller or its authorized agent or subcontractor, and for options and add-on products purchased subsequent to delivery and installation of Products purchased under this Agreement, delivery shall be complete upon transfer of possession to common carrier, F.O.B. Destination whereupon title to and all risk of loss, damage to or destruction of the Products shall pass to Purchaser.

(b) For Products that require installation by Seller or its authorized agent or subcontractor, delivery shall be complete upon delivery of the Products to Purchaser's designated site, F.O.B. Destination; title to and all risk of loss, damage to or destruction of such Products shall pass to Purchaser upon completion of the installation by Seller or its authorized agent or subcontractor.

(c) All freight charges and other transportation, packing and insurance costs, license fees, custom duties and other similar charges shall be the sole responsibility of the Seller, unless included in the purchase price or otherwise agreed to in writing by Purchaser. In the event of any loss or damage to any of the Products during shipment, Seller and Purchaser shall cooperate in making a claim against the carrier.

7. SECURITY INTEREST/FILING

7.1 From the F.O.B. point, Seller shall have a purchase money security interest in the Products (and all accessories and replacements thereto and all proceeds thereof) until payment in full by Purchaser and satisfaction of all other obligations of Purchaser hereunder. Purchaser hereby (i) authorizes Seller to file (and Purchaser shall promptly execute, if requested by Seller) and any Uniform Commercial Code financing statements with respect to the Products and/or this Agreement. Purchaser also agrees that an original or a photocopy of this Agreement (including any addenda, attachments and amendments hereto) may be filed by Seller as a Uniform Commercial Code financing statement. Purchaser further represents and covenants that (a) it will keep the Products in good order and repair until the purchase price has been paid in full, (b) it will promptly pay all taxes and assessments upon the Products or the use thereof, (c) it will not attempt to transfer any interest in the Products until the purchase price has been paid in full, and (d) it is solvent and financially capable of paying the full purchase price for the Products.

8. CHANGES, CANCELLATION, AND RETURN

8.1 Orders accepted by Seller are not subject to change except upon written agreement.

8.2 Orders accepted by Seller are cancellable by Purchaser 30 days prior to shipment, upon Seller's written consent and payment by Purchaser of a cancellation charge equal to 10% of the price of the affected Products, plus any shipping, insurance, inspection and refurbishment charges; the cost of providing any training, education, site evaluation or other services; and any return,

cancellation or restocking fees with respect to any Third Party Products ordered by Seller on behalf of Purchaser. Seller may retain any payments received from Purchaser up to the amount of the cancellation charge. In no event can an order be cancelled by Purchaser or Products be returned to Seller after shipment has been made.

8.3 Seller shall have the right to change the manufacture and/or design of its Products if, in the judgment of Seller, such change does not alter the general function of the Products.

9. FORCE MAJEURE

9.1 Seller will make every effort to complete shipment, and installation where indicated, but shall not be liable for any loss or damage for delay in delivery, inability to install or any other failure to perform due to causes beyond its reasonable control including, but not limited to, acts of government or compliance with any governmental rules or regulations, acts of God or the public, war, civil commotion, blockades, embargoes, calamities, floods, fires, earthquakes, explosions, storms, strikes, lockouts, labor disputes, or unavailability of labor, raw materials, power or supplies. Should such a delay occur, Seller may reasonably extend delivery or production schedules or, at its option, cancel the order in whole or part without liability other than to return any unearned deposit or prepayment.

10. WARRANTY

10.1 Seller warrants that the Products manufactured by Seller and sold hereunder shall be free from defects in material or workmanship under normal use and service for the warranty period. The final assembled Products shall be new although they may include certain used, reworked or refurbished parts and components (e.g., circuit boards) that comply with performance and reliability specifications and controls Seller's obligation under this warranty is limited, at Seller's option, to the repair or replacement of the Product or any part thereof. Unless otherwise set forth in the Product Warranty attached hereto and incorporated herein by reference, the warranty period shall commence upon the earlier of the date that the Products have been installed in accordance with 12.6 hereof, which date shall be confirmed in writing by Seller, or first patient use, and shall continue for 12 consecutive months. Seller makes no warranty for any Products made by persons other than Seller or its affiliates, and Purchaser's sole warranty therefor, if any, is the original manufacturer's warranty, which Seller agrees to pass on to Purchaser, as applicable. The warranty provided by Seller under this Section 10 extends only to the original Purchaser, unless the Purchaser obtains the Seller's prior written consent with respect to any sale or other transfer of the Equipment during the term of the warranty.

10.2 No warranty extended by Seller shall apply to any Products which have been damaged by accident, misuse, abuse, negligence, improper application or alteration or by a force majeure occurrence as described in Section 9 hereof or by the Purchaser's failure to operate the Products in accordance with the manufacturer's instructions or to maintain the recommended operating environment and line conditions; which are defective due to unauthorized attempts to repair, relocate, maintain, service, add to or modify the

TERMS AND CONDITIONS OF SALE

Products by the Purchaser or any third party or due to the attachment and/or use of non-Seller supplied equipment without Seller's prior written approval; which failed due to causes from within non-Seller supplied equipment; which have been damaged from the use of operating supplies or consumable parts not approved by Seller. In addition, no warranty extended by Seller shall apply to any transducer failure due to events such as cracking from high impact drops, cable rupture from rolling equipment over the cable, or delamination from cleaning with inappropriate solutions. Seller's obligation under this warranty is limited to the repair or replacement, at Seller's option, of defective parts. Seller may effectuate such repair at Purchaser's facility, and Purchaser shall furnish Seller safe and sufficient access for such repair. Repair or replacement may be with parts or products that are new, used or refurbished. Repairs or replacements shall not interrupt, extend or prolong the term of the warranty. Purchaser shall, upon Seller's request, return the noncomplying Product or part to Seller with all transportation charges prepaid, but shall not return any Product or part to Seller without Seller's prior written authorization. -Purchaser shall pay Seller its normal charges for service and parts for any inspection, repair or replacement that is not, in Seller's sole judgment, required by noncompliance with the warranty set forth in Section 10.1.- Seller's warranty does not apply to consumable materials, disposables, supplies, accessories and collateral equipment, except as specifically stated in writing or as otherwise set forth in the Product Warranty attached hereto and incorporated herein by reference, nor to products or parts thereof supplied by Purchaser.

10.3 This warranty is made on condition that immediate written notice of any noncompliance be given to Seller and Seller's inspection reveals that the Purchaser's claim is valid under the terms of the warranty (i.e., that the noncompliance is due to traceable defects in original materials and/or workmanship).

10.4-Purchaser shall provide Seller with full and free access to the Products, network cabling and communication equipment as is reasonably necessary for Seller to provide warranty service. This access includes establishing and maintaining connectivity to the Products via VPN IPsec Tunneling (non-client) Peer-to-Peer connection, modem line, internet connection, broadband internet connection or other secure remote access reasonably required by Seller, in order for Seller to provide warranty service, including remote diagnostics, monitoring and repair services.

10.5 Warranty service will be provided without charge during Seller's regular working hours (8:30-5:00), Monday through Friday, except Seller's recognized holidays. If Purchaser requires that service be performed other than during these times, such service will be made available 24x7, at an additional charge, at Seller's then current rates. The obligations of Seller described in this section are Seller's only obligations and Purchaser's sole and exclusive remedy for a breach of product warranty.

10.6 SELLER MAKES NO WARRANTY OTHER THAN THE ONE SET FORTH HEREIN AND IN THE ATTACHED PRODUCT WARRANTY

COVERING THE APPLICABLE PRODUCT CATEGORY. SUCH WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED, INCLUDING BUT NOT LIMITED TO ANY EXPRESS OR IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR PARTICULAR PURPOSES, AND SUCH CONSTITUTES THE ONLY WARRANTY MADE WITH RESPECT TO THE PRODUCTS AND ANY DEFECT, DEFICIENCY OR NONCONFORMITY IN ANY PRODUCT, SERVICE OR OTHER ITEM FURNISHED UNDER THIS AGREEMENT.

10.7 In the event of any inconsistencies between the terms of this Section 10 and the terms of the attached Product Warranty, the terms of the attached Product Warranty shall prevail.

11. LIMITATION OF LIABILITY

11.1 In no event shall Seller's liability hereunder exceed the actual loss or damage sustained by Purchaser, up to the purchase price of the Products.

11.2 SELLER SHALL NOT BE LIABLE FOR ANY LOSS OF USE, REVENUE OR ANTICIPATED PROFITS, COST OF SUBSTITUTE PRODUCTS OR SERVICES, LOSS OF STORED, TRANSMITTED OR RECORDED DATA, OR FOR ANY INCIDENTAL, UNFORESEEN, SPECIAL, PUNITIVE OR CONSEQUENTIAL DAMAGES ARISING OUT OF OR IN CONNECTION WITH THIS AGREEMENT OR THE SALE OR USE OF THE PRODUCTS. This provision does not affect third party claims for personal injury arising as a result of Seller's negligence or product defect. **THE FOREGOING IS A SEPARATE, ESSENTIAL TERM OF THIS AGREEMENT AND SHALL BE EFFECTIVE UPON THE FAILURE OF ANY REMEDY, EXCLUSIVE OR NOT.**

12. INSTALLATION - ADDITIONAL CHARGES

12.1 General. Unless otherwise expressly stipulated in writing, the Products covered hereby shall be installed by and at the expense of Seller except that Seller shall not provide rigging or site preparation services unless otherwise agreed to in writing by Seller for an additional charge. Seller will not install accessory items such as cabinets, illuminators, darkroom equipment or processors for X-Ray and CT equipment, unless otherwise agreed to in writing by Seller.

12.2 Installation by Seller. If Seller specifies it will install the Products, the following applies: subject to fulfillment of the obligations set forth in 12.4 below, Seller shall install the Products covered hereby and connect same to the requisite safety switches and power lines to be installed by Purchaser. Except as otherwise specified below, if such installation and connection are performed by Seller technical personnel, prices shown include the cost thereof, provided that the installation and connection can be performed within the Continental United States or Puerto Rico and during normal business hours. Any overtime charges or other special expenses needs to be pre approved by Purchaser before work begins and shall be additional charges to the prices shown.

12.3 Trade Unions. In the event that a trade union, or unions, or other local labor conditions prevent Seller from performing the above work with its own employees or contractors then Purchaser shall either make all required arrangements with

the trade union, or unions, to permit Seller completion of said work or shall provide the personnel, at Purchaser's sole cost and expense. Moreover, any additional cost incurred by Seller and related to such labor disputes shall be paid by the Purchaser and Seller's obligations under such circumstances will be limited to providing engineering supervision of installation and connection of Seller equipment to existing wiring.

12.4 Purchaser's Obligations. Purchaser shall, at its expense, provide all proper and necessary labor and materials for plumbing service, carpentry work, conduit wiring, and other preparations required for such installation and connection. All such labor and materials shall be completed and available at the time of delivery of the Products by Seller. Additionally, the Purchaser shall provide free access to the premises of installation and, if necessary, safe and secure space thereon for storage of Products and equipment prior to installation by Seller. Purchaser shall be responsible, at its sole cost and expense, for obtaining all permits, licenses and approvals required by any federal, state or local authorities in connection with the installation and operation of the Products, including but not limited to any certificate of need and zoning variances. Purchaser shall provide a suitable environment for the Products and shall ensure, at its sole cost and expense, that its premises are free of asbestos, hazardous conditions and any concealed, unknown or dangerous conditions and that all site requirements are met. Seller shall delay its work until Purchaser has completed the removal of the asbestos or other hazardous materials or has taken any other precautions and completed any other work required by applicable regulations. Purchaser shall reimburse Seller for any increased costs and expenses incurred by Seller that are the result of or are caused by any such delay. In the event that Seller is requested to supervise the installation of the Products, it remains the Purchaser's responsibility to comply with local regulations. Seller is not an architect and all drawings furnished by Seller are not construction drawings.

12.5 Regulatory Reporting. In the event that any regulatory activity is performed by other than Seller authorized personnel, Purchaser shall be responsible for fulfilling any and all reporting requirements.

12.6 Completion of Installation. Installation shall be complete upon the conclusion of final calibration and checkout under Seller standard procedures to verify that the Products meet applicable written performance specifications. Notwithstanding the foregoing, first use of the Products by Purchaser, its agents or employees for any purpose after delivery shall constitute completion of installation.

13. PATENT, TRADEMARK AND OTHER INFRINGEMENT CLAIMS

13.1 Infringement by Seller. Seller warrants that the Products manufactured by Seller and sold hereunder do not infringe any U.S. patent or copyright. If Purchaser receives a claim that any such Product, or parts thereof, infringe upon the rights of others under any U.S. patent or copyright, Purchaser shall notify Seller immediately in writing. As to all infringement claims relating to Products or

TERMS AND CONDITIONS OF SALE

parts manufactured by Seller or one of its affiliates:

(a) Purchaser shall give Seller information, assistance and exclusive authority to evaluate, defend and settle such claims.

(b) Seller shall then, at its own expense, defend or settle such claims, procure for the Purchaser the right to use the Products, or remove or modify them to avoid infringement. If none of these alternatives is available on terms reasonable to Seller, then Purchaser shall return the Products to Seller and Seller shall refund to Purchaser the purchase price paid by the Purchaser less reasonable depreciation for Purchaser's use of the Products. The foregoing states Seller's entire obligation and liability, and the Purchaser's sole remedy, for claims of infringement.

14. INFRINGEMENT. Seller shall, at its own expense, indemnify, defend and hold harmless Purchaser against any claim, suit or proceeding brought against Purchaser which is based upon a claim, that the Product furnished under this Agreement infringes on any U.S. patent, and Seller shall pay all damages and expenses arising from such claim finally awarded by a court therein. If Product(s) supplied under this Agreement is held by a court to constitute infringement and its use is enjoined, Seller shall, at its own expense, and only with Purchaser's written approval, either procure for Purchaser the right to continue such use; or replace same with equal or better non-infringing goods; or modify it so it becomes non-infringing, provided such modification does not diminish the product's functionality. If Seller is unable to provide a non-infringing fully functional product, then Seller will provide a refund of the purchase price as depreciated on a straight-line, five (5) year basis.

Notwithstanding the foregoing, Seller shall have no liability under this provision to the extent that any such claim is based on (i) use of the product in combination with equipment or software not supplied hereunder where the product would not itself be infringing; (ii) compliance with designs, plans or specifications provided by any of the Indemnitees; (iii) use of the product in an application or environment for which it was not designed or not contemplated in the related documentation; (iv) use of other than the most recent release of the product provided by Seller; or (v) modification of the product by anyone other than Seller, its employees or agents. Notwithstanding the foregoing, in no event shall the Seller be liable for any loss of use, revenue or anticipated profits, loss of stored, transmitted or recorded data, or for any incidental, unforeseen, special, punitive or consequential damages.

15. DESIGNS AND TRADE SECRETS; LICENSE; CONFIDENTIALITY

15.1 Any drawings, data, designs, software programs or other technical information supplied by Seller to Purchaser in connection with the sale of the Products are not included in the sale of the Products to Purchaser, shall remain Seller's property and shall at all times be held in confidence by Purchaser. Such information shall not be reproduced or disclosed to others without Seller's prior written consent.

15.2 For all goods purchased hereunder which utilize software for their operation, such "Applications Software" shall be licensed to Purchaser under the terms of Seller's Software

License Schedule as attached hereto.

15.3 Diagnostic/Maintenance Software is not included under 14.2 above, is available only as a special option under a separate Diagnostic Materials License Agreement and may be subject to a separate licensing fee.

15.4 Seller and Purchaser shall maintain the confidentiality of any information provided or disclosed to the other party relating to the business, customers and/or patients of the disclosing party, as well as this Agreement and its terms (including the pricing and other financial terms under which the Purchaser will be purchasing the Products hereunder). Each party shall use reasonable care to protect the confidentiality of the information disclosed, but no less than the degree of care it would use to protect its own confidential information, and shall only disclose the other party's confidential information to its employees and agents having a need to know this information. The obligations of confidentiality set forth herein shall not apply to any information in the public domain at the time of disclosure or that is required to be disclosed by court order or by law.

16. ENGINEERING CHANGES

16.1 Seller makes no representation that engineering changes which may be announced in the future will be suitable for use on, or in connection with, the Products.

17. ASSIGNMENT

17.1 Neither party may assign any rights or obligations under this Agreement without the written consent of the other and any attempt to do so shall be void, except that Seller may assign this Agreement without consent to any subsidiary or affiliated company, and may delegate to authorized subcontracts or service suppliers any work to be performed under this Agreement so long as Seller remains liable for the performance of its obligations under this Agreement. This Agreement shall inure to and be binding upon the parties and their respective successors, permitted assigns and legal representatives. Seller shall have no obligations under this Agreement to any assignee of Purchaser that is not approved by Seller in advance.

18. DAMAGES, COSTS AND FEES

18.1 In the event that any dispute or difference is brought arising from or relating to this Agreement or the breach, termination or validity thereof, the prevailing party shall NOT be entitled to recover from the other party any punitive damages. The prevailing party shall be entitled to recover from the other party all reasonable attorneys' fees incurred, together with such other expenses, costs and disbursements as may be allowed by law.

19. MODIFICATION

19.1 This Agreement may not be changed, modified or amended except in writing signed by duly authorized representatives of the parties.

20. GOVERNING LAW; WAIVER OF JURY TRIAL

20.1 This Agreement shall be governed by the laws of the State of North Carolina.

20.2 EACH OF THE PARTIES EXPRESSLY WAIVES ALL RIGHTS TO A JURY TRIAL IN CONNECTION WITH ANY DISPUTE UNDER THIS AGREEMENT.

21. COST REPORTING

21.1 Purchaser agrees that it will fully and accurately account for and report in all cost reports

and otherwise fully and accurately disclose to federal and state health care program payors and fully and accurately reflect where and as appropriate to the applicable reimbursement methodology, all services and other items, including any and all discounts, received from Seller under this Agreement, in compliance with all applicable laws, rules and regulations, including but not limited to the Social Security Act and implementing regulations relating to Medicare, Medicaid and other federal and state health care reimbursement programs.

22. INTEGRATION

22.1 These terms and conditions, including any attachments or other documents incorporated by reference herein, constitute the entire agreement and the complete and exclusive statement of agreement with respect to the subject matter hereof, and supersede any and all prior agreements, understandings and communications between the parties with respect to the Products.

23. SEVERABILITY; HEADINGS

23.1 No provision of this Agreement which may be deemed unenforceable will in any way invalidate any other portion or provision of this Agreement. Section headings are for convenience only and will have no substantive effect.

24. WAIVER

24.1 No failure and no delay in exercising, on the part of any party, any right under this Agreement will operate as a waiver thereof, nor will any single or partial exercise of any right preclude the further exercise of any other right.

25. NOTICES

25.1 Any notice or other communication under this Agreement shall be deemed properly given if given in writing and delivered in person or mailed, properly addressed and stamped with the required postage, to the intended recipient at its address specified on the face hereof. Either party may from time to time change such address by giving the other party notice of such change in accordance with this section.

26. RIGHTS CUMULATIVE

26.1 The rights and remedies afforded to Seller under this Agreement are in addition to, and do not in anyway limit, any other rights or remedies afforded to Seller by any other agreement, by law or otherwise.

27. END USER CERTIFICATION

27.1 Purchaser represents, warrants and covenants that it is acquiring the Products for its own end use and not for reselling, leasing or transferring to a third party (except for lease-back financings).

28. COMPLIANCE PLAN

28.1 Seller acknowledges that it has been informed that Purchaser University Health System has implemented a compliance program for a number of purposes including, but not limited to, ensuring that the provision of, billing for, and care at Health System is in compliance with applicable Federal and State laws ("Compliance Plan") and all Health System operations are conducted in accordance with applicable laws and regulations. Seller acknowledges that it has received a copy of the PURCHASER Code of Conduct. Seller further states that 1) It has not been convicted of a criminal offense related to healthcare, 2) It is not currently under sanction, exclusion, or investigation

TERMS AND CONDITIONS OF SALE

(civil or criminal) by a federal or state enforcement, regulatory, administrative, or licensing agency or otherwise ineligible for federal or state program participation, and 3) It is not currently listed on the General Service Administration List of Parties Excluded from the Federal Procurement and Non-Procurement Programs. Seller will notify Purchaser in the event of any investigation, civil or criminal, of any of the above.

Seller will notify Purchaser promptly of any allegation of unethical behavior, wrongdoing, or violation of Purchaser Policies or procedures or federal or state laws or regulations.

Purchaser University Health System, upon Seller's request, shall make available such other information relating to its Compliance Program as is appropriate to assist Seller in adhering to the policies set forth in the Code of Conduct developed as part of the Compliance Plan. Seller is invited to attend educational sessions of Health System related to its Compliance Program. The failure of Seller to conduct its activities in accordance with the Code of Conduct and Compliance Program shall constitute a material breach of this Agreement and Health System shall have the absolute right to terminate this Agreement immediately. Seller agrees to verify that its employees and subcontractors meet all the requirements outlined above prior to their performing services under this Agreement.

29. PUBLICITY

Neither party shall in any way or in any form publicize or advertise in any manner the existence of this Agreement without the express prior written approval of the other party. However, nothing herein shall preclude Seller from listing Purchaser on its routine client list for matters of reference.

30. USE OF PURCHASER NAME, LOGOS, AND TRADEMARKS

Each parties' name and logos are the exclusive property of that party and may only be used for approved purposes. Each party agrees not use any of the other parties' logos outside of the sole purpose indicated in this Agreement without the expressed written consent of that party. Neither party may change or alter either parties' logo in any way.

31. INDEPENDENT CONTRACTOR

Seller agrees that the relationship between Purchaser and Seller is neither an employer-employee relationship nor an agency relationship. Rather, Seller is retained as and shall perform as an independent Seller to Purchaser. As an independent Seller, Seller will be solely responsible for determining the means and methods for performing the described services set forth on the face of this Agreement.

**AMENDMENT
TO
SIEMENS MEDICAL SOLUTIONS USA, INC. & DUKE UNIVERSITY AND HEALTH
SYSTEM
TERMS AND CONDITIONS OF SALE**

THIS AMENDMENT OF SIEMENS (this "Amendment") is entered into as of the 29 day of September, 2009 by and between SIEMENS MEDICAL SOLUTIONS USA, INC., ("SIEMENS"), and DUKE UNIVERSITY & HEALTH SYSTEM ("BUYER").

WITNESSETH:

BUYER and SIEMENS are parties to the Siemens Terms and Conditions of Sale dated as of 31 August, 2009, as amended. BUYER and SIEMENS wish to now amend the Siemens Terms and Conditions of Sale

NOW THEREFORE, in consideration of the promises set forth in the Terms and Conditions of Sale and for other good and valuable consideration, the receipt of which is hereby acknowledged, and intending to be legally bound, the parties hereto agree as follows:

1.1 In the event there is a conflict between the terms and conditions of this Amendment and Siemens Terms and Conditions of Sale, the terms and conditions of this Amendment shall control. The Siemens Terms and Conditions of Sale and this Amendment shall hereinafter be referred to collectively as the "Agreement."


1.2 Add a new Section 32, General Indemnification, as follows:


"Each party hereto (as the "Indemnifying Party") agrees to indemnify and hold harmless the other party (as the "Indemnified Party") and its directors, officers, members, employees and agents from and against any losses, judgments, claims, costs, expenses (including reasonable attorneys' fees), liabilities, or damages (collectively "Losses"), for bodily injuries, death, and/or tangible property damage, asserted against the Indemnified Party by a third party and resulting from the Indemnifying Party's breach of its obligations under this Agreement or the negligent act or omission of the Indemnifying Party or its directors, officers, employees or agents in connection with this Agreement. The indemnification obligations of this Section shall survive termination of this Agreement."

1.3 Except as expressly amended hereby, the Agreement (Siemens Terms and Conditions of Sale and this Amendment) shall remain in full force and effect and the parties hereby ratify and confirm their rights, duties and obligations.


IN WITNESS WHEREOF, the parties hereto have entered into this Amendment the day and year first above written.

SIEMENS MEDICAL SOLUTIONS USA, INC.:

By: 
Name: David B. Hutton
Title: Zone Controller

By: 
Name: Michael Wendt
Title: Zone GM&VP

DUKE UNIVERSITY & HEALTH SYSTEM

By: 
Name: Brooke Benson
Title: Director Clinical Resource Mgmt

SIEMENS

Siemens Medical Solutions USA, Inc.
51 Valley Stream Parkway, Malvern, PA 19355
Fax: (866) 309-6967

SIEMENS REPRESENTATIVE
Stuart Waddey - (919) 605-9227

Software License Schedule to the Siemens Medical Solutions USA, Inc. General Terms and Conditions

1. DEFINITIONS: The following definitions apply to this Schedule:

"Agreement" shall mean the attached (i) Quotation for Products and/or Services including the Terms and Conditions of Sale and applicable schedules; and/or (ii) Software License Agreement describing the software licensed herein and the specific system for which the license is issued.

"Licensor" shall mean Siemens Medical Solutions USA, Inc.

"Licensee" shall mean the end-user to whom Licensor provides Software or Documentation for its internal use under the Agreement.

"Software" shall mean the software described in the attached Agreement, including the following as contained therein: (i) software programs consisting of a series of statements or instructions to be used directly or indirectly in a programmable controller or computer to bring about a certain result and (ii) databases consisting of systemized collections of data to be used or referenced directly or indirectly by a programmed controller or computer. Notwithstanding the foregoing, "Software" does not include "firmware" as such term is conventionally understood. Diagnostic/Maintenance Software also is not included within the scope of the Software licensed under this Schedule, and is available only as a special option under a separate Diagnostic Materials License Agreement and may be subject to a separate licensing fee.

"Documentation" shall mean the documents and other supporting materials which are intended to support the use of an associated product, including (but not limited to) instructions, descriptions, flow charts, logic diagrams and listings of the Software, in text or graphic form, on machine readable or printed media.

"Designated Unit" shall mean a single control unit or computer identified on the first page of the Agreement, on which Software licensed hereunder may be used by Licensee.

2. SCOPE: The following terms and conditions shall apply to all Software and Documentation provided by Licensor to Licensee under the Agreement (whether included with other products listed in the Agreement or listed separately in the Agreement), together with any updates or revisions thereto which Licensor may provide to Licensee, and all copies thereof, except any Software and/or Documentation licensed directly by Licensor's supplier under a separate end-user license agreement accompanying the Software or the Documentation, in which case Licensee agrees to be bound by that license agreement as a condition to using the Software and/or Documentation. Except as expressly provided herein, and provided that in no event shall the warranties or other obligations of Licensor with respect to such Software or Documentation exceed those set forth in this Schedule, this Schedule shall be subject to the liability limitations and exclusions and other terms and conditions set forth in the Agreement. **ANY USE OF THE SOFTWARE, INCLUDING BUT NOT LIMITED TO USE ON THE DESIGNATED UNIT, WILL CONSTITUTE LICENSEE'S AGREEMENT TO THIS SOFTWARE LICENSE SCHEDULE (OR RATIFICATION OF ANY PREVIOUS CONSENT).**

3. SOFTWARE AND DOCUMENTATION LICENSE: Subject to the payment of any applicable annual license fee(s), whether stated separately or included in the purchase price of another product, and to Licensee's acceptance of all of the obligations set forth herein and to the fulfillment of those obligations, Licensor or, if applicable, its licensor or supplier, hereby grants to Licensee a paid-up, nonexclusive and nontransferable (except as expressly provided in this Schedule) limited license to use the Software provided by Licensor under the Agreement solely for Licensee's own use on the Designated Unit and to use the Documentation in support of Licensee's authorized use of the Software, for the purpose of operating the Designated Unit in accordance with the instructions set forth in the user's manual supplied with the Designated Unit and for no other purpose whatsoever. A separate license is required for each Designated Unit on which the Software is to be used. Licensee may obtain from Licensor one copy of the Software licensed hereunder for backup and archival purposes only as is necessary to support Licensee's own authorized use of the Software, provided that Licensee includes on or in all copies (in any form) all copyright, trade secret or other proprietary notices contained on or in the Software as provided by Licensor. Additional copies of the Documentation may be licensed from Licensor at its then applicable charges. Licensee may make the Software and Documentation (including any copies) available only to its employees and other persons on Licensee's premises to whom such disclosure is necessary to enable Licensee to use the Software or Documentation within the scope of the license provided in this Schedule. If the Software is

supplied to any unit or agency of the United States Government other than the Department of Defense, the Software and Documentation are classified as "restricted computer software" and the Government's rights in the Software and Documentation shall be as provided in paragraph (c) (2) of the Commercial Computer Software-Restricted Rights clause in FAR 52.227-19 and any successor laws, rules or regulations thereto. If the Software is supplied to the United States Department of Defense, the Software is classified as "commercial computer software" and the Government is furnished the Software and Documentation with "restricted rights" as defined in paragraph (c) (1) of the Rights in Technical Data and Computer Software clause in DFARS 252.227-7013 and any successor laws, rules or regulations thereto.

4. PROPRIETARY PROTECTION AND CONFIDENTIALITY: Ownership of and title to the Software and Documentation and all copies, in any form, licensed under this Schedule are and will remain in Licensor or its suppliers at all times. Licensee shall not (i) remove any copyright, trade secret or other proprietary right notices contained on or in the Software or Documentation as provided by Licensor, (ii) reproduce or modify any Software or Documentation or copy thereof, (iii) reverse assemble, reverse engineer or decompile any Software, or copy thereof, in whole or in part (except and only to the extent that such activity is expressly permitted by applicable law notwithstanding this limitation), (iv) sell, transfer or otherwise make available to others the Software or Documentation, or any copy thereof, except as expressly permitted by this Schedule, or (v) apply any techniques to derive any trade secrets embodied in the Software or Documentation. Licensee shall take all appropriate actions to ensure that: (i) the Software does not leave the Designated Unit's equipment location as set forth above, (ii) the Software is not copied by Licensee or any third parties, and (iii) the Software is not used in any equipment other than the Designated Unit. Licensee shall secure and protect the Software and Documentation and copies thereof from disclosure and shall take such actions with its employees and other persons who are permitted access to the Software or Documentation or copies as may be necessary to satisfy Licensee's obligations hereunder. Prior to disposing of any computer medium, computer memory or data storage apparatus, Licensee shall ensure that all copies of Software and Documentation have been erased therefrom or otherwise destroyed. In the event that Licensee becomes aware that any Software or Documentation or copies are being used in a manner not permitted by the license, Licensee shall immediately notify Licensor in writing of such fact and if the person or persons so using the Software or Documentation are employed or otherwise subject to Licensee's direction and control, Licensee shall use reasonable efforts to terminate such impermissible use. Licensee will fully cooperate with Licensor so as to enable Licensor to enforce its proprietary and property rights in the Software. Licensee agrees that, subject to Licensee's reasonable security procedures, Licensor shall have immediate access to the Software at all times and that Licensor may take immediate possession thereof upon termination or expiration of the associated license or this Schedule. Licensee's obligations under this paragraph shall survive any termination of a license, the Schedule or the Agreement.

5. UPDATES AND REVISIONS: During the warranty period or under a separate service contract or software update subscription, revised or updated versions of the Software licensed under this Schedule may be made available, at Licensor's option, to Licensee to use or to test while Licensee continues use of a previous version. Licensee has the right to decide whether to install any such revised or updated versions or to continue use of the previous version after giving due regard to the United States Food and Drug Administration rules and regulations. However, Licensee shall pay Licensor for any services necessitated by any modifications of the Software by Licensee or by Licensee's failure to utilize the current non-investigational version of the Software provided by Licensor. Software updates that provide new features or capabilities or that require hardware changes will be offered to Licensee at purchase prices established by Licensor. Licensor retains the sole right to determine whether an update represents an enhancement of a previously purchased capability or a new capability for which the Licensee will be charged. In addition, some updates may require Applications Training performed by Licensor's personnel that will be offered at Licensor's prevailing rates. Licensor retains the sole right to determine whether an update requires such training.

SIEMENS

Siemens Medical Solutions USA, Inc.
51 Valley Stream Parkway, Malvern, PA 19355
Fax: (866) 309-6967

SIEMENS REPRESENTATIVE
Stuart Waddey - (919) 605-9227

6. DELIVERY, RISK OF LOSS AND TITLE: Notwithstanding the provisions of Section 6 of the attached Terms and Conditions of Sale, if any, the Software and Documentation licensed hereunder shall be delivered on or about the delivery date stated in the Agreement unless a separate delivery date is agreed upon. If Software or Documentation licensed hereunder is lost or damaged during shipment from Licensor, Licensor will replace it at no charge to Licensee. If any Software or Documentation supplied by Licensor and licensed hereunder is lost or damaged while in the possession of Licensee, Licensor will replace it at Licensor's then current applicable charges, if any, for materials, processing and distribution. Notwithstanding the provisions of Section 6 of the attached Terms and Conditions of Sale, if any, the Software and Documentation, in any form, and all copies made by Licensee, including partial copies, and all computer media provided by Licensor are and remain the property of Licensor or its supplier. Licensee has no right, title or interest in the Software, the Documentation, or any computer media provided by Licensor, or copies, except as stated herein, and ownership of any such Software, Documentation and computer media shall at all times remain with Licensor or its suppliers.

7. LICENSE TRANSFER: The Software and Documentation, and the license hereunder, may not be assigned, transferred or sublicensed except as hereinafter provided. Upon the sale or lease of the Designated Unit to a third party, Licensee may transfer to such third party, with Licensor's written consent and in accordance with Licensor's then current policies and charges, the license to use the Software and Documentation hereunder, together with the Software, the Documentation, the computer media provided by Licensor, and all copies provided that: (i) Licensee notifies Licensor in writing of the name and address of such third party; (ii) such third party agrees in a written instrument delivered to Licensor to the terms of this Schedule; and (iii) Licensee does not retain any copies of the Software or Documentation in any form.

8. WARRANTIES: Licensor warrants that for the warranty period provided by Licensor under the attached Terms and Conditions of Sale, if any, the Software shall conform in all material respects to Licensor's published specifications as contained in the applicable supporting Documentation. This paragraph replaces Paragraphs 10.1 and 10.4 of any such Terms and Conditions of Sale with respect to the Software and Documentation. Such Documentation may be updated by Licensor from time to time and such updates may constitute a change in specification. Licensee acknowledges that the Software is of such complexity that it may have inherent or latent defects. As Licensee's sole remedy under the warranty, Licensor will provide services, during the warranty period, to correct documented Software errors which Licensor's analysis indicates are caused by a defect in the unmodified version of the Software as provided by Licensor. Licensor does not warrant that the Software will meet Licensee's requirements, or will operate in combinations which may be selected for use by Licensee, or that the operation of the Software will be uninterrupted or error free. Licensee is responsible for determining the appropriate use of and establishing the limitations of the Software and its associated Documentation as well as the results obtained by use thereof.

LICENSOR MAKES NO WARRANTY WITH RESPECT TO THE SOFTWARE AND DOCUMENTATION OTHER THAN THOSE SET FORTH IN THIS SECTION. THE WARRANTY HEREIN IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED, INCLUDING BUT NOT LIMITED TO ANY EXPRESS OR IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, WHICH ARE HEREBY DISCLAIMED, AND CONSTITUTES THE ONLY WARRANTY MADE WITH RESPECT TO THE SOFTWARE AND DOCUMENTATION.

9. LICENSE TERM AND TERMINATION: The license for the Software and Documentation is effective on the shipment date of the Software and Documentation (F.O.B. shipping point or F.A.S., as the case may be) and continues until Licensee's possession of the Software and all copies ceases (except in connection with a transfer of the license as permitted by this Schedule) or until otherwise terminated as provided herein. Licensee may terminate the license for the Software and Documentation at any time after discontinuance of use of the Software and Documentation and all copies, upon written notice to Licensor. If Licensee (i) fails to comply with its obligations herein and does not cure such failure within ten (10) days after receipt of notice from Licensor, or (ii) attempts to assign the Agreement or this Schedule or any rights or obligations hereunder without Licensor's prior written consent, then Licensor may terminate the license hereunder and require the immediate discontinuance of all use of the Software and Documentation and all copies thereof in any form, including modified versions and updated works. Within five (5) days after the termination of the license, Licensee shall, at Licensor's option either: (i) return to Licensor the Software and Documentation, and all copies, in any form, including updated versions, along with any computer media provided by Licensor; or (ii) destroy the affected Software and Documentation, and all copies, in any form, including updated versions, and certify such return or destruction in writing to Licensor.

Created: 2/14/2013 9:53:00 AM
PRO 1-5OJS40

Siemens Medical Solutions USA, Inc. Confidential

10. MISCELLANEOUS: Since the unauthorized use of the Software and/or Documentation may leave Licensor without an adequate remedy at law, Licensee agrees that injunctive or other equitable relief will be appropriate to restrain such use, threatened or actual. Licensee further agrees that to the extent applicable, (i) any of Licensor's suppliers of Software and/or Documentation is a direct and intended beneficiary of this Schedule and may enforce it directly against Licensee with respect to the Software and/or Documentation provided by such supplier, and that (ii) **NO SUPPLIER OF LICENSOR SHALL BE LIABLE FOR ANY GENERAL, SPECIAL, DIRECT, INDIRECT, CONSEQUENTIAL, INCIDENTAL OR OTHER DAMAGES ARISING OUT OF ANY SUBLICENSE OF THE SOFTWARE AND/OR DOCUMENTATION. THIS LIMITATION ON LIABILITY SHALL APPLY EVEN IF ANY REMEDY FAILS OF ITS ESSENTIAL PURPOSE.**

11. ADDITIONAL PROVISIONS RELATING TO THIRD-PARTY SOFTWARE: If the Software includes software licensed by Licensor from third parties, the following additional provisions shall apply:

(a) If Software is provided by Licensor on separate media and labeled "Recovery Media," Licensee may use the Recovery Media solely to restore or reinstall the Software and/or Documentation originally installed on the Designated Unit.

(b) Licensee is licensed to use the Software to provide only the limited functionality (specific tasks or processes) for which the Designated Unit has been designed and marketed by Licensor. This license specifically prohibits any other use of the software programs or functions, or inclusion of additional software programs or functions that do not directly support the limited functionality, on the Designated Unit. If Licensee uses the Designated Unit to access or utilize the services or functionality of Microsoft Windows Server products (such as Microsoft Windows NT Server 4.0 (all editions) or Microsoft Windows 2000 Server (all editions)), or uses the Designated Unit to permit workstation or computing devices to access or utilize the services or functionality of Microsoft Windows Server products, Licensee may be required to obtain a Client Access License for the Designated Unit and/or each such workstation or computing device. Licensee should refer to the end user license agreement for its Microsoft Windows Server product for additional information.

(c) The Software may contain support for programs written in Java. Java technology is not fault tolerant and is not designed, manufactured, or intended for use or resale as online control equipment in hazardous environments requiring fail-safe performance, such as in the operation of nuclear facilities, aircraft navigation or communication systems, air traffic control, direct life support machines, or weapons systems, in which the failure of Java technology could lead directly to death, personal injury, or severe physical or environmental damage. Sun Microsystems, Inc. has contractually obligated Licensor's supplier to make this disclaimer.

(d) The Software may permit Licensor, its supplier(s), or their respective affiliates to provide or make available to Licensee Software updates, supplements, add-on components, or Internet-based services components of the Software after the date Licensee obtains its initial copy of the Software ("Supplemental Components").

- If Licensor provides or makes available to Licensee Supplemental components and no other end-user software licensing agreement terms are provided along with the Supplemental Components, then the terms of this Software License Schedule shall apply.

- If a supplier of Licensor or affiliates of such a supplier make available Supplemental Components, and no other end-user software licensing agreement terms are provided, then the terms of this Schedule shall apply, except that the supplier or affiliate entity providing the Supplemental Component(s) shall be the licensor of the Supplemental Component(s). Licensor, its supplier(s), and their respective affiliates reserve the right to discontinue any Internet-based services provided to Licensee or made available to Licensee through the use of the Software.

(e) **The Software and Documentation supplied by Licensor's suppliers are provided by such suppliers "AS IS" and with all faults. SUCH SUPPLIERS DO NOT BEAR ANY OF THE RISK AS TO SATISFACTORY QUALITY, PERFORMANCE, ACCURACY, OR EFFORT (INCLUDING LACK OF NEGLIGENCE) WITH RESPECT TO SUCH SOFTWARE AND DOCUMENTATION. ALSO, THERE IS NO WARRANTY BY SUCH SUPPLIERS AGAINST INTERFERENCE WITH LICENSEE'S ENJOYMENT OF THE SOFTWARE OR AGAINST INFRINGEMENT. IF LICENSEE HAS RECEIVED ANY WARRANTIES REGARDING THE DESIGNATED UNIT OR THE SOFTWARE, THOSE WARRANTIES DO NOT ORIGINATE FROM, AND ARE NOT BINDING ON, LICENSOR'S SUPPLIERS.**

(f) Licensee acknowledges that portions of the Software are of U.S. origin. Licensee agrees to comply with all applicable international and national laws that apply to the Software, including the U.S. Export Administration Regulations, as well as applicable end-user, end-use and destination restrictions issued by U.S. and other governments. For additional

Page 16 of 42

SIEMENS

Siemens Medical Solutions USA, Inc.
51 Valley Stream Parkway, Malvern, PA 19355
Fax: (866) 309-6967

SIEMENS REPRESENTATIVE
Stuart Wadley - (919) 605-9227

Information on exporting software supplied by Microsoft, see
<http://www.microsoft.com/exporting/>.

Revised 03/15/05

SIEMENS

Siemens Medical Solutions USA, Inc.
51 Valley Stream Parkway, Malvern, PA 19355
Fax: (866) 309-6967

SIEMENS REPRESENTATIVE
Stuart Wadley - (919) 605-9227

MR Warranty Information

| <u>Product</u> (New Systems and "Proven Excellence" Refurbished Systems Only) | <u>Period of Warranty</u> ¹ | <u>Coverage</u> |
|---|--|----------------------------------|
| MR System (not including consumables) | 12 month | Full Warranty (parts & labor) |
| <u>Post Warranty (after expiration of system warranty) – Replacement parts only!</u> | | |
| Magnet | 12 month | Parts only |
| Spare Parts | 6 month | Parts only |
| Consumables | Not Covered | |

Note: Optional extended warranty coverage can be obtained by purchase of a service agreement.

¹ Period of warranty commences from the date of first use or completion of installation, whichever occurs first. In the event the completion of installation is delayed for reasons beyond Siemens' control, the stated warranty period shall commence 60 days after delivery of equipment.

SIEMENS

Siemens Medical Solutions USA, Inc.
51 Valley Stream Parkway, Malvern, PA 19355
Fax: (866) 309-6967

SIEMENS REPRESENTATIVE
Stuart Waddey - (919) 605-9227

Detailed Technical Specifications

MAGNETOM Aera - USA

| Part No. / Product | Description |
|--|---|
| <p>14416900 MAGNETOM Aera - System</p> | <p>Aera ex Erlangen ENS_14416900</p> <p>MAGNETOM Aera - the first 1.5T Tim+Dot system - integrates the next generation Tim - Tim 4G and the Siemens unique Dot Engines (Day Optimizing Throughput Engines) enabling workflow efficiency combined with higher diagnostic confidence due to consistent results.</p> <p>The system includes:</p> <p>Tim 4G+Dot</p> <p>Tim 4G provides increased patient comfort and optimized workflow efficiency. Only one patient setup, no repositioning, and no changing of coils. Ultra-light-weighted coils with high density of coil elements for maximized patient comfort and increased SNR. Feet-first positioning for almost all examinations possible reduces claustrophobia.</p> <p>Tim 4G is 4G flexibility, accuracy and speed and brings image quality and acquisition speed to a new level.</p> <p>Dot helps to take away the complexity in MRI scanning through patient personalization, user guidance and process automation. Optimized scan strategies are provided and can be selected based on patient condition, which allows for high quality exams even when conditions change. Integrated decision points allow the user to easily add or remove one or a group of protocols with one click. Step by step real-time on board guidance guides novice users even through the most complicated exams. Process automation allows optimal timing for breathing, scanning, and planning. Dot can be easily customized to follow the individual standards of care.</p> <p>Dot is personalized, guided and automated and designed to improve workflow efficiency and image consistency.</p> <p>MAGNETOM Aera with its 70 cm Open Bore design and a system length of only 145 cm gives a patient friendly appearance that can significantly help patients with anxiety or claustrophobia.</p> <p>Magnet:</p> <ul style="list-style-type: none"> - Ultra-short 137 cm long (145 cm with covers), whole-body superconductive 1.5T magnet with active shielding (AS) technology with counter coils - External Interference Shielding (E.I.S.) - Excellent homogeneity enabled by TrueForm magnet design which allows for a cylindrically optimized homogeneity volume resulting in higher image quality (50 x 50 x 45 cm³ DEV, typ. 3.6 ppm based on the 24-plane plot method) - The magnet has a helium capacity of approximately 1,280 liters and a typical Helium boil-off rate of 0 l/yr during typical, undisturbed clinical operation depending on the sequences used and examination time, and provided the system is serviced in regular intervals. - It has an integrated magnet cooling system. <p>Gradient system:</p> <ul style="list-style-type: none"> - Actively shielded water-cooled world-class gradient system - True Form Gradient Design - All axes force compensated <p>DirectRF - RF Transmit/Receive System:</p> <ul style="list-style-type: none"> - Fully integrated Transmit and Receive path in the magnet housing including extremely compact water-cooled solid state amplifier with 26.1 kW peak power - High dynamic range - Immediate feedback loop for real-time sequence adaptation - Integrated no tune transmit/receive Body Coil - The revolutionary Tim 4G technology allows connecting up to 204 coil elements simultaneously enabling higher SNR and iPAT in all directions. No repositioning of patients is needed even for large Field of View |

SIEMENS

Siemens Medical Solutions USA, Inc.
 51 Valley Stream Parkway, Malvern, PA 19355
 Fax: (866) 309-6967

SIEMENS REPRESENTATIVE
 Stuart Waddey - (919) 605-9227

| Part No. / Product | Description |
|---|--|
| <p>(Continued) 14416900 MAGNETOM Aera - System</p> | <p>examinations.</p> <ul style="list-style-type: none"> - Dual-Density Signal Transfer enables ultra-high density coil design by integrating key RF components into the local coil. <p>Tim 4G Coils: The new Tim 4G coil technology with Dual-Density Signal Transfer, DirectConnect and SlideConnect Technology combines key imaging benefits: Excellent image quality, high patient comfort, and unmatched flexibility.</p> <p>The Tim 4G coils are designed for highest image quality combined with easy handling. The high element density of the coils increases SNR and reduces examination times. DirectConnect and SlideConnect™ technology reduce patient set up time significantly. The coils are designed with the patient in mind. Light weighted coils and open design ensure highest patient comfort which results in better patient cooperation and image quality. No coil changing with multi-exam studies saves patient setup- and table time. AutoCoilSelect enables dynamic, automatic, or interactive selection of the coil elements within the Field of View and speeding the exam preparation at the host. All coils are time-saving "no-tune" coils. A comprehensive set of pads for comfortable and stable patient positioning together with safety straps are included.</p> <ul style="list-style-type: none"> - Head/Neck 20 The 20-channel coil with its 20 integrated pre-amplifiers ensures excellent signal-to-noise ratio. The unique DirectConnect technology allows users connecting the 20 coil elements of the Head/Neck20 without cables. The patient friendly open design allows for maximum patient comfort which is supported in addition by a look-out mirror for claustrophobic patients. The high channel coil is iPAT compatible in all directions. <p>The open and light design of the upper coil part increases patient comfort and is removable for easy patient handling. The lower coil part may remain on the table for most of the examinations can be used without the upper part. The Head/Neck 20 and Spine 32 are smoothly integrated into the patient table, thus enabling high flexibility in imaging and fewer coil changes and easy handling when switching patients. The Head/Neck 20 coil is equipped with two removable cushioned head stabilizers for stable and comfortable patient positioning.</p> <p>The Head/ Neck 20 can be used for applications like head examinations, neck examinations, MR Angiography, combined head/neck examinations or for imaging of the TMJ (temporomandibular joints).</p> <p>Typically combined with the Spine 32 and Body 18 or Peripheral Angio 36 but also other combinations eg with flexible coils like the Flex Large 4 are possible.</p> <ul style="list-style-type: none"> - Body 18 The 18-channel coil with its 18 integrated pre-amplifiers ensures maximum signal-to-noise ratio. The 18 coil elements of the Body 18 with only one SlideConnect Plug allows for fast and easy patient preparation resulting in less table time. Fast acquisition times enabled by iPAT in all directions. The light-weighted coil ensures highest patient comfort. <p>Body 18 operates in an integrated fashion with the Spine 32 as a 30 channel body coil</p> <p>Body 18 can be combined with further Body 18 coils for larger coverage and positioned in different orientations (0°, 90°, 180°, 270°) for patient specific adaptations</p> <p>The Body 18 is typically used in combination with the Spine 32 for examinations of the thorax, abdomen, pelvis or hip and operates as a 30 channel body coil (3 rings 10 elements). The Body 18 can also be used for cardiac or vascular applications. Through its perfect combinability with the Spine 32, further Body 18 (optional), the Peripheral Angio 36 (optional), but also the Head/Neck20 and all flexible coils (e.g. Flex Large 4, Flex Small 4) it contributes for a broad range of indications up to whole-body imaging.</p> <ul style="list-style-type: none"> - Spine 32 The 32-channel coil with its 32 integrated pre-amplifiers ensures maximum signal-to-noise ratio. The unique DirectConnect technology allows connecting the 32 coil elements of the Spine 32 without the need to plug in any cable. The patient friendly ergonomic design allows for maximum patient comfort. The high element coil is iPAT compatible in all directions. <p>Smoothly integrated into the patient table the Spine 32 may remain on the patient table for nearly all exams.</p> |

SIEMENS

Siemens Medical Solutions USA, Inc.
 51 Valley Stream Parkway, Malvern, PA 19355
 Fax: (866) 309-6967

SIEMENS REPRESENTATIVE
 Stuart Waddey - (919) 605-9227

| Part No. / Product | Description |
|--------------------|--|
| | The Spine 32 is typically combined with Body 18, Head/Neck 20, Peripheral Angio 36 or Flex Large 4, Flex |

| Part No. / Product | Description |
|---|---|
| <p>(Continued) 14416900 MAGNETOM Aera - System</p> | <p>Small 4.</p> <ul style="list-style-type: none"> - Flex Large 4/ Flex Small 4 Light-weighted, very flexible, iPAT compatible, 4-element no-tune receiver coils which are made of soft and smooth material. The coils can be wrapped around or used flat. <p>Both coils can be connected via Flex Coil interface. One Flex Coil interface is already delivered as standard.</p> <p>The coils can be used for different examinations ranging from examinations of the extremities to abdominal examinations.</p> <p>TimTable</p> <ul style="list-style-type: none"> - The maximum scan range of the Tim Table is 140 cm. A scan range of 205 cm can be achieved with the Tim Whole Body suite (optional) - The maximum patient weight of 250 kg (550 lbs) is valid for horizontal and vertical movements, which ensures maximized patient comfort for obese patients. - The patient table can be lowered to a minimum height of 52 cm from the floor, for easier patient positioning and better accessibility for geriatric, pediatric or immobile patients. An infusion stand is integrated to ensure fast patient set up also for critical patients. - Multiple Tim4G coils can be connected at once for efficient and patient friendly examinations. - The Tim Table can be moved with two clicks into the isocenter - one click to the upmost position and one click into the isocenter. <p>Dot (Day Optimizing Throughput) Engine Dot multiplies the power of Tim resulting in greater image consistency and diagnostic confidence</p> <p>Dot Control Centers and Dot Display</p> <ul style="list-style-type: none"> - The ergonomically designed Dot Control Centers are integrated left and right into the front covers for controlling table movement and interaction with the Dot Display. The Dot Control Centers are well illuminated for easy visual recognition. - Automated table move up to upmost position, to center position or Home position facilitate smooth patient preparation and will reduce table time - Variable (6 levels) ventilation and lighting inside the magnet bore or volume adjustments are possible for increased patient comfort - The Dot Display provides on board guidance for patient set up where it's needed - directly at the scanner. Information such as Patient name or exam type or required patient position, guidance for ECG set up and immediate visualization of physiological curves will be provided for convenient operation. - Almost all table control functions, including ventilation and illumination of the magnet bore, can be also controlled from the operator console for convenient operation. <p>Dot Technology Dot makes it easy to get the best possible results for virtually any type of patient. Dot gives uniquely tailored, optimized scans configurable to patient condition or clinical question. Dot provides patient personalization, user guidance and process automation and is of course configurable by the user to adapt to the different clinical needs and standards of care.</p> <p>Brain Dot Engine The Brain Dot Engine simplifies general brain examinations with guided and automated workflows customized to the site specific standards of care. The Brain Dot Engine supports the user in achieving reproducible image quality with increased ease of use and time efficient exams. The brain workflow can be personalized to the individual patient condition and clinical need. Several predefined strategies are included, which can be easily selected with one click. They can be changed at any time during the brain workflow.</p> |

SIEMENS

Siemens Medical Solutions USA, Inc.
 51 Valley Stream Parkway, Malvern, PA 19355
 Fax: (866) 309-6967

SIEMENS REPRESENTATIVE
 Stuart Waddey - (919) 605-9227

| Part No. / Product | Description |
|--------------------|--|
| | <p>Protocols tailored for use of contrast media are integrated.</p> <ul style="list-style-type: none"> - Standard: Standard examination with 2D protocols - Resolution focus: Examination with 3D protocols (with e.g. SPACE) for detailed views - Speed focus: Examination with fast 2D protocols (with e.g. HASTE) for further speeding up the exam - Limited patient capabilities: Examination with <i>syngo</i> BLADE protocols - to minimize and correct for the effects of motion automatically |

| Part No. / Product | Description |
|---|---|
| <p><i>(Continued)</i> 14416900 MAGNETOM Aera - System</p> | <p>Step-by-step user guidance is seamlessly integrated. Example images and guidance text are displayed for each individual step of the scanning workflow. Both - images and text - are easily configurable by the user.</p> <p>Easy positioning of the patient with AutoPosition. The patient is automatically placed at the isocenter without any laser marking required.</p> <p>AutoAlign Head allows automatically slice positioning and aligns on the anatomically derived sagittal, coronal, and axial slices of the localizer. The operator-free alignment and anatomical marking are consistent, independently of patient age, head position, or disease.</p> <p>Automatic real-time calculation of trace-weighted images and ADC maps with Inline DiffusionTechnology.</p> <p>Easy rerun or repeat with functionality allows for reduced table time even in case of patients with pain or claustrophobia. An image inside the examination UI can be selected and a rerun of the corresponding series can be triggered with identical sequences or parameters. Alternatively an exam can be repeated with a changed strategy.</p> <p>The Brain Dot Engine as all Dot Engines can be modified by the user to their individual standard of care.</p> <p>Tim Application Suite The Tim Application Suite offers a complete range of clinically optimized sequences, protocols and workflow functionalities for all body regions. Excellent head-to-toe imaging can be accomplished with the sequences and features included in this application suite. To enable this comprehensive application range, nine dedicated application packages have been included.</p> <ul style="list-style-type: none"> - Neuro Suite - Angio Suite - Cardiac Suite - Body Suite - Onco Suite - Breast Suite - Ortho Suite - Pediatric Suite - Scientific Suite <p>Neuro Suite Comprehensive head and spine examinations can be performed with dedicated programs. High resolution protocols and fast protocols for uncooperative patients are provided. The Neuro Suite also includes protocols for diffusion imaging, perfusion imaging, and fMRI. It includes for example:</p> <ul style="list-style-type: none"> - EPI sequences and protocols for diffusion, perfusion and fMRI for advanced neurological applications. Diffusion weighted imaging is possible with up to 16 b-values in the orthogonal directions. Dynamic Analysis software (included in standard configuration) enables calculation of: <ul style="list-style-type: none"> - ADC maps - t-test maps from the EPI images for fMRI - Time-to-Peak maps for perfusion analysis. - Whole spine protocols acquire in multiple steps via software controlled table movement in a single click. |

SIEMENS

Siemens Medical Solutions USA, Inc.
51 Valley Stream Parkway, Malvern, PA 19355
Fax: (866) 309-6967

SIEMENS REPRESENTATIVE
Stuart Wadley - (919) 605-9227

| Part No. / Product | Description |
|--------------------|--|
| | <ul style="list-style-type: none"> - 3D isotropic resolution volume imaging using T1 3D MPRAGE / 3D FLASH, SPACE DarkFluid, T2 SPACE and 3D TSE - T2-weighted high resolution 3D Restore protocols optimized for inner ear examinations - Whole-spine protocols in multiple steps with software controlled table movement - 2D and 3D MEDIC protocols for T2-weighted imaging, particularly for C-spine examinations in axial orientation where reproducibility is difficult due to CSF pulsations and blood flow artifacts - 3D Myelograms with 3D HASTE and 3D True-FISP for anatomical details - Dynamic sacro-iliac joint imaging after contrast administration using a fast T1-weighted FLASH 2D sequence - Spine diffusion protocols to differentiate osteoporosis versus tumor infiltration and post-radiotherapy changes versus residual tumor with PSIF sequence |

| Part No. / Product | Description |
|---|--|
| <p><i>(Continued)</i> 14416900 MAGNETOM Aera - System</p> | <ul style="list-style-type: none"> - Precision filter for high spatial accuracy e.g. for neuro intra-operative imaging and stereotactic planning - 3D CISS (Constructive Interference in Steady State) for excellent visualization of fine structures such as cranial nerves. High resolution imaging of inner ear and spine - AutoAlign Head LS providing a fast, easy, standardized, and reproducible patient scanning supporting reading by delivering a higher and more standardized image quality <p>Angio Suite Excellent MR Angiography can be performed to visualize arteries and veins with or without contrast agent. <i>Contrast-enhanced MRA</i></p> <ul style="list-style-type: none"> - 3D contrast-enhanced MRA protocols for e.g. single step, dynamic, peripheral, whole body MRA with the shortest TR and TE. The strong gradients make it possible to separate the arterial phase from the venous phase. - TestBolus workflow for optimized bolus timing and superb image quality. - CareBolus functionality for accurate determination of the bolus arrival time and the "Stop and Continue" of the 3D ce-MRA protocol after the 2D bolus control scan. - Dynamic ce-MRA for 3D imaging over time. <p><i>Non-contrast-MRA and venography</i></p> <ul style="list-style-type: none"> - 2D and 3D Time-of-Flight (ToF) protocols for MRA for the Circle of Willis, carotids, neck vessels, and breath-hold protocols for abdominal vessels - Triggered 2D ToF sequences for non-contrast MRA, particularly of the abdomen and the extremities - 2D/3D Phase-Contrast - MR venography with 2D/3D Time-of-Flight (ToF) and Phase-Contrast - TONE (Tilted Optimized Non-saturation Excitation) and MTC (Magnetization Transfer Contrast) techniques for improved Contrast-to-Noise Ratio (CNR) <p><i>Image processing tools</i></p> <ul style="list-style-type: none"> - MPR, MIP, MinIP, and 3D SSD - Inline MIP for immediate results - Inline subtraction of pre- and post-contrast measurements - Inline standard deviation maps of Phase-Contrast measurements for delineation of arteries and veins <p>Cardiac Suite The cardiac suite covers comprehensive 2D routine cardiac applications, ranging from morphology and ventricular function to tissue characterization. Featuring <i>syngo</i> BEAT 2D in conjunction with iPAT and T-PAT techniques. <i>Cardiac views</i></p> <ul style="list-style-type: none"> - Fast acquisition of the basic cardiac orientations for further examination planning - Cardiac scouting provides users with a step-by-step procedure for the visualization and planning of typical cardiac views, e.g. based on TrueFISP or Dark Blood TurboFLASH: short axis, 4-chamber and 2-chamber views. <p><i>syngo BEAT</i></p> <ul style="list-style-type: none"> - Unique tool for fast and easy cardiovascular MR imaging |

SIEMENS

Siemens Medical Solutions USA, Inc.
 51 Valley Stream Parkway, Malvern, PA 19355
 Fax: (866) 309-6967

SIEMENS REPRESENTATIVE
 Stuart Waddey - (919) 605-9227

| Part No. / Product | Description |
|--------------------|---|
| | <ul style="list-style-type: none"> - E.g. 1 click change from FLASH to TrueFISP for easy contrast optimization - 1-click to switch arrhythmia rejection on / off - 1-click change from Cartesian to radial sampling to increase effective image resolution (e.g. in pediatric patients) and avoid folding artifacts in large patients - <i>Visualization of structural cardiovascular pathologies with CMR - syngo BEAT</i> - Breath-hold and free breathing techniques for strong contrast between the blood and vascular structures. Dark Blood TSE and HASTE imaging are available for the structural evaluation of the cardiothoracic anatomy, including vessels or heart valves. Cine techniques (FLASH & TrueFISP) for high-resolution valve evaluation - Multiple contrasts such as T1- and T2-weighted imaging for use in diseases such as myocarditis (inflammation / hyperaemia), ARVD (fibrous-fatty degeneration) or acute myocardial infarction (edema) - Dark-blood TSE with motion compensation for high-quality vessel wall imaging in small or large vessels - <i>Tools for rapid evaluation of left or right ventricular function</i> - Acquisition of a stack of short-axis slices (standard segmented FLASH, or advanced segmented TrueFISP) - Automatic adjustment of the acquisition window to the current heart rate - Use of the Inline ECG for graphical ECG triggering setup - Retrospective gating with cine sequences (TrueFISP, FLASH) |

| Part No. / Product | Description |
|---|---|
| <p><i>(Continued)</i> 14416900 MAGNETOM Aera - System</p> | <ul style="list-style-type: none"> - Protocols for whole-heart coverage - iPAT integration for highest temporal and spatial resolution - Real-time imaging in case the patient is not able to hold his breath - <i>Dynamic imaging and tissue characterization with syngo BEAT</i> - Protocols for high-contrast and high-resolution tissue characterization - Protocols for stress and rest imaging with TrueFISP or TurboFLASH contrast support the acquisition of multiple slices with high resolution and arbitrarily adjustable slice orientation for each slice - T-PAT with mSENSE and GRAPPA for advanced parallel imaging provides fast high-resolution dynamic imaging - Segmented IR TrueFISP / FLASH with T1 scout for optimization of tissue contrast - Advanced tissue characterization with 2D phase-sensitive IR (PSIR) sequences TrueFISP and FLASH contrast. Magnitude and phase-sensitive images with one acquisition - Simple: no adjustment of inversion time (TI) necessary with PSIR technique - Ungated single-shot PSIR imaging for tissue characterization under difficult conditions: free-breathing technique that can be applied even in case of arrhythmia <p>Physiological Measurement Unit (PMU) - Wireless Physio Control</p> <ul style="list-style-type: none"> - Synchronizes the measurement with the physiological cycles (triggering to minimize motion artifacts caused by cardiac and respiratory movements) - Wireless Sensors - Wireless Vector ECG / respiration and pulse sensors for physiologically synchronized imaging, rechargeable battery-powered - for optimized patient handling - Physiological Signals Display - ECG (3 channels) - Pulse - Respiration - External Trigger Input Display <p>ECG Triggering:</p> <ul style="list-style-type: none"> - Acquisition of multiple slices, e.g. of the heart, at different phases of the cardiac cycle - Excellent image quality by synchronizing data acquisition with cardiac motion - Peripheral Pulse Triggering: - Reduces flow artifacts caused by pulsatile blood flow |

SIEMENS

Siemens Medical Solutions USA, Inc.
51 Valley Stream Parkway, Malvern, PA 19355
Fax: (866) 309-6967

SIEMENS REPRESENTATIVE
Stuart Waddey - (919) 605-9227

| Part No. / Product | Description |
|--------------------|--|
| | <ul style="list-style-type: none"> - Excellent image quality by synchronizing data acquisition to the pulsatile blood flow - Respiratory Triggering: - Excellent image quality by synchronizing data acquisition with the respiratory motion - External Triggering: - Interface for trigger input from external sources (e.g. Patient Monitoring System) inside the examination room - Interface for trigger input from external sources (e.g. pulse generator, trigger sources for fMRI) outside the examination room - Optical trigger output for fMRI - Retrospective gating for ECG, peripheral pulse, and external trigger input <p>Breast Suite MR imaging has proven a very high sensitivity for breast lesions and is the gold standard for the examination of silicone implants. Extremely high spatial and temporal resolution can be achieved in very short measuring times by using IPAT with GRAPPA. Excellent soft tissue differentiation, customized protocols (e.g. with fat saturation or water excitation or silicone excitation), as well as flexible multiplanar visualization allow for fast, simple and reproducible evaluation of MR breast examinations. This package includes:</p> <ul style="list-style-type: none"> - Quantitative evaluation and fast analysis of the data with colorized Wash-in, Wash-out, Time-To-Peak, Positive-Enhancement-Integral, MIPtime and combination maps with Inline technology or for offline calculation |

| Part No. / Product | Description |
|--|--|
| <p><i>(Continued)</i> 14416900 MAGNETOM Aera - System</p> | <ul style="list-style-type: none"> - High-resolution 2D protocols for morphology evaluation - High-resolution 3D protocols covering both breasts simultaneously - Protocols to support interventions (fine needle and vacuum biopsies, wire localization) - Protocols for evaluating breasts with silicone implants - Automatic and manual frequency adjustment, taking into account the silicone signal - Detection of the silicone signal either to suppress the silicone signal, if the surrounding tissue is to be evaluated, or to suppress the tissue signal in order to detect an implant leakage - SPAIR - robust fat sat (robust fat suppression using an adiabatic frequency selective inversion pulse) - DIXON - 2-point Dixon with 3D VIBE, the following contrasts can be obtained: in-phase, opposed phase, fat and water image. - IPAT with GRAPPA for maximum resolution in short time - Inline subtraction and MIP display - Offline subtraction, MPR and MIP display - <i>syngo</i> REVEAL: diffusion imaging for breast exams - IPAT Extension allows bilateral 3D sagittal breast imaging with Fat Sat or Water excitation <p>The Breast Suite also includes: <i>syngo</i> VIEWS (Volume Imaging with Enhanced Water Signal)</p> <ul style="list-style-type: none"> - bilateral - both breasts are examined simultaneously - axial - the milk ducts are directly displayed - fat-saturated or water-excited - fat complicates clinical evaluation and is suppressed - near-isotropic 3D measurement - the same voxel size in all three directions for reconstruction in any slice direction - submillimeter voxel - highest resolution for precise evaluation <p>Body Suite Body Suite covers your needs for clinical body applications. Ultrafast high resolution 2D and 3D protocols are provided for abdomen, pelvis, MR Colonography, MRCP, dynamic kidney, and MR Urography applications. Siemens unique 2D PACE technique makes body imaging easy allowing for multi-breath hold examinations as well as free breathing during the scans. Motion artifacts are greatly reduced with 2D PACE Inline technology.</p> |

SIEMENS

Siemens Medical Solutions USA, Inc.
 51 Valley Stream Parkway, Malvern, PA 19355
 Fax: (866) 309-6967

SIEMENS REPRESENTATIVE
 Stuart Waddey - (919) 605-9227

| Part No. / Product | Description |
|--------------------|--|
| | <p>This package includes:</p> <ul style="list-style-type: none"> - Free breathing 2D PACE applications with 2D/3D HASTE (RESTORE) and 2D/3D TSE (RESTORE) - Optimized fast single shot HASTE protocols and high-resolution 3D RESTORE protocols based on SPACE and TSE for MRCP and MR Urography examinations <p>ABDOMEN:</p> <p>2D:</p> <ul style="list-style-type: none"> - T1w (FLASH) breath-hold scans +/- Fat Sat (SPAIR, Q-FatSat, in-/opp-phase) - T2w (HASTE, TSE/BLADE, EPI) breath-hold scans +/- Fat Sat (SPAIR, FatSat, STIR) - T1w (TFL) triggered scans (2D PACE free breathing) in-/opp-phase - T2w (HASTE, TSE/BLADE, EPI) triggered scans (2D PACE free breathing) +/- Fat Sat (SPAIR, FatSat, STIR) as well as HASTE- and TSE-multi-echo - Optimized fast single shot HASTE protocols and high-resolution 3D RESTORE protocols based on SPACE and TSE for MRCP and MR urography examinations <p>3D:</p> <ul style="list-style-type: none"> - Dixon (VIBE 2pt-Dixon) breath-hold scans, following contrasts can be obtained: in-phase, opposed phase, fat and water image. - Dynamic (VIBE + Q-FatSat) protocols for best visualization of focal lesions with high spatial and temporal resolution - Colonography bright lumen with T2-weighted TrueFISP and dark lumen with T1-weighted VIBE <p>PELVIS:</p> <ul style="list-style-type: none"> - High-resolution T1w, T2w pelvic imaging (prostate, cervix) - Isotropic T2w SPACE 3D protocols for tumor search in the pelvis - Dynamic volume examinations with 3D VIBE - <i>syngo</i> REVEAL: diffusion imaging for liver and whole body exams |

| Part No. / Product | Description |
|--|---|
| <p><i>(Continued)</i> 14416900 MAGNETOM Aera - System</p> | <p>Onco Suite MR imaging has an excellent advantage of soft tissue contrast, multi-planar capabilities and the possibility of selectively suppressing specific tissue e.g. fat or water. This helps visualize pathologies, particularly metastases. The Onco Suite features a collection of sequences as well as protocols and evaluation tools that guide through a detailed screening of clinical indications, such as in hepatic neoplasms.</p> <p>This package includes:</p> <ul style="list-style-type: none"> - STIR TSE and HASTE, FLASH in-phase and opposed-phase protocols with a high sensitivity to metastases visualization - Dynamic imaging protocols for assessment of the kinetic behavior for lesion visualization and characterization - Quantitative evaluation and fast analysis of the data with colorized Wash-in, Wash-out, Time-To-Peak, Positive-Enhancement-Integral, MIPTIME and combination maps with Inline technology or for offline calculation - Display and analysis of the temporal behavior in selected regions of interest with the included MeanCurve post-processing application. This includes the capability of using additional datasets as a guide for defining regions of interest even faster and easier than before. - <i>syngo</i> REVEAL: diffusion imaging for liver and whole body exams <p>Dedicated prostate protocols for detection, localization, and staging of tumors and recurrences</p> <ul style="list-style-type: none"> - <i>syngo</i> REVEAL (diffusion-weighted imaging) - Protocols with high temporal resolution allow time course evaluation based on pharmacokinetic modeling <p>OrthoSuite Ortho Suite is a comprehensive collection of protocols for joint and spine imaging. MR imaging is especially suitable for avascular necrosis and internal derangements. The protocols included in this Suite can also be applied for imaging of tumors and infections.</p> <p>This package includes:</p> <ul style="list-style-type: none"> - 2D TSE protocols for PD, T1 and T2-weighted contrast with high in-plane resolution and thin slices |

SIEMENS

Siemens Medical Solutions USA, Inc.
51 Valley Stream Parkway, Malvern, PA 19355
Fax: (866) 309-6967

SIEMENS REPRESENTATIVE
Stuart Waddey - (919) 605-9227

| Part No. / Product | Description |
|--------------------|--|
| | <ul style="list-style-type: none"> - 3D MEDIC, 3D TrueFISP protocols with water excitation for T2-weighted imaging with high in-plane resolution and thin slices - High resolution 3D VIBE protocol for MR arthrography (knee, shoulder and hip) - 3D MEDIC, 3D TrueFISP, 3D VIBE protocols with water excitation having high isotropic resolution, optimized for 3D post-processing - PD SPACE with fat saturation and T2 SPACE with high isotropic resolution optimized for 3D post-processing - Whole spine single-step or multi-step protocols - Excellent fat suppression in off-center positions, e.g. in the shoulder due to high magnet homogeneity - Dynamic TMJ and ilio-sacral joint protocol - Susceptibility-insensitive protocols for imaging in the presence of a prosthesis - Multi-Echo SE sequence with up to 32 echoes for the calculation of T2 time maps (calculation included in the Scientific Suite) - High resolution 3D DESS (Double Echo Steady State): T2 / T1-weighted imaging for excellent fluid-cartilage differentiation <p>Pediatric Suite The parameters for pediatric imaging vary significantly in comparison to the parameters for adults. The reasons are developing tissues, body size, faster heart rates and restricted compliance with breath-hold commands. Protocols can be adapted for imaging infants.</p> <p>Scientific Suite Scientific Suite supports the scientifically oriented user with an easy access to application-specific data for further processing and advanced image computation methods.</p> <ul style="list-style-type: none"> - Support of USB memory sticks - Access to the file system by means of a secure and convenient browser - Anonymization of patient data - Easy generation of AVIs and screenshots for integration into presentations and training videos - Export function for tables, statistics and signal-time-courses in a communal format (MeanCurve, Spectroscopy, DTI evaluation) |

| Part No. / Product | Description |
|--|---|
| <p><i>(Continued)</i> 14416900 MAGNETOM Aera - System</p> | <ul style="list-style-type: none"> - Advanced image computation methods such as T2 and T1 time calculation, addition, subtraction, multiplication, division, and integration of images <p>The sequences, features and techniques for acquisition and reconstruction included in the Tim Application Suite are described in detail below.</p> <p>Sequences Spin Echo family of sequences:</p> <ul style="list-style-type: none"> - Spin Echo (SE) - Single, Double, and Multi Echo (up to 32 echoes); Inversion Recovery (IR) - 2D / 3D Turbo Spin Echo (TSE) - Restore technique for shorter TR times while maintaining excellent T2 contrast; TurboIR: Inversion Recovery for STIR, DarkFluid T1 and T2, TrueIR; Echo Sharing for dual-contrast TSE - 2D / 3D HASTE (Half-Fourier Acquisition with Single Shot Turbo Spin Echo) - Inversion Recovery for STIR and DarkFluid contrast - SPACE for 3D imaging with high isotropic resolution with T1, T2, PD, and DarkFluid Contrast <p>Gradient Echo family of sequences:</p> <ul style="list-style-type: none"> - 2D / 3D FLASH (spoiled GRE) - dual echo for in- / opposed phase imaging 3D VIBE (Volume Interpolated Breathhold Examination) - quick fat saturation; double echo for in-phase / opposed phase 3D imaging; DynaVIBE: Inline 3D elastic motion correction for multi phase data sets of the abdomen; Inline Breast Evaluation - 2D / 3D MEDIC (Multi Echo Data Image Combination) for high resolution T2 weighted orthopedic imaging and excellent contrast |

SIEMENS

Siemens Medical Solutions USA, Inc.
51 Valley Stream Parkway, Malvern, PA 19355
Fax: (866) 309-6967

SIEMENS REPRESENTATIVE
Stuart Waddey - (919) 605-9227

| Part No. / Product | Description |
|--------------------|--|
| | <ul style="list-style-type: none"> - 2D / 3D TurboFLASH - 3D MPRAGE; single shot T1 weighted imaging e.g. for abdominal imaging during free breathing - 3D GRE for field mapping - 2D / 3D FISP (Fast Imaging with Steady State Precession) - 2D / 3D PSIF - PSIF Diffusion - Echo Planar Imaging (EPI) - diffusion-weighted; single shot SE and FID e.g. for BOLD imaging and Perfusion-weighted imaging; 2D / 3D Segmented EPI (SE and FID) - ce-MRA sequence with Inline subtraction and Inline MIP - 2D / 3D Time-of-Flight (ToF) Angiography - single slab and multi slab; triggered and segmented - 2D / 3D Phase Contrast Angiography - <i>syngo</i> BEAT Tool - TrueFISP segmented; 2D FLASH segmented; - Magnetization-prepared TrueFISP (IR, SR, FS); IR TI scout; Retrogating <p>Standard Fat/Water Imaging</p> <ul style="list-style-type: none"> - Fat and Water Saturation. Additional frequency selective RF pulses used to suppress bright signal from fatty tissue. Two selectable modes: weak, strong - Quick FatSat - SPAIR: robust fat suppression for body imaging using a frequency selective inversion pulse - Fat / Water Excitation. Spectral selective RF pulses for exclusive fat / water excitation - Dixon technique for fat and water separation - based on VIBE (2 point Dixon) <p>Standard Techniques</p> <ul style="list-style-type: none"> - True Inversion Recovery to obtain strong T1-weighted contrast - Dark Blood inversion recovery technique that nulls fluid blood signal - Saturation Recovery for 2D TurboFLASH, gradient echo, and T1-weighted 3D TurboFLASH with short scan time (e.g. MPRAGE) - Freely adjustable receiver bandwidth, permitting studies with increased signal-to-noise ratio - Freely adjustable flip angle. Optimized RF pulses for image contrast enhancement and increased signal-to-noise ratio - MTC (Magnetization Transfer Contrast). Off-resonance RF pulses to suppress signal from certain tissues, thus enhancing the contrast. Used e.g. in MRA - Argus viewer for reviewing cine studies• - Report Viewer for DICOM structured reports including report editing |

| Part No. / Product | Description |
|---|---|
| <p><i>(Continued)</i> 14416900 MAGNETOM Aera - System</p> | <ul style="list-style-type: none"> - Dynamic Analysis for addition, subtraction, division, standard deviation, calculations of ADC maps, T1 and T2 values, TTP, t-Test, etc. - Image Filter - 3D post-processing MPR, MIP, MinIP, SSD - Flexible film formats and paper print - Data storage of images and cine AVI files on CD / DVD with DICOM viewer as the viewing tool for hand out to the patients or referrals - Selectable centric elliptical phase reordering via the user interface - Multiple Direction Diffusion Weighting (MDDW) - perform diffusion tensor imaging with multiple diffusion weightings and up to 12 directions for generating data sets. - Inversion Recovery to nullify the signal of fat, fluid or any other tissue <p>Standard techniques for Flow Artifact reductions</p> <ul style="list-style-type: none"> - LOTA (LongTerm Data Averaging) technique to reduce motion and flow artifacts - Pre-saturation techniques using RF saturation pulses to suppress flow and motion artifacts |

SIEMENS

Siemens Medical Solutions USA, Inc.
 51 Valley Stream Parkway, Malvern, PA 19355
 Fax: (866) 309-6967

SIEMENS REPRESENTATIVE
 Stuart Waddey - (919) 605-9227

| Part No. / Product | Description |
|--------------------|---|
| | <ul style="list-style-type: none"> - Tracking SAT bands maintain constant saturation of venous and/or arterial blood flow e.g. for 2D/3D sequential MRA - TONE (Tilted Optimized Non-saturating Excitation - variable excitation flip angel to compensate inflow saturation effects in 3D MRA - selectable on desired flow direction and speed - Gradient Motion rephasing permitting effective reduction of flow artifacts <p>Standard Motion Correction</p> <ul style="list-style-type: none"> - <i>syngo</i> Blade - improves image quality by minimizing and correcting for the effects of motion during an MR sequence acquisition. e.g. head, spine, orthopedic imaging and the abdomen - 1D PACE (Prospective Acquisition CorrEction) allows examination of patients with free breathing - 2D PACE (Precise Motion Correction) detects and corrects respiratory motion eg of the heart or liver <p>MAGNETOM Aera runs <i>syngo</i> MR software. <i>syngo</i>® is the unique software platform for medical applications. Parallel working and one-click exams are efficiently supported and increase productivity. Parallel scanning and reconstruction are standard.</p> <p>The unique Phoenix technique is the easiest way to exchange protocol data. It supports intelligent extraction of sequence parameters from images acquired on a MAGNETOM Aera system.</p> <p>Inline technologies, scan@center or AutoVoiceCommands speed up the workflow further.</p> <p>The context-sensitive "Online Help" function and <i>syngo</i> Scan Assistant offer support and propose solutions to MR-specific questions and parameter conflicts.</p> <p>Studies can be easily networked and managed using the standard DICOM 3.0 protocol for efficient support of workflow. The following standard functions are supported: Send/Receive, Query/Retrieve, and Basic Print for DICOM-compatible laser cameras (camera is not included in the basic unit), DICOM Worklist, DICOM Storage Commitment (SC) DICOM Modality Perform Procedure Step (MPPS), DICOM Structured Report (SR), DICOM Study Split.</p> <p>Patient Communication</p> <ul style="list-style-type: none"> - The intercom system includes an ergonomically designed patient communication unit for desktop positioning on the <i>syngo</i> Acquisition Workplace and pneumatic headphones for the patient. - Active Noise Cancellation allows for increased user comfort in the control room combined with comprehensive patient supervision. - It controls emergency table stop, volume control of speaker and headphones in the examination room, volume control of speaker in the control room, response to the patient's activation of the assistance-call button and provides a connection to an external audio system (external audio system is not included in the basic unit) for music playback. <p>Computer system</p> <p>The high performance host computer and the new high performance measurement and reconstruction system are ideally suited for even the most demanding applications. The PC-based computer system uses the intuitive <i>syngo</i> MR user interface. The computer system includes the following components: High-performance measurement and reconstruction system</p> |

| Part No. / Product | Description |
|---|---|
| <p>(Continued) 14416900 MAGNETOM Aera - System</p> | <ul style="list-style-type: none"> - Two Intel Quadcore Processor ≥ E 5540 - clock rate of ≥ 2 × 2.53 GHz - Main memory (RAM) of 48 GB, - Hard disk for raw data ≥ 300 GB - Hard disk for system software ≥ 100 GB - Parallel Scanning and Reconstruction of up to 8 data sets - Reconstruction speed <ul style="list-style-type: none"> - 12.195 recons per second (256 x 256 FFT, full FoV) - 37.914 recons per second (256 x 256 FFT, 25 % recFoV) <p>High-performance host computer</p> |

SIEMENS

Siemens Medical Solutions USA, Inc.
 51 Valley Stream Parkway, Malvern, PA 19355
 Fax: (866) 309-6967

SIEMENS REPRESENTATIVE
 Stuart Wadley - (919) 605-9227

| Part No. / Product | Description |
|--|--|
| | <ul style="list-style-type: none"> - Intel Xeon processor ≥ W3520 QuadCore - clock rate ≥ 2.66 GHz - Main Memory (RAM) ≥ 4 GB - three hard disks <ul style="list-style-type: none"> - system SW ≥ 146 GB SAS - data base ≥ 146 GB SAS - images ≥ 146 GB SAS - DVD-R writer for CD-R (approx. 4000 images 256² DICOM Standard, ISO 9660) and DVD-R (approx. 25 000 images 256² DICOM Standard, ISO 9660) storage of DICOM data or other data like AVI files <ul style="list-style-type: none"> - DVD-ROM drive - Electronic mouse. - The combination of host computer and the measurement and reconstruction system offers a truly powerful imaging system designed for large image matrix sizes of up to 1024 x 1024. The unrestricted multitasking capability allows time-saving parallel scanning and reconstruction. - High-resolution 19" color LCD flat screen monitor with 1280 x 1024 pixel display, integrated gamma correction for optimum display of radiographic grayscale images and automatic backlight control for long term brightness stability. <p>Installation:</p> <ul style="list-style-type: none"> - The relatively lightweight design of the MAGNETOM Aera in most cases eliminates the need for structural building reinforcements and thus facilitates installation in upper floors. - The compact integrated design allows for short installation times and reduces the required space to less than 30 sqm (323 sq. ft.) for the entire installation. The minimum room height clearance is only 2.40 m (7' 10"). - MAGNETOM Aera allows siting of the system without a dedicated computer room - no additional cooling or floor requirements. - MAGNETOM Aera combines state-of-the-art performance with peace of mind. High system availability is ensured by the expert, highly trained Siemens MR service engineers; - Your Siemens service contract (not included in the basic unit) offers a comprehensive range of benefits such as Uptime Remote Diagnostics for improved productivity and maximum uptime. |
| <p>14416902 Tim [204x48] XQ Gradients #Ae</p> | <p>Tim [204x48] performance level Tim 4G offers DirectRF a completely redesigned RF architecture. This new all.digital-in/ digital-out design integrates all RF transmit and receive components at the magnet, eliminating analog cables for true signal purity. This compact and efficient design enables a dynamic feedback control for temporal stability and power linearity. The all-new innovative coil architecture packs more coil elements in a smaller space. Therefore up to 204 coil elements can be simultaneously connected. The newly designed ultra high density array is an essential part supplementing Tim4G. Combined with the 48 independent RF channels advanced IPAT capabilities and SNR are enabled. An additional benefit of multiple coil elements and receiver channels is improved performance in multi-directional, i.e. three dimensional, high-speed, high-resolution IPAT in the head-feet, anterior-posterior or left-right directions.</p> <p>This option includes also Advanced High Order Shim.</p> <p>XQ gradients</p> |

| Part No. / Product | Description |
|--|--|
| <p><i>(Continued)</i> 14416902 Tim [204x48] XQ Gradients #Ae</p> | <p>Siemens XQ gradients provide actively shielded, water cooled world-class gradients. All axes are force-compensated.</p> <p>The XQ gradients have:</p> <ul style="list-style-type: none"> - Maximum gradient amplitude of 45 mT/m, per axis, i.e. 78 mT/m vector summation gradient performance, - max. slew rate 200 T/m/s per axis, i.e. 346 T/m/s vector summation, |

SIEMENS

Siemens Medical Solutions USA, Inc.
51 Valley Stream Parkway, Malvern, PA 19355
Fax: (866) 309-6967

SIEMENS REPRESENTATIVE
Stuart Waddey - (919) 605-9227

| Part No. / Product | Description |
|---|--|
| | <ul style="list-style-type: none"> - minimal rise time 225 μs, from 0 to 45 mT/m amplitude - Max. output voltage for each of the gradient axes 2250 V - Max. output current for each of the gradient axes 900 A - Separate cooling channels that simultaneously cool primary and secondary coils allow the application of extremely gradient intensive techniques in a new class of performance. - 100% duty cycle for fast and demanding techniques such as ultra-short TE MRA in continuous operation, thin slice single breath-hold liver studies and EPI imaging techniques (all optional in appropriate clinical packages). - Variable Field-of-View selection from 0.5 cm to 50 cm (up to 45 cm in z direction) for optimal coverage and highest spatial resolution in diagnostic. The minimum slice thickness in 2D and 3D is 0.1 mm and 0.05 mm, respectively. - Acquisition of sagittal, transverse, coronal, single oblique and double oblique slices with highest resolution. - The extremely compact water-cooled gradient amplifier features a modular expandable design with excellent linearity and pulse reproducibility. It is digitally controlled and has very low switching losses due to ultrafast solid state technology. |
| 14416906 Tim Dockable Table #Ae | <p>The new MAGNETOM Aera table with its light appealing design allows for a fast patient preparation and maximized patient comfort. It provides unobstructed foot space for attending staff and direct access to the patient. The patient table can be lowered to a minimum height of 56 cm (18.5") from the floor, for easier moving of immobile patients and better access for geriatric, pediatric patients or immobile patients. The Tim Table can be moved with two clicks into the isocenter - one click to the upmost position and one click into the isocenter. The tabletop travels beyond the rear end of the system, enabling additional patient access.</p> <p>Multiple Tim4G coils can be connected at once for efficient patient set up and patient friendly examinations. The seamless integration of multiple Tim4G coils is possible via 4 SlideConnect and 4 DirectConnect connector slots, which are embedded in the table. This allows for comprehensive examinations without the need of repositioning.</p> <p>The Tim Dockable Table is easily adjustable for height even in the undocked state. A minimum height of 61 cm allows for easy wheelchair access or easy patient movement to the hospital bed. The integrated infusion stand and arm rests allow for fast patient set up anywhere and also for critical patients</p> |
| 08464872 PC Keyboard US english #Tim | <p>The keys of the numerical key panel are assigned to syngo-specific functions and labeled with the corresponding syngo icons. The keyboard supports the country specific special characters.</p> |
| 14416914 Pure White Design #T+D | <p>The unique color and material selection enhances the visual appeal of the new system design, thereby creating an enticing, patient-friendly impression. The Dot Control Centers and the unique Dot Display are neatly integrated into this main face plate. The aesthetically pleasing and ergonomically designed control elements of the Dot Control Centers are well illuminated for easy visual recognition. In particular, the table cover and the asymmetric left deco area cover have also been designed to promote a modern visual appearance. This combination of ingenuity and practical design as presented with "Pure White" design with its brilliant white and the silver trim simply makes the MAGNETOM an overall visually appealing system and creates a patient-friendly environment.</p> |
| 08857828 UPS Cable #Tim | <p>Power cable to connect the 3 KVA Powerware 9125 small UPS system (pn PWR9125H3000) to the ACC cabinet of the MAGNETOM Avanto/ Espree/ Tim Trio for backing up the host computer and imager.</p> <p>Configuration includes connection box.</p> |

| Part No. / Product | Description |
|--------------------|-----------------------------------|
| <i>(Continued)</i> | The standard cable length is 9 m. |

SIEMENS

Siemens Medical Solutions USA, Inc.
51 Valley Stream Parkway, Malvern, PA 19355
Fax: (866) 309-6967

SIEMENS REPRESENTATIVE
Stuart Waddey - (919) 605-9227

| Part No. / Product | Description |
|--|---|
| 08857828 UPS Cable #Tim | |
| 14413662 UPS Powerware PW9130G-3000T- XLEU | Voltage range: 180 - 276 V Input frequency: 50 / 60 Hz Output voltage: 230 VAC Dimensions (H x W x D): UPS 346 x 214 x 412 mm incl. UPS bracket set Weight: approx. 36 kg |
| 14407258 MR Workplace Table 1.2m | The table design matches the MED-wide uniform design with silver-finished rim, use of friendly colors matching the Siemens color pattern for MAGNETOM and SOMATOM. - Width 120 cm - Depth 80 cm - Height 72 cm |
| 14407261 MR Workplace Container, 50cm | The table design matches the MED-wide uniform design with silver-finished rim, use of friendly colors matching the Siemens color pattern for MAGNETOM and SOMATOM. Table height 72 cm, matching the <i>syngo</i> Acquisition Workplace and <i>syngo</i> MR Workplace console table, for installation in the operator room either directly to the left or right of the <i>syngo</i> Acquisition Workplace or <i>syngo</i> MR Workplace console table or separately. - Width 50 cm - Depth 80 cm - Height 72 cm Alternatively this casing is also suited for the Recon image processor (except for the MR systems with the Tim generation: there the Recon image processor is always placed inside the electronics cabinet). |
| 4MR5142869 Armrest #MR | An MR-compatible arm rest that supports the patient's arm on the magnet patient table when starting intravenous lines. The board is removed after the IV is inserted. This product has been tested and verified for compatibility with the following Siemens' products: MAGNETOM Trio, Verio, Espree, Essenza, Avanto, Symphony, Area Skyra and Biograph mMR. Compatibility with other products cannot be assured and may void service contracts and/or system warranties. |
| CHILINST_AVT Chiller Start-up and Warranty for TIM | Start up and initial set up service performed by the chiller manufacturer or designated service representative. This service does not include the piping and other prerequisite siting, of the waterchiller, which are the responsibility of the customer. 12 months warranty and performed by the chiller manufacturer. |
| MRWE MR Wall sign - English | Highly durable 1mm PVC wall signs with high-tack, double-back tape. Sticks to most any surface. English. 12" x 18". |
| 14402592 Inline Composing syngo #Tim | The Inline Composing option includes the following functions: - Inline calculation of full-format images of the spine, the central nervous system or the vessel tree, for example, combined from multiple overlapping steps. - Dedicated composing algorithms, optimized for the generation of anatomical or angiographic full-format images. - Data sets with different FoV, resolution, matrix and slice thickness can be combined. - Generation of full-format images from inline-computed MIPs. - Different inline functions can be combined; e.g. in case of multiple-step angios, Inline subtraction, Inline MIP and Inline Composing can be performed fully automatically. - Full-format acquisitions from Inline Composing are ideal for further measurement planning on large FoV, e.g. with the Tim Planning Suite (optional, urgently recommended). <i>Prerequisite: Software syngo MR B13.</i> |
| 14402593 | - Easy planning on a FoV of any desired size (up to 205 cm). |

SIEMENS

Siemens Medical Solutions USA, Inc.
 51 Valley Stream Parkway, Malvern, PA 19355
 Fax: (866) 309-6967

SIEMENS REPRESENTATIVE
 Stuart Waddey - (919) 605-9227

| Part No. / Product | Description |
|---|--|
| Tim Planning Suite | <ul style="list-style-type: none"> - Planning of multiple steps simultaneously, e.g. on a whole-body image, with only one Set-n-Go protocol - which includes several steps. - Tim Planning Suite UI: Dedicated user interface and exclusive tools for effective and smooth working on a large FoV. - Multiple slice groups with their overlap are displayed together and can be easily arranged. - All steps can have independent sets of parameters. - All steps are displayed together with a single mouse click. - Easy positioning of all steps, for example, through Align FoV. - Full support of . - Full support of Phoenix, thus maximum reproducibility, for example, for follow-up studies, multi-centric studies or exchange of experiences across different institutions. - Dedicated protocols are provided for the Tim Planning Suite, for example, for orthopedic, oncological or angiographic indications. - Inline Composing for optimized workflow for the generation of full-format images of anatomic or angiographic data sets is a prerequisite. Efficient measurement planning on these full-format images with Tim Planning Suite. - It is highly recommendable to order application training! <p><i>Prerequisite: Software syngo MR B13</i></p> |
| 14416923 Abdomen Dot Engine #T+D | <p>Abdomen Dot Engine Guidance view</p> <ul style="list-style-type: none"> - Step-by-step user guidance is seamlessly integrated. - Example images and guidance text displayed for each step of scanning workflow. - Both images and text are easily configurable by the user <p>Patient View</p> <ul style="list-style-type: none"> - Easily tailored to the individual patient. - Several pre-defined, integrated Dot Exam Strategies are included - Single click update of queue and the complete scan set-up. - Integrated contrast media protocols (Vibe Dynamic) <p>Parameter View</p> <ul style="list-style-type: none"> - A new view that displays the essential parameters - Can be opened at any time during an examination <p>Automatic sequence scaling</p> <ul style="list-style-type: none"> - Auto FoV: optimal FoV is proposed, based on the localizer images. - AutoNavigator: based on automatic breathing pattern detection and scaling of triggered scans. - Breath-hold adaptations <p>Dot Exam Strategies Personalize to the individual patient condition and clinical need.</p> <ul style="list-style-type: none"> - Predefined strategies: <ul style="list-style-type: none"> - Standard with breath-hold - Standard with PACE triggering - Limited patient capabilities using syngo BLADE and PACE triggering. <p>Dot Decisions Seamlessly integrated into scanning workflow:</p> <ul style="list-style-type: none"> - Select the queue and the appropriate protocol or set of protocols are automatically added. - Abdomen Dot Engine integrates MRCP and Diffusion decision points. <p>Timeline setup and monitoring Convenient visual overview of multi-phase breath-hold examinations and CM enhancement curve visualization.</p> <p>Auto Voice Commands</p> <ul style="list-style-type: none"> - Played automatically |

SIEMENS

Siemens Medical Solutions USA, Inc.
 51 Valley Stream Parkway, Malvern, PA 19355
 Fax: (866) 309-6967

SIEMENS REPRESENTATIVE
 Stuart Waddey - (919) 605-9227

| Part No. / Product | Description |
|--|---|
| <p><i>(Continued)</i> 14416923 Abdomen Dot Engine #T+D</p> | <ul style="list-style-type: none"> - Facilitate timing of scanning, breathing and contrast media. - The user controls breath-hold or pauses are actually played - Ability to add pauses between automatic breath-holds. <p>Auto Bolus Detection</p> <ul style="list-style-type: none"> - Automatically initiates the dynamic upper abdomen examination based on bolus detection. - The user can override this function. <p>Inline radial range calculation for MRCP</p> <ul style="list-style-type: none"> - MRCP is measured - Inline Radial Ranges are automatically generated. <p>Inline Subtraction Automatically subtracts the native (non-contrast) measurement from the arterial, portal-venous and late phase.</p> <p>Inline Registration The system automatically performs a registration / alignment of the anatomy for the different dynamic phases, of interest when examining nodular enhancing pathologies.</p> <p>Customization Existing Dot Engines can be modified by the user to their individual standard of care.</p> <ul style="list-style-type: none"> - Add / remove protocol steps - Change guidance content (images and text) - Change or add Dot Exam Strategies and Decision Points - Modify the Parameter View |
| <p>08464740 Flow Quantification #Tim</p> | <p>Flow Quantification enables the acquisition of flow encoded images and the evaluation of blood as well as of cerebro-spinal fluid (CSF).</p> <p>Sequences include:</p> <ul style="list-style-type: none"> - ECG triggered 2D phase contrast with iPAT support - Retrospective reconstruction algorithms for full R-R interval coverage - Maxwell Term Compensation |
| <p>07365419 Argus Flow</p> | <p>The combination of automated contouring and easy-to-use editing tools, provides users with a rapid way to quantify flow parameters.</p> <p>Argus Flow includes:</p> <ul style="list-style-type: none"> - Calculation of flow and velocity parameters(e.g. mean and max velocity, mean, cumulative, prograde, retrograde flow) for large and small vessels. - Semi-automatic detection of regions of interest over time - Color-coded display of velocity values - Calculation of flow and velocity parameters (e.g. peak velocity, average velocity, flow, integral flow) - Graphical and tabular display of the results (e.g. flow-time curves) - Integration of the results in Argus structured report and storage in DICOM format for documentation. |
| <p>14418563 Neuro Perfusion Evaluation,USA #T+D</p> | <p>Post-processing features:</p> <ul style="list-style-type: none"> - Flexible selection of the Arterial Input Function (AIF) by the user. - Pixelwise calculation of the hemodynamic parameters relative Mean Transit Time (relMTT), relative Cerebral Blood Volume (relCBV), relative Cerebral Blood Flow (relCBF), corrected relative Cerebral and Blood Flow (relCBF) for compensation of blood brain barrier leakage. - Pixelwise calculation of maximum signal loss due to contrast agent enhancement (Percentage of Baseline at Peak, PBP) and of the time to the maximum signal loss (Time-To-Peak, TTP). |

SIEMENS

Siemens Medical Solutions USA, Inc.
51 Valley Stream Parkway, Malvern, PA 19355
Fax: (866) 309-6967

SIEMENS REPRESENTATIVE
Stuart Waddey - (919) 605-9227

| Part No. / Product | Description |
|--------------------|--|
| | - Display of the global signal time course (averaged over all slices) to assess the quality of the exam. |

| Part No. / Product | Description |
|--|---|
| <i>(Continued)</i> 14418563 Neuro Perfusion Evaluation, USA #T+D | <ul style="list-style-type: none"> - Predefined post-processing protocols available, user definable post-processing protocol are possible. <p>Visualization features:</p> <ul style="list-style-type: none"> - Colored display of $relMTT$-, $relCBV$-, $relCBF$-, $relCBFcor$, PBP- and TTP-maps. - Zoom, pan, annotate. - Colored images can be saved as DICOM images. |

| | |
|------------------------------|---|
| 14416944 DTI Package #T+D | <p>Diffusion Tensor Imaging Diffusion Tensor Imaging allows for a complete description of the diffusion properties of the brain within the scope of the tensor diffusion model, both for anisotropic and isotropic diffusion. Efficient diffusion direction schemes are pre-defined to allow for optimal diffusion directional resolution. Schemes with up to 256 directions can be selected. Inline technology enables automatic and immediate calculation of the diffusion tensor, including grey-scale and colored "fractional anisotropy" (FA) map derived from it.</p> <p>Details:</p> <ul style="list-style-type: none"> - Measurements with up to 256 different directions and with up to 16 different b-values - Inline calculation of tensor, grey-scale and colored FA map, ADC map and trace-weighted image - Support of parallel imaging (IPAT) - Clinical protocols with full head coverage, incl. inline calculation of tensor, FA, ADC and trace-weighted images in 4 minutes. <p>DTI Tractography syngo <i>syngo</i> DTI Tractography is optimized for the clinical use by providing advanced 3D visualization of white matter tracts in the context of 2D or 3D anatomical datasets and DTI datasets. DTI data sets can be explored fast and intuitively using the interactive QuickTracking. QuickTracking instantaneously displays the tract originating from the mouse pointer position while moving over the DTI data set. This also allows identifying qualified regions to place seeding ROIs. Seed points can be set to assess connectivity by tracking with single ROI and with multiple ROIs. Furthermore they can be placed in fused views displaying the anatomical reference and e.g. the colored FA map simultaneously.</p> <p>Texture Diffusion, a highly versatile in-plane visualization of white matter tracts, allows to display and read DTI Tractography results on PACS reading stations and in the OR.</p> <p>At the same time the package provides the scientific user with the flexibility to configure the tracking algorithm and to change display settings for the tracts. Tract and seeding ROI statistics are included to support publications (e.g. mean/max FA value, min/mean/max ADC value).</p> <p>All views can be exported as DICOM images or bitmaps. Tract and seeding ROI statistics can be exported as html files.</p> <p>DTI Evaluation Clinical applications are supported by a dedicated DTI evaluation mode to support diagnostics of white matter diseases (e.g. multiple sclerosis and brain maturation disorders). Based on the tensor, in addition to the already inline-calculated parameter maps, further maps characterizing the anisotropy of diffusion properties can be calculated and stored. Multiple diffusion parameter maps (e.g. Fractional Anisotropy, ADC, $b=0$) and an anatomical image are displayed next to each other in the same slice position for comparison. The images can be evaluated together based on ROIs and the results can be documented in a table. The display options include 2D and 3D tensor graphics, colour-coded images and overlay images on the anatomical images.</p> <p>In addition, the package offers the scientific user full flexibility of 2- and 3-dimensional visualization of the diffusion tensor with measures of isotropic and anisotropic (fractional and relative) diffusion, Eigen vectors (E1, E2, E3) of the diffusion tensor and shape-descriptive measures of the diffusion tensor (linear, planar, spherical).</p> |
|------------------------------|---|

| | |
|---------------------------------------|---|
| 14409110 Arterial Spin Labeling 2D | The measurement of relative CBF is done by combination of an ASL preparation module with the 2D multi-slice EPI sequence with full IPAT compatibility. Technically, the ASL module consists of a spatially selective inversion pulse combined with different types of saturation pulses (slice presaturation, label slab saturation) and can be |
|---------------------------------------|---|

SIEMENS

Siemens Medical Solutions USA, Inc.
51 Valley Stream Parkway, Malvern, PA 19355
Fax: (866) 309-6967

SIEMENS REPRESENTATIVE
Stuart Waddey - (919) 605-9227

| Part No. / Product | Description |
|------------------------------------|--|
| | classified under PASL (Pulsed Arterial Spin Labeling). Output quantities are a series of labeled and unlabeled images, a series of inline calculated perfusion weighted color maps and a series of relative cerebral blood flow (relCBF) color maps. Prospective motion correction and spatial filtering can be applied for inline calculation to improve the image quality. |
| 14402527 SWI #Tim | Despite a strong sensitivity for local magnetic field inhomogeneities Susceptibility Weighted Imaging (SWI) as a 3D technology keeps up the signal near large susceptibility leaps due to very thin slices and high resolution in the slice |

| Part No. / Product | Description |
|--|--|
| <i>(Continued)</i> 14402527 SWI #Tim | <p>(high image quality e.g. in the area of the forebrain near the frontal sinus). Moreover, the phase information of the MR signal is integrated in the image display. In order to further increase sensitivity for localized microscopic magnetic field inhomogeneities, large-area magnetic field inhomogeneities (e.g. caused by susceptibility leaps near the sinus) are specifically suppressed in the phase images. This allows even small amounts of deoxygenated hemoglobin (e.g. in cerebral veins) or from products of hemoglobin decomposition (e.g. from hemorrhages) to be displayed. Interesting measuring times for the ultra-high-resolution 3D protocols are achieved through parallel imaging with iPAT (GRAPPA).</p> <p>The Susceptibility Weighted Imaging package includes:</p> <ul style="list-style-type: none"> - SWI measuring sequence, iPAT compatible - optimized measuring protocols for the head - inline-postprocessing for automatic calculation of relevant images within the scope of image reconstruction: <ul style="list-style-type: none"> - calculation of susceptibility-weighted images - venous angiography: MIP of a thin slice block <p>SWI has been optimized for clinical use to support diagnostics with cerebrovascular diseases (e.g. cerebral insult), venous malformation, brain trauma and tumors.</p> <p><i>Prerequisite: Software syngo MR B13</i></p> |

| | |
|---|---|
| 14416941 Spectroscopy Package #T+D | <p>The Single Voxel Spectroscopy option is used to measure proton spectra from single voxels. The spectra may show alterations in brain metabolism e.g. in brain tumors, in degenerative changes of the brain and in metabolic diseases. The possibility of automatic adjustment, measurement and evaluation permits near automatic spectroscopy measurements. The whole procedure, including the evaluation of the spectra using the mandatory spectroscopy evaluation option, takes approx. 6 minutes and can be done by doctors or technologists.</p> <p>The 2D Chemical Shift Imaging option is used to measure 2D proton spectroscopic data to generate metabolite images e.g. in brain tumors, metabolic diseases of the brain and degenerative changes in brain metabolism. The whole procedure, including the generation of metabolite images using the spectroscopy evaluation takes approximately 8 minutes.</p> <p>The 3D Chemical Shift Imaging option is used to measure 3D proton spectroscopic data and allows for the evaluation of the spectra in measured volumes and the generation of metabolite images and spectral maps, e.g. in cases of brain tumors, metabolic diseases of the brain and degenerative changes in brain metabolism. The whole procedure, including the generation of metabolite images using the spectroscopy evaluation takes approximately 10-16 minutes</p> <p>Optimized protocols for 3D CSI in the prostate are also included.</p> <p>The evaluation software is fully integrated in <i>syngo MR</i>. Evaluation protocols adapted to the scan protocols carry out a complete and automatic evaluation of the measured data.</p> <p>The following functions are included:</p> <ul style="list-style-type: none"> - Subsequent water suppression with optional phase correction - Apodization - Zero filling - Fourier transformation - Base line correction |
|---|---|

SIEMENS

Siemens Medical Solutions USA, Inc.
51 Valley Stream Parkway, Malvern, PA 19355
Fax: (866) 309-6967

SIEMENS REPRESENTATIVE
Stuart Waddey - (919) 605-9227

| Part No. / Product | Description |
|--------------------|---|
| | <ul style="list-style-type: none"> - Automatic or manual phase correction - Curve fitting and peak labeling - Summaries in tabular form of the essential results specifying the metabolites, their position, integrals and signal ratios in relation to a selectable reference. - Capability of exporting spectroscopy header information and data into a documented external format. - Automated peak normalization to tissue, water or reference. <p>For CSI the following functions are included:</p> <ul style="list-style-type: none"> - Spectra of selected voxels are automatically calculated, corrected for possible B0 deviations and displayed. - Spectral fit is automatically optimized for each voxel. - CSI data can be represented as spectral maps and colored metabolite images that can be superposed onto |

| Part No. / Product | Description |
|--|---|
| <i>(Continued)</i> 14416941 Spectroscopy Package #T+D | anatomical images. |
| 14405328 TWIST syngo #Tim | <p><i>syngo</i> TWIST provides:</p> <ul style="list-style-type: none"> - Visualization of contrast agent dynamics in the vessel system of interest with maximum flexibility. - Needs only a low amount of contrast agent. - Imaging in all body regions, e.g. carotids, pulmonary and peripheral vessels with brilliant spatial and temporal resolution. - Clear separation of the arterial and venous phase. - High speed acquisition by intelligent k-space strategies and use of iPAT, powered by Tim. - <i>syngo</i> TWIST provides fat suppression using water selective excitation. - Inline technologies, such as subtraction and MIP are provided for optimal workflow. - In case of very high spatial resolution <i>syngo</i> TWIST may even replace conventional static MR angio. Moreover, <i>syngo</i> TWIST does not require any bolus timing - just inject and go. |
| 14416955 Body 18 #Ae | <p>The Body 18 has a 18-element design with 18 integrated preamplifiers that are arranged in 3 clusters of 6 coil elements each. The Body 18 will be typically used together with the Spine 32 with which it operates in an integrated fashion as a 30-element coil, resulting in 3 rings of 10 elements each for highest SNR and fast imaging. It can be positioned in different orientations and addresses the requirement range for the examinations of obese patient to pediatric patients. The light weight coil improves patient comfort and can be easily connected via SlideConnect technology. No tuning of the fully iPAT-compatible Body 18 is necessary allowing for efficient and patient friendly set-up.</p> <p>For examinations where larger anatomical coverage is required, several Body 18 coils can be used simultaneously. Up to four Body 18 can be used simultaneously, but typically two Body 18 will be used for coverage of the entire abdomen or in the case of large patients.</p> <p>The Body 18 is typically used in combination with the Spine 32 for examinations of the thorax, abdomen, pelvis or hip and is also well suited for cardiac or vascular applications. In addition, the Body 18 can be combined with the Spine 32, further Body 18 (optional), the Peripheral Angio 36 (optional), but also the Head/Neck20 and the 4-channel flex coils (e.g. Flex Large 4, Flex Small 4) it contributes for all large-Field-of-View applications up to whole-body imaging.</p> <p>The dimensions of the Body 18 are 385 mm x 590 mm x 65 mm (L x W x H). Its weight is about 2 kg (4.5 lbs), whereas the patient feels as little weight as 1kg (2,25 lbs).</p> |
| 14405341 Mapit syngo #Tim | <p>Features:</p> <ul style="list-style-type: none"> - 3D VIBE sequence for Inline T1 mapping |

SIEMENS

Siemens Medical Solutions USA, Inc.
51 Valley Stream Parkway, Malvern, PA 19355
Fax: (866) 309-6967

SIEMENS REPRESENTATIVE
Stuart Waddey - (919) 605-9227

| Part No. / Product | Description |
|---|--|
| | <ul style="list-style-type: none"> - Multiecho spin echo sequence for Inline T2 mapping - 3D Multiecho gradient echo sequence for Inline T2* mapping - iPAT compatibility - Protocols for Inline parametric mapping <p>Using iPAT the 3D sequences provide isotropic imaging extremely high resolution while maintaining clinical measurement times. These data sets allow for the multi planar reconstruction of all planes. 3D is necessary to properly visualize the whole articular cartilage since it typically has a complex shape. In addition the accuracy of isotropic high resolution 3D data sets is superior because partial volume effects between e.g. synovial fluid and cartilage are minimized.</p> <p>For the visualization of the parametric maps in the anatomical context the maps can be displayed as a colored overlay onto anatomical images using the optional package "syngo Image Fusion"</p> |
| 07820082 Inline Perfusion #Tim | <p>Inline Technology – Processing Instead of Post-processing.</p> <p>Inline Technology helps to streamline the clinical workflow by automating post-processing steps before image viewing. This facilitates getting clinical results immediately. This package integrates Inline technology with perfusion imaging. Automatic real-time calculation of Global Bolus Plot (GBP), Percentage of Baseline at map (PBP) and Time-to-Peak map (TTP) with Inline technology is possible.</p> |

| Part No. / Product | Description |
|---|---|
| <i>(Continued)</i> 07820082 Inline Perfusion #Tim | An optimized EPI sequence for perfusion-diagnostics is included in the standard Tim Application Suite. With this package real-time calculations are done of anatomical images and, in addition, of a global bolus plot and a Time-to-Peak map for visualizing the time dependence of tissue perfusion. |
| 14416958 Peripheral Angio 36 #Ae | <p>The Peripheral Angio 36 has a 36-element design with 36 integrated preamplifiers distributed over 6 planes with 6 elements each.</p> <p>A uniquely designed non-ferromagnetic coil cart for safe coil storage is included. The PA Matrix Coil is also shipped with a set of positioning cushions for proper handling.</p> <p>No tuning of the fully iPAT-compatible Peripheral Angio 36 is required.</p> <p>With a length of about 1m both legs are covered from the iliac artery level down to the foot arch vessels using multiple, flexible wings. For the visualization of the abdominal aorta and the iliac bifurcation it can be combined with the Body 18 and Spine 32. For larger body coverage eg whole body with up to 205 cm possible coverage, it can be combined with Head/Neck20 or a further Body18 to allow for large Field of View examinations with high patient comfort. Patient set up is done once and no repositioning is necessary</p> <p>For peripheral Angiography the PA Matrix coil will be typically used in feet-first position, but also head-first positioning for whole-body examinations is possible (optional Tim Whole Body Suite required).</p> <p>The dimensions of the Peripheral Angio 36 are: 860 mm x 300 - 640 mm x 280 mm</p> |
| 14416961 Hand/Wrist 16 #Ae | <p>The 16-element coil with 16 integrated pre-amplifiers excels in highest resolution imaging with exceptional signal/noise ratio, while taking full advantage of iPAT in all directions.</p> <p>Hand/Wrist 16 is ergonomically designed and adapted to the shape of the hand/wrist region. The coil features a hinged design of the upper part and slidable attachment to the base plate. Together with the included stabilization pads the coil allows easy, fast and comfortable patient positioning.</p> |
| 14416962 Foot/Ankle 16 #Ae | <p>The 16-element coil with 16 integrated pre-amplifiers excels in highest resolution imaging with exceptional signal/noise ratio, while taking full advantage of iPAT in all directions.</p> <p>Foot/Ankle 16 is ergonomically designed and features a boot-like coil design. Together with the included</p> |

SIEMENS

Siemens Medical Solutions USA, Inc.
51 Valley Stream Parkway, Malvern, PA 19355
Fax: (866) 309-6967

SIEMENS REPRESENTATIVE
Stuart Waddey - (919) 605-9227

| Part No. / Product | Description |
|------------------------------------|--|
| | stabilization pads the coil allows easy, fast and comfortable patient positioning. |
| 14407269 Loop Coil, large #Tim | The coil has a diameter of approx. 11 cm and can be positioned easily. |
| 14407270 Loop Coil, medium #Tim | The coil has a diameter of approx. 7 cm and can be positioned easily. |
| 14407271 Loop Coil, small #Tim | The coil has a diameter of approx. 4 cm and can be positioned easily. |
| 14405224 Composing syngo #Tim | <p>The option features:</p> <ul style="list-style-type: none"> - Display and storage of full-format images, e.g. of the spine, the central nervous system or the vessel tree (starting from <i>syngo</i> MR B13), combined from multiple overlapping stages. - Dedicated composing algorithms, optimized for the generation of anatomical or angiographic (starting from <i>syngo</i> MR B13) full-format images. - Data sets with different FoV, resolution, matrix and slice thickness can be combined (starting from <i>syngo</i> MR B13). - Generation of full-format images from inline MIPs (starting from <i>syngo</i> MR B13). - Original, detail and reconstructed images can be displayed in different layouts. - Comparison of two reconstructed images for evaluation and diagnosis is thus made possible. - Filming in different layouts is supported. - Measurements of basic functions via reconstructed images is then possible. - Measurements of extended orthopedic functions: scoliotic angle, kyphotic angle, vertical distance measurement and differences in width of the intervertebral |

| Part No. / Product | Description |
|--|--|
| (Continued) 14405224 Composing syngo #Tim | <p>spaces.</p> <p><i>Prerequisite: SW syngo MR B13.</i></p> |
| 14409198 Native syngo #Tim | <p><i>syngo</i> NATIVE offers:</p> <ul style="list-style-type: none"> - Non-contrast enhanced MRA - Separate imaging of arteries and veins - Visualization of - e.g. - renal arteries or peripheral vessels <p>The <i>syngo</i> NATIVE package comprises:</p> <ul style="list-style-type: none"> - <i>syngo</i> NATIVE TrueFISP - <i>syngo</i> NATIVE SPACE |
| 14430403 Tx/Rx 15-channel Knee Coil DDST #Ae | <p>Thanks to its 15-channel design this coil is perfectly suited for high-resolution images with excellent SNR. With the arrangement of the antennas in three rings of 5 elements each, the coil is specially designed for parallel imaging with high acceleration factors.</p> <p>The coil is positioned on a laterally movable support and therefore allows for comfortable patient positioning of both legs for off-center examinations. SlideConnect Technology allows for fast and easy patient preparation, resulting in less table time. Furthermore, the upper part can be removed for easier patient positioning. Additional cushions allow for optimum patient immobilization.</p> <p>The integrated transmission function makes volume-sensitive excitation with greatly reduced RF power possible on the one hand and, on the other, prevents aliasing artifacts (e.g. due to the other knee).</p> |
| 14416968 CP Extremity Coil #Ae | <p>The coil is placed on a laterally movable holder and is capable to allow off-center scanning with comfortable positioning of the other leg and has special fixation aids with automatic inflation. The coil may be placed on top of the Spine Matrix Coil.</p> <p>The upper part of the coil can be removed for easy patient positioning and has an opening for examinations of the</p> |

SIEMENS

Siemens Medical Solutions USA, Inc.
51 Valley Stream Parkway, Malvern, PA 19355
Fax: (866) 309-6967

SIEMENS REPRESENTATIVE
Stuart Wadley - (919) 605-9227

| Part No. / Product | Description |
|--|--|
| | ankle. Because of the circular polarization this coil is suited for highest resolution imaging with excellent signal/noise ratio. The integrated transmit functionality allows volume selective excitation with significantly reduced RF-power, and avoids the occurrence of aliasing artifacts (e.g. from the other knee). The inner diameter of the CP Extremity Coil is: 195 mm (min.) |
| 14416960 Shoulder 16 Coil Kit #Ae | The iPAT compatible Shoulder 16 Large and Shoulder 16 Small are ergonomically designed and adapted to the shape of the shoulder. The different sizes obtain maximum image quality for different body sizes: - 165 mm (6.5 in) diameter for small and medium sized shoulders - 200 mm (7.9 in) diameter for large shoulders The coils can be used either for left or right shoulders. It features sliding attachments to the base plate and can easily be adjusted for comfortable positioning. The coils excels in highest resolution imaging with exceptional signal/noise ratio. |
| 14416952 Coil Storage Cart #T+D | The cart may be rolled to convenient locations in the examination room and can be opened up to work like a shelf. The coil storage cart has multiple drawers and trays as well as many other storage spaces for coils, cushions and miscellaneous items. Its dimensions are: Width 140 cm (4' 7") when closed and 280 cm (9' 12") when opened, depth 54 cm (1'9") and height 121 cm (3'12"). |
| 14416975 Eco Chiller 60kW | Chiller KKT ECO 133 - L Function: Supplies dedicated primary chilled water in cases where no chilled water supply is available on site. Air-cooled version, for outdoor installation up to a maximum distance of 25 m for connection to the IFP, incl. 50 m FOC for control. The cooling capacity of the chiller is 60 kW, the chilled water temperature is 20°C, the water flow is 130 l/min. |

| Part No. / Product | Description |
|--|--|
| <i>(Continued)</i> 14416975 Eco Chiller 60kW | Ambient temperature: -20 to +48°C Connection rating: 28 kW Voltage: 3/PE 400 V to 480 V / 50/60 Hz Fuse rate: 80 A Power consumption: 66 A Dimensions: 2000 mm x 1100 mm x 2100 mm (height x width x depth). Weight: 760 kg Noise level at a distance of 10 m at outside temperatures of: 21°C 47 dB(A) 32°C 52 dB(A) 48°C 58 dB(A) IFP (Interface Panel) Main functions of the IFP: - Interface function between the KKT chiller and the MR cabinet. - Water supply for MREF, MBB, CBB and TX box. Additional devices such as integrated differential pressure control, a pressure gage, and a filter are used in order to guarantee the precise functioning of the cooling circuit, especially for the cold head compressor (MREF). The connection must be made locally with 2" lines up to a maximum distance of 25 m. Dimensions: 800 mm x 1150 mm x 210 mm (height x width x depth). Weight: 67 kg |

SIEMENS

Siemens Medical Solutions USA, Inc.
51 Valley Stream Parkway, Malvern, PA 19355
Fax: (866) 309-6967

SIEMENS REPRESENTATIVE
Stuart Waddey - (919) 605-9227

| Part No. / Product | Description |
|--|--|
| 14418521 syngo Expert-i #T+D | <p>The option is integrated in the <i>syngo</i> user interface thus enables easy access to the user interface of the <i>syngo</i> Acquisition Workplace for planning and processing support purposes. The access is protected by appropriate security mechanisms (active enabling prior to every connection through the user present on site, password protection), in order to prevent unwanted connections.</p> <p>The client software can be operated on any commercial PC with the following specification:</p> <ul style="list-style-type: none"> - Operating system: Windows NT or XP - 32 bit graphics board - 850 MHz CPU |
| 14416969 Flex Coil Interface 1.5T | <p>The interface is not permanently mounted and therefore allows free positioning of the flexible coils as required by the examination procedure.</p> |
| M3SSMREPICBC Spectris Solaris EP Injector iCBC (Optional) | <p>Includes Spectris Solaris EP injector and Integrated Continuous Battery Charger (ICBC).</p> <ul style="list-style-type: none"> - Optimized color touch screen with few keystrokes. - Six user-programmable phases for added flexibility. - Independent Keep Vein Open (KVO) allows more time to focus on patient. - Large 115 mL saline syringe allows for longer KVO and multiple flushes. - Design of low pressure tubing eliminates dead space in the "T" connection that can waste contrast. - The clear barrel design with molded FluidDots help detect the presence of air in a syringe. - Pressure Limit Setting control software enables user to select from one to six preset maximum pressure limits, ranging from 100-300 psi, and to view current pressure during injection next to the pre-selected maximum value on the Solaris display. <p>Installation, applications and one year warranty provided by Medrad.</p> <p>Not for mobile use, refer to Siemens part number M3SSMR300EPM for the Solaris injector used in a mobile environment.</p> <p>This product has been tested and verified for compatibility with the following Siemens' products: MAGNETOM Trio, Espree, Essenza, Verio, Avanto, Symphony, Aera, Skyra and Biograph mMR. Compatibility with other products cannot be guaranteed and use with any other products may void service contracts and/or system warranties.</p> |
| 14413612 Tissue 4D syngo #Tim (Optional) | <p>Visualization features:</p> <ul style="list-style-type: none"> - 4D visualization (3D and over time) - Color display of parametric cards (Ktrans, Kep, Ve, Vp, IAUC) - Additional visualization of 2D or 3D morphological dataset <p>Post-processing features:</p> <ul style="list-style-type: none"> - Elastic 3D motion correction - Fully automatic calculation of subtracted images <p>Standard curve evaluation:</p> <ul style="list-style-type: none"> - Calculation and display of enrichment curves <p>Pharmacokinetic model:</p> <ul style="list-style-type: none"> - Pharmacokinetic calculation on a pixel-by-pixel basis using a 2-compartment model - Calculation is based on the Toft model. Various model functions are available. - Manual segmentation and calculation on the result images. <p>The following result images can be saved as DICOM images:</p> <ol style="list-style-type: none"> 1. 3D motion-corrected, dynamic images 2. Colored images 3. Possibility for exporting results in the relevant layout format. |

SIEMENS

Siemens Medical Solutions USA, Inc.
51 Valley Stream Parkway, Malvern, PA 19355
Fax: (866) 309-6967

SIEMENS REPRESENTATIVE
Stuart Wadley - (919) 605-9227

| Part No. / Product | Description |
|---|--|
| 14416929 Advanced Cardiac Package #T+D (Optional) | <p>Combining the unique advantages of Tim and <i>syngo</i> BEAT with iPAT and powerful gradients, it allows performing cardiac MR examinations without compromise in image resolution or acquisition speed. <i>syngo</i> BEAT is a unique tool for fast and easy cardiovascular MR imaging. It provides 1-click switch from cine imaging to tagging for wall motion evaluation and 1-click switch from 2D to 3D imaging. <i>syngo</i> BEAT automatically adjusts all parameters associated with the changes.</p> <p>Cardiac and Vessel Morphology</p> <ul style="list-style-type: none"> - Multi echo technique for e.g. thalassemia assessment - 3D aortopathy imaging with free breathing (SPACE) <p>Global or Regional Wall Motion Analysis with <i>syngo</i> BEAT</p> <ul style="list-style-type: none"> - 3D cine acquisition for full CT-like heart coverage - 2D segmented FLASH for visualization of the regional wall motion using various tagging techniques (grid or stripes) <p>Dynamic myocardial imaging with <i>syngo</i> BEAT</p> <ul style="list-style-type: none"> - Ultra-fast, high-SNR sequence for dynamic imaging with GRE EPI contrast for stress and rest exams <p>Tissue characterization with <i>syngo</i> BEAT</p> <ul style="list-style-type: none"> - Robust myocardial tissue characterization with 3D PSIR (phase-sensitive inversion recovery), e.g. after myocardial infarction or for differentiation of cardiomyopathies - Fast and complete coverage of the myocardium with IR 3D FLASH and TrueFISP <p>Coronary imaging with <i>syngo</i> BEAT</p> <ul style="list-style-type: none"> - 3D Whole-Heart non-contrast Coronary MRA - 3D Whole-Heart MRA with advanced free-breathing navigator compensating diaphragm shifts during the acquisition (motion-adaptive respiratory gating) |
| 14407334 Argus 4D Ventr.Function syngo #Tim (Optional) | <p>This package includes Argus Function as well as Argus 4D Ventricular Function.</p> <p>Argus Function:</p> <ul style="list-style-type: none"> - Automatic, semi-automatic, or manual segmentation of the left and semi-automatic or manual segmentation of the right ventricle. - Volumetric analysis and wall thickness analysis. |

| Part No. / Product | Description |
|---|---|
| <p><i>(Continued)</i></p> 14407334 Argus 4D Ventr.Function syngo #Tim (Optional) | <ul style="list-style-type: none"> - Output of parametric results, volume-time curves and bull's-eye plots. - DICOM Structured Reporting. <p>Argus 4D Ventricular Function:</p> <ul style="list-style-type: none"> - Calculation of volumetric cardiac data of a given patient very quickly and easily. - Parametric results and volume-time curves are calculated upon automatic creation and adaptation of a 4D model of the left ventricle. - The resulting 4D model of the patient's heart can be visualized superimposed to anatomical images as a reference. |
| 14416908 Tim Whole Body Suite #T+D (Optional) | <p>Tim and the Tim Whole Body Suite enable for true whole body MR scanning for head-to-toe imaging. Whole body imaging with highest image quality without patient repositioning and without the need to change a single coil, not even once, this means whole body imaging without compromise.</p> <p>The Tim Whole Body Suite features:</p> <ul style="list-style-type: none"> - The all-new Tim Table or Tim Dockable Table enable a full Field-of-View with coverage up to 205 cm (6' 9"). The table top has the same length as the standard system without whole body capabilities. Additional free space is required at the rear part of the magnet to ensure, that the table movement is not limited by the rear |

SIEMENS

Siemens Medical Solutions USA, Inc.
 51 Valley Stream Parkway, Malvern, PA 19355
 Fax: (866) 309-6967

SIEMENS REPRESENTATIVE
 Stuart Wadley - (919) 605-9227

| Part No. / Product | Description |
|--|---|
| | <p>wall.</p> <ul style="list-style-type: none"> - Table movement to its full extent can be remotely controlled from the operator console either by the operator or by sequence protocols. - Protocols and programs for whole body MR angiography and morphology e.g. for metastasis visualization and preventive care examinations. - Whole body MR Angiography is possible with high speed, high resolution and high image contrast on the entire volume combining high speed gradients and iPAT. - The large FoV of 205 cm supports the assessment of metastases distribution in the body with sequences such as TIRM (Turbo Inversion Recovery). |
| <p>14426310 Angio Dot Engine (Optional)</p> | <p><i>Guidance View</i></p> <ul style="list-style-type: none"> - Step-by-step user guidance is seamlessly integrated. - Example images and guidance text are displayed for each individual step of the scanning workflow. - Both images and text are easily configurable by the user <p><i>Test bolus</i></p> <ul style="list-style-type: none"> - Automatic detection of arterial / venous timing window <p><i>Feedback of bolus timing information</i></p> <ul style="list-style-type: none"> - Timing information is fed back into planning steps and parameters are adapted automatically <p><i>Auto Voice Commands</i></p> <ul style="list-style-type: none"> - Integrated into the scanning workflow. - The system plays them automatically at the right point in time. - This ensures optimal timing of scanning, breathing and contrast media. - The user can monitor which breath hold or pauses are actually played, and could add pauses between the automatic breath hold commands if necessary <p><i>Customization</i></p> <p>Existing Dot Engines can be modified by the user to their individual standard of care.</p> <ul style="list-style-type: none"> - Add/remove protocol steps - Change guidance content (images and text) - Change or add Dot Exam Strategies and Decision Points - Modify the Parameter View <p>Application Packages: <i>syngo</i> Inline Composing</p> <ul style="list-style-type: none"> - Automatic anatomical or angiographic composing of multiple adjacent coronal or sagittal images |

| Part No. / Product | Description |
|--|---|
| <p><i>(Continued)</i> 14426310 Angio Dot Engine (Optional)</p> | <ul style="list-style-type: none"> - Composed images can be automatically loaded into Graphical Slice Positioning for scan planning purposes <p>Tim Planning Suite With the Tim Planning Suite, multiple regions in the entire body can be examined in a minimum of time through measurement planning on a single FoV of any desired size.</p> |
| <p>07365484 Image Fusion syngo (Optional)</p> | <p>CT, MR, NM, or PET images are accepted as input for image fusion. Studies can be done with the same modality or with different modalities.</p> <p>Registration Algorithms:</p> <ul style="list-style-type: none"> - Easy-to-use visual alignment with 6 degrees of freedom (3x translation, 3x rotation) - Landmark based registration with convenient landmark editor for point-based registration using anatomical |

SIEMENS


Siemens Medical Solutions USA, Inc.
51 Valley Stream Parkway, Malvern, PA 19355
Fax: (866) 309-6967

SIEMENS REPRESENTATIVE
Stuart Wadley - (919) 605-9227

| Part No. / Product | Description |
|--------------------|---|
| | <p>landmarks</p> <ul style="list-style-type: none">- Automatic registration- Storage of transformation matrix after registration for later retrieval with datasets <p>Visualization Techniques:</p> <ul style="list-style-type: none">- Side by side visualization of both datasets with correlated pointer and correlated scrolling with dog ears- 2D alpha-blending in monochrome or pseudo-color with adjustable balance between the two superimposed data sets. <p>Storage of fused results as secondary capture images.</p> |

SIEMENS

MAGNETOM AERA 1.5T TYPICAL ROOM PLAN




The intended use for this Cut Sheet is to communicate the spatial requirements as well as the basic architectural, electrical, structural, and mechanical requirements for this piece of imaging equipment. The information provided in this document is for reference only, during the pre-planning stage, and therefore does not contain any site specific detailed requirements. This information is subject to change without notice. Federal, state and/or local requirements may impact the final placement of the components. It is the customer's responsibility to ensure that the final layout and placement of the equipment complies with all applicable requirements.

SIEMENS

FOR REFERENCE ONLY,
NOT FOR CONSTRUCTION.

MAGNETOM AERA 1.5T TYPICAL ROOM PLAN



TYPICAL PLAN

SCALE: 1/8" = 1'-0"

| EQUIPMENT LEGEND | | | | | | | | |
|------------------|---|------------|-----------------|------------------|---------------------|----------|--------|---------------------|
| NO | DESCRIPTION | SMS SYM | WEIGHT (LBS) | BTU/HR TO AIR | DIMENSIONS (INCHES) | | | REMARKS |
| | | | | | W | D | H | |
| ① | MRC OPERATING CONSOLE AND KEYBOARD | Ⓜ | 132 | --- | 45 11/16 | 35 1/4 | 28 3/8 | |
| ② | COLOR MONITOR FOR MRC | Ⓜ | 22 | 239 | 18 5/16 | 16 15/16 | 4 3/4 | ON CONSOLE/COUNTER |
| ③ | HOST PC MRC | Ⓜ | 49 | 2,389 | 11 | 27 | 18 1/8 | |
| ④ | CONTAINER FOR HOST 500 | Ⓜ | 238 | --- | 19 5/8 | 31 1/2 | 28 3/8 | |
| ⑤ | ALARM BOX | Ⓜ | 2 | --- | 9 | 4 | 9 | |
| ⑥ | 1.5T MAGNET WITH COVERS AND PATIENT TABLE | Ⓜ | 10,093 | 7,506 | 91 | 170 | 86 | |
| ⑦ | RF-FILTER PLATE | Ⓜ | 285 | 853 | 46 1/2 | 21 3/4 | 21 1/2 | |
| ⑧ | MAGNET STOP | Ⓜ | 1 | --- | 3 | 5 | 3 | |
| ⑨ | ELECTRONICS CABINET (GPA/EPC CABINET) | Ⓜ | 3,307 | 13,649 | 61 1/2 | 26 | 77 1/2 | |
| ⑩ | SEP CABINET | Ⓜ | 750 | 3,415 | 25 5/8 | 25 5/8 | 73 5/8 | |
| ⑪ | POWERWARE 9130 UPS WITH EBM (OPTION) | Ⓜ | 186 | 1,257* | 16 7/8 | 12 7/8 | 16 1/4 | *1,755 ON BATTERIES |

SIEMENS

FOR REFERENCE ONLY,
NOT FOR CONSTRUCTION.

MAGNETOM AERA 1.5T SPECIFICATIONS

| POWER REQUIREMENTS | |
|--|--------------------|
| VOLTAGE RANGE: 480 VAC $\pm 10\%$ FOR ALL LINE AND LOAD CONDITIONS. VOLTAGE BALANCE: 2% MAXIMUM DIFFERENCE BETWEEN PHASES | |
| FREQUENCY: | 60 Hz ± 1.0 Hz |
| LINE IMPEDENCE: | 95 mOHMS |
| STAND BY POWER: | 7.2 kVA |
| HIGHEST AVERAGE POWER | 56 kVA |
| CONNECTION VALUE (LESS THAN 5 MINUTES) | 110 kVA |
| MOMENTARY POWER | 140 kVA |
| RECOMMENDED TRANSFORMER | 150 kVA |
| MR SYSTEM OVERCURRENT PROTECTION | 175 A |
| RECOMMENDED UPS | TBD |
| UPS SYSTEM OVERCURRENT PROTECTION | TBD |
| MAX. ALLOWABLE VOLTAGE DROP AT MAX. POWER | 6.0% |

| NOISE LEVELS | |
|---|--------------------------|
| SYSTEM ROOM | NOISE LEVEL / dB(A) |
| CONTROL ROOM | <55 |
| EXAMINATION ROOM | 86.1 dB(A) AVERAGE VALUE |
| EQUIPMENT ROOM | <65 |
| IT IS THE CUSTOMER'S RESPONSIBILITY TO ENSURE THAT ALL LOCAL/ STATE/OSHA NOISE REGULATIONS ARE ADHERED TO. ADDITIONAL NOISE DATA MAY BE PROVIDED BY SIEMENS PROJECT MANAGER UPON REQUEST. | |

| POWER REQUIREMENTS |
|--|
| DEMAND AND CAPACITY REQUIREMENTS NOTES |
| 1) IF EQUIPMENT UPGRADE IS ANTICIPATED, INSTALLING ELECTRICAL POWER TO MEET THE REQUIREMENTS OF THE HIGHER POWER GRADIENT PACKAGE AT THE TIME OF INITIAL INSTALLATION WILL REDUCE THE COST TO UPGRADE THE ELECTRICAL SYSTEM LATER. |
| 2) RECOMMENDED TRANSFORMER SIZE (SYSTEM WITHOUT UPS) IS BASED ON INDUSTRY STANDARD ISOLATION TRANSFORMER KVA RATINGS. SOURCE IMPEDANCE FEEDING THE MAGNETOM SYSTEM, INCLUDING ANY ISOLATION TRANSFORMERS, MUST MEET EQUIPMENT REQUIREMENTS AS LISTED HERE. SIEMENS RECOMMENDS A TRANSFORMER WITH COPPER WINDINGS, AN ELECTRO-STATIC SHIELD, AND A LOW IMPEDANCE (<3%) TO ENSURE THAT SOURCE IMPEDANCE REQUIREMENTS ARE MET. |
| 3) OVERCURRENT PROTECTION IS SPECIFIED FOR SYSTEMS WITHOUT AN UNINTERRUPTIBLE POWER SUPPLY (UPS). ADDITION OF A UPS REQUIRES A HIGHER CAPACITY MAINS CONNECTION (DEPENDENT UPON UPS MODEL AND SIZE). MAXIMUM FAULT CURRENT IS DEPENDENT UPON THE IMPEDANCE OF THE FACILITY ELECTRICAL SYSTEM. CUSTOMER'S ARCHITECT OR ELECTRICAL CONTRACTOR TO SPECIFY AIC RATING OF OVERCURRENT PROTECTION BASED ON FACILITY IMPEDANCE CHARACTERISTICS. |
| 4) MOMENTARY POWER IS BASED ON A MAXIMUM RMS VALUE FOR A PERIOD NOT TO EXCEED FIVE (5) SECONDS, AS DEFINED IN NEC 517.2. STAND-BY AND AVERAGE CURRENT ARE SUBSTANTIALLY LOWER. |
| 5) THE CONDUCTOR SIZE SHOULD BE SELECTED TO MEET THE VOLTAGE DROP REQUIREMENTS, TAKING INTO CONSIDERATION THE MAINS CAPACITY, RUN LENGTH, AND ANY ADDITIONAL TRANSFORMERS USED TO OBTAIN THE PROPER EQUIPMENT VOLTAGE LEVEL. NEMA STANDARD XR-9-1989 (R1994,R2000) PROVIDES GENERAL GUIDELINES FOR SIZING CONDUCTORS, TRANSFORMERS, AND ELECTRICAL SYSTEMS FOR MEDICAL IMAGING SYSTEMS. |
| 6) LONG-TIME POWER IS BASED ON THE HIGHEST AVERAGE RMS VALUES FOR A PERIOD EXCEEDING 5 MINUTES DURING CLINICAL SYSTEM OPERATION, AS DEFINED IN NEC 517.2. |
| 7) A CIRCUIT BREAKER WITH A HIGH INRUSH RATING (>8x RATED CURRENT) IS REQUIRED TO PERMIT SWITCH-ON OF THE UPS SYSTEM WITHOUT SPURIOUS TRIPPING. CIRCUIT BREAKERS WITH AN ADJUSTABLE MAGNETIC TRIP (SIEMENS FD6 SERIES OR SIMILAR) ARE HIGHLY RECOMMENDED. |

| CEILING HEIGHTS |
|------------------------------|
| EXAM ROOM 7'-11" MINIMUM |
| CONTROL ROOM 6'-11" MINIMUM |
| EQUIPMENT ROOM 7'-3" MINIMUM |

| REMOTE SYSTEM DIAGNOSTICS |
|--|
| SIEMENS REMOTE SERVICES (SRS) REQUIRES A CONNECTION BETWEEN THE SRS REMOTE SERVER AND SIEMENS SYSTEMS VIA REMOTE LOCAL AREA NETWORK ACCESS, TO ENSURE THE UPTIME OF YOUR SYSTEM. |
| THIS SERVICE REQUIRES ONE OF THE FOLLOWING CONNECTION METHODS: |
| 1. (PREFERRED) VPN - WHERE THE CUSTOMER HAS AVAILABLE A VPN CAPABLE FIREWALL OR OTHER VPN APPLIANCE. |
| 2. (OPTIONAL) *SRS ROUTER* - CONNECTED TO ANALOG PHONE LINE VIA *ANALOG MODEM*, ETHERNET CONNECTION TO CUSTOMER'S LAN, AND A POWER OUTLET. |
| NOTE: = *SUPPLIED BY SIEMENS* |

| FOR MORE INFORMATION |
|--|
| FOR MORE DETAILED PLANNING REQUIREMENTS FOR THIS SYSTEM, SEE THE TYPICAL FINAL DRAWING SET NUMBER: 10023 |

SIEMENS

FOR REFERENCE ONLY,
NOT FOR CONSTRUCTION.

MAGNETOM AERA 1.5T SPECIFICATIONS

CHILLED WATER SUPPLY

A CHILLED WATER SUPPLY IS REQUIRED TO THE MRI SYSTEM 24 HOURS A DAY, YEAR ROUND FOR THE COLD HEAD AND GRADIENT SYSTEMS. THIS CAN BE PROVIDED BY A CENTRAL CHILLED WATER SUPPLY OR A SEPARATE STAND ALONE CHILLER THAT MEETS THE STATED REQUIREMENTS. THE CHILLED WATER CAN ALSO BE SUPPLIED BY A DEDICATED KRAUS ECO CHILLER AND INTERFACE PANEL.

WITHOUT THE USE OF A DEDICATED KRAUS CHILLER, A SEP (SYSTEM SEPARATOR CABINET), MUST BE INCLUDED WITH THE SIEMENS ORDER. THE PIPE SIZE BETWEEN THE KRAUS CHILLER AND INTERFACE PANEL, OR BETWEEN THE WATER SUPPLY AND SEP MUST BE 2 INCH UP TO 82 FEET, 2-1/2 INCH UP TO 148 FEET, CONSULT FOR LONGER PIPE. PERMISSIBLE MATERIALS THAT CAN BE USED FOR THE PIPING ARE: STAINLESS STEEL (V2A, V4A), NON-FERROUS METAL (COPPER, BRASS), SYNTHETIC MATERIAL, PLASTICS, BRAZING SOLDER, HARD SOLDER, OR FITTING SOLDER TYPE 3 AND 4. THERE ARE MATERIALS THAT MAY CAUSE DAMAGE TO THE COOLING SYSTEM AND CANNOT BE USED, THESE MATERIALS ARE ALUMINUM, IRON, CARBON STEEL, ZINC, ZINC PLATED STEEL, OR STANDARD STEEL PIPES.

THESE REQUIREMENTS ARE REQUIRED FOR NEW INSTALLATIONS, IF EXISTING WATER PIPES COMPLY WITH SIEMENS WATER SPECIFICATIONS, THEY DO NOT NEED TO BE REPLACED.

NORMAL TAP WATER MUST BE AVAILABLE FOR FILLING THE SECONDARY WATER CIRCUIT. THERE SHALL BE A HOSE BIB LOCATED WITHIN 65' OF THE SEP, IFF, ACC OR THE KRAUS CHILLER.

THE SUPPLY AND RETURN CHILLED WATER PIPES MUST BE LABELED. THE LOCATION OF THE LABELS MUST BE AT ALL CONNECTION AND REFILLING POINTS AND MUST CONTAIN FLOW DIRECTION AND CONTENTS.

ENVIRONMENTAL REQUIREMENTS

1) AIR CONDITIONING IS TO PROVIDE A TEMPERATURE OF 70°F ±5°F IN THE EXAM ROOM, 70°F±10°F IN THE EQUIPMENT & CONTROL AREAS, RELATIVE HUMIDITY OF 40-60% (NON-CONDENSING) IS REQUIRED EXAMINATION ROOM AND 40-80% (NON-CONDENSING) IN ALL OTHER AREAS WHERE SIEMENS EQUIPMENT IS INSTALLED. THESE CONDITIONS ARE TO BE MET AT ALL TIMES; 24 HOURS A DAY, 7 DAYS A WEEK.

2) A DEDICATED AIR CONDITIONING AND HUMIDIFICATION SYSTEM IS RECOMMENDED FOR THE EXAM ROOM. A MINIMUM FRESH AIR EXCHANGE RATE OF 6 TIMES PER HOUR FOR THE EXAM ROOM IS REQUIRED. AIR SUPPLY AND RETURN ABOVE THE FINISHED CEILING IN THE EXAM ROOM IS RECOMMENDED. EACH ROOM SHOULD HAVE A DEDICATED CONTROL AND SENSOR TO MONITOR AND ADJUST THE AIR.

3) THE HEAT INTO THE EXAM ROOM IS LESS THAN 10,236 BTU/HR. THE HEAT INTO THE EQUIPMENT ROOM IS 17,060 MAXIMUM, BUT IS TYPICALLY 8,530 BTU/HR. THIS HEAT DISSIPATION IS FROM THE SIEMENS EQUIPMENT ONLY. AUXILIARY SUPPORT EQUIPMENT (i.e UPS) AND LIGHTING MUST BE CONSIDERED FOR TOTAL HEAT LOADS.

4) IT IS IMPORTANT FOR FRESH AIR INTAKE SYSTEMS TO EXHAUST AIR DIRECTLY OUT OF THE BUILDING. THE EXHAUST AIR MUST NOT BE DEFLECTED INTO ANOTHER ROOM. THE MAGNET ROOM EXHAUST AIR SHOULD BE INSTALLED AT LEAST 6'-6" ABOVE FINISHED FLOOR.

5) THE AIR INTAKE OF THE AIR CONDITIONING SYSTEM MUST NOT BE LOCATED IN THE VICINITY OF THE QUENCH VENT EXHAUST.

6) IF THE INPUT DRAWS UPON AIR FROM OUTSIDE THE BUILDING, IT IS RECOMMENDED TO INSTALL AN ON-SITE FILTER TO REMOVE DUST PARTICLES GREATER THAN 10 MICRONS.

CHILLED WATER REQUIREMENTS

WATER REQUIREMENTS TO BE MEASURED AT THE SEP CABINET.

| | |
|----------------------------------|-------------------------------|
| FLOW RATE: | 23.78-29.05 GPM |
| WATER TEMPERATURE: | 48°F ±4°F |
| BTU DISCHARGE TO THE WATER | 204,729 BTU/HR |
| WATER PRESSURE | MAXIMUM 87 PSI |
| LOSS OF PRESSURE FOR SEP CABINET | 14.5 PSI MAXIMUM |
| CHILLED WATER ACIDITY RANGE | 6 pH TO 8 pH |
| CHILLED WATER HARDNESS | <250 ppm CALCIUM CARBONATE |
| CHLORINE GAS CONCENTRATION | <200 ppm |
| FILTRATION | 500 μm |

FOR INSTALLATION OF A KRAUS KSC 215 CHILLER, IT IS THE RESPONSIBILITY OF THE CUSTOMER/MECHANICAL CONTRACTOR TO PROVIDE A MIXTURE OF WATER WITH 35%-38% ETHYLENE GLYCOL PRIOR TO CHILLER START UP. DO NOT USE PROPYLENE GLYCOL OR AUTOMOTIVE ANTI-FREEZE.

THE AMOUNT OF THE MIXTURE MUST FILL THE CHILLER, MR SYSTEM AND PIPING (SUPPLY AND RETURN), SEE EXAMPLES BELOW.

(1) GALLON OF UNDILUTED GLYCOL, OR (2) GALLONS OF WATER/GLYCOL MIXTURE MUST REMAIN ON SITE FOR USE AFTER START UP.

MIXTURE VOLUME INCLUDING SUPPLY & RETURN+15 GAL. CHILLER & MR

| PIPE DIAMETER | TOTAL LENGTH | MIXTURE VOLUME | GLYCOL NEEDED |
|---------------|--------------|----------------|---------------|
| 2" | 100' | 31.3 GALLONS | 11.9 GALLONS |
| 2" | 200' | 47.6 GALLONS | 18.1 GALLONS |
| 2.5" | 100' | 40.5 GALLONS | 15.4 GALLONS |
| 2.5" | 200' | 66.0 GALLONS | 25.1 GALLONS |

MIXTURE VOLUME = $3.14 \times (\text{PIPE RADIUS})^2 \times \text{PIPE LENGTH} + 15$ GALLONS.
GLYCOL AMOUNT = 35-38% OF MIXTURE VOLUME.

QUENCH VENT NOTES

LIQUID AND GASEOUS HELIUM ARE USED IN THE OPERATION OF A SUPERCONDUCTING MRI SYSTEM. THE MECHANICAL CONTRACTOR SHALL PROVIDE A VENT, ACCORDING TO SIEMENS SPECIFICATIONS, TO EXHAUST GASEOUS HELIUM FROM THE MAGNET TO OUTSIDE THE BUILDING. PLEASE SEE THE SIEMENS TYPICAL DRAWINGS FOR DETAILS.

SIEMENS

FOR REFERENCE ONLY,
NOT FOR CONSTRUCTION.

MAGNETOM AERA 1.5T SPECIFICATIONS

PROTECTING THE ENVIRONMENT

PROTECTING THE IMMEDIATE ENVIRONMENT FROM THE EFFECT OF THE MAGNETIC FIELD REQUIRES CONSIDERATION. INFORMATION STORED ON MAGNETIC DATA CARRIERS SUCH AS DISKS, TAPES, AND CREDIT CARDS MAY BE ERASED IF IN CLOSE PROXIMITY. CAUTION WITH REGARD TO HEART PACEMAKERS MUST BE EXERCISED. MOST PACEMAKER UNITS EMPLOY A REED RELAY WHICH MAY CHANGE OPERATING MODE WHEN EXPOSED TO AN EXTERNAL MAGNETIC FIELD. THEREFORE, PACEMAKER USERS MUST BE KEPT AT A SPECIFIED DISTANCE FROM THE MAGNET WHICH IS DETERMINED BY THE MAGNETIC FIELD STRENGTH.

PROTECTING THE MAGNETIC FIELD

THE SIEMENS MAGNETOM UTILIZES A SUPERCONDUCTIVE MAGNET WITH AN EXTREMELY HOMOGENEOUS FIELD WITHIN THE MAGNET TO PROVIDE DISTORTION-FREE IMAGING. THE PRESENCE OF FERROMAGNETIC MATERIAL WITHIN THE VICINITY OF THE MAGNET CAN ADVERSELY AFFECT THE UNIFORMITY OF THE USEFUL MAGNETIC FIELD. THIS APPLIES TO STATIONARY FERROUS MATERIAL (STRUCTURAL STEEL) WHICH IS TO BE MINIMIZED. STATIONARY STEEL COMPENSATION MAY BE ACHIEVED BY MAGNET POSITIONING AND SELECTIVE USE OF SHIMS. FIELD DISTORTION ENCOUNTERED BY MOVING FERROMAGNETIC OBJECTS IS MORE DIFFICULT TO COMPENSATE AND MAY REQUIRE THE USE OF MAGNETIC SHIELDING.

MAGNETIC FRINGE FIELDS

MAGNETIC FIELDS MAY AFFECT THE FUNCTION OF DEVICES IN THE VICINITY OF THE MAGNET. THESE DEVICES MUST BE OUTSIDE CERTAIN MAGNETIC FIELDS. THE DISTANCES LISTED ARE FROM THE MAGNET ISOCENTER AND DO NOT CONSIDER ANY MAGNETIC ROOM SHIELDING.

| X/Y AND Z AXIS | DEVICES |
|---------------------------|--|
| 6'-1" / 9'-2" 3.0mT | SMALL MOTORS, WATCHES, CAMERAS, CREDIT CARDS, MAGNETIC DATA CARRIERS (SHORT-TERM EXPOSURE) |
| 7'-3" / 11'-6" 1.0mT | COMPUTERS, MAGNETIC DISK DRIVES, OSCILLOSCOPES, PROCESSORS |
| 8'-3" / 13'-2" 0.5mT | CARDIAC PACEMAKERS, X-RAY TUBES, INSULIN PUMPS, B/W MONITORS, MAGNETIC DATA CARRIERS (LONG-TERM STORAGE) |
| 9'-9" / 16'-1" 0.2mT | SIEMENS CT SCANNERS |
| 10'-4" / 17'-1" 0.15mT | COLOR MONITORS, SIEMENS LINEAR ACCELERATORS |
| 13'-1" / 22'-3" 0.05mT | X-RAY IMAGE INTENSIFIERS, GAMMA CAMERAS, PET/CYCLOTRON, ELECTRON MICROSCOPES, LINEAR ACCELERATORS |

THE OWNER/USER IS TO VERIFY THE LOCATION OF THE 0.5mT FIELD AND ENSURE THAT IT IS MAINTAINED AS A RESTRICTED AREA.

MAGNET SITING REQUIREMENTS

IT MUST BE ENSURED THAT THE MAGNET IS LOCATED SO THAT THE STABILITY AND HOMOGENEITY OF THE MAGNETIC FIELD ARE NOT ADVERSELY AFFECTED BY EXTRANEOUS FIELDS AND STATIC OR DYNAMIC FERROMAGNETIC OBJECTS.

| X/Y AND Z AXIS | SOURCE OF INTERFERENCE |
|-----------------|--|
| 3'-6" | STEEL REINFORCEMENT RODS IN FLOOR - MAXIMUM 20 LBS/SQ. FT. |
| 18'-1" / 21'-4" | STRETCHERS UP TO 110 LBS. |
| 13'-1" | A/C CHILLERS |
| 19'-9" / 23'-0" | TRANSPORT DEVICES UP TO 440 LBS. |
| 21'-4" / 26'-3" | VEHICLES UP TO 2,000 LBS. |
| 23'-0" / 31'-3" | ELEVATORS, TRUCKS UP TO 10,000 LBS. |
| 39'-4"/26'-2" | AC TRANSFORMERS LESS THAN 100 KVA |
| 41'-0"/32'-9" | AC TRANSFORMERS LESS THAN 250 KVA |
| 42'-7"/39'-4" | AC TRANSFORMERS LESS THAN 650 KVA |
| 45'-11"/49'-3" | AC TRANSFORMERS LESS THAN 1600 KVA |
| 9'-10"/6'-6" | AC CABLES, MOTORS LESS THAN 100 AMPS |
| 22'-11"/9'-10" | AC CABLES, MOTORS LESS THAN 250 AMPS |
| 131'-2" | ELECTRIC RAILWAY SYSTEMS |

FOR IRON OBJECTS LOCATED UP TO 45' FROM THE Z AXIS, THE DISTANCES FOR THE Z AXIS MUST BE USED. REDUCTION IS POSSIBLE WITH STEEL SHIELDING.

MAXIMUM CABLE LENGTH

THERE ARE 3 DIFFERENT LENGTHS OF CABLE THAT ARE AVAILABLE FOR THE MRI SYSTEM DIFFERENTIATED BY MAXIMUM LENGTHS FROM THE MAGNET TO THE FILTER PANEL (INSIDE) AND FROM THE FILTER PANEL TO THE ELECTRONICS (OUTSIDE).

| INSIDE | OUTSIDE |
|--------|---------|
| 20' | 4' |
| 20' | 32' |
| 20' | 39' |

THE VERTICAL DISTANCE FOR CABLE TRAVEL FROM THE FILTER PANEL TO THE CABLE TRAY, AND FROM THE CABLE TRAY TO THE MAGNET MUST BE CONSIDERED.

THE MAXIMUM DISTANCE FROM THE ACC CABINET TO THE CONTROL CONSOLE IS 75 FEET.

SIEMENS

FOR REFERENCE ONLY,
NOT FOR CONSTRUCTION.

MAGNETOM AERA 1.5T SPECIFICATIONS

RF SHIELDING

THE EXAMINATION AREA MUST BE SHIELDED TO PROVIDE A REDUCTION OF RADIO FREQUENCY WAVES EMANATING FROM EXTERNAL TRANSMITTERS. THE REQUIRED ATTENUATION IS 90dB IN THE FREQUENCY RANGE OF 15-128 MHz. IF CO-SITING TWO SYSTEMS EACH ROOM SHOULD BE 100 dB. THE RF SHIELD MUST BE TESTED BEFORE AND AFTER MAGNET PLACEMENT IN THE RF ROOM AND AFTER THE SIEMENS RF FILTER PANEL IS INSTALLED.

THE RF-SHIELDING MUST BE INSULATED FROM ALL GROUNDS SUCH THAT THE ONLY GROUND IS THE SINGLE POINT GROUND ON THE OUTSIDE OF THE RF-ROOM WALL. RESISTANCE \geq 100 OHMS.

ALL ELECTRICAL LINES INTO THE RF ROOM MUST BE ROUTED THROUGH RF FILTERS (PROVIDED BY RF SHIELDING SUPPLIER).

ALL ELECTRICALLY NON-CONDUCTIVE SUPPLY LINES (E.G. FIBER OPTIC CABLES, OR HOSES) INTO THE RF ROOM MUST BE ROUTED THROUGH RF SEALED WAVEGUIDES (PROVIDED BY RF SHIELDING SUPPLIER).

FOR PRESSURE EQUALIZATION PURPOSES THE RF DOOR SHOULD OPEN TO THE OUTSIDE OF THE RF ROOM. AS AN ALTERNATIVE A 24"x24" OPENING IN THE RF ROOM FOR PRESSURE EQUALIZATION IS REQUIRED.

BUILDING VIBRATIONS

VIBRATION OF THE SITE HAS THE ABILITY TO AFFECT THE STABILITY AND HOMOGENEITY OF THE MAGNETIC FIELD. THEREFORE EXTERNAL VIBRATIONS OR SHOCKS AFFECTING THE MAGNET MAY DEGRADE IMAGE QUALITY. IN THE THREE SPATIAL ORIENTATIONS THE BUILDING MUST NOT EXCEED ACCELERATION OF 0.001m/s or -80dB(g) $g=9.81$ m/s

THE REQUIREMENT FOR a_{max} IS MEASURED AS MAXIMUM RMS VALUE PER FREQUENCY COMPONENT <0.5 Hz IN THE FOURIER TRANSFORMATION OF THE RECORDED SIGNAL (SPECTRUM).

THE VIBRATION LEVEL OF CONTINUOUS VIBRATIONS (CAUSED BY AIR CONDITIONER, COMPRESSOR, ETC.) AT THE LOCATION OF THE MAGNET MUST NOT EXCEED THE SPECIFIED VALUES.

FOR ALL NON-CONTINUOUS TRANSIENT VIBRATIONS THE FIGURES SHOULD BE MULTIPLIED BY 4 (OR 12dB).

CONTACT SIEMENS PROJECT MANAGER FOR MORE DETAILS.

TRANSPORTING REQUIREMENTS

LARGEST ITEM - MAGNET - 9,566 LBS.

MAGNET DIMENSIONS WITH TRANSPORT DEVICE:
7'-7" HIGH x 7'-7" WIDE x 5'-2" DEEP (6'-0" WITH SUPPORT)
THE ROOF HATCH/DELIVERY OPENING SHOULD BE 4" LARGER.

TO TRANSPORT THE GPA/EPC CABINET (3,307 POUNDS) A MINIMUM ROOM HEIGHT OF 6'-9" IS REQUIRED. 6'-3" WITH WHEELS REMOVED. 6'-1" WITH WHEELS AND MAINS CONNECTION REMOVED.